

# 4FRONT™

The 4Front™ is a quiet, more responsive front-wheel drive power chair that features automotive-grade suspension with unprecedented comfort and ride quality. The 4Front's CASE (Caster Angle Sensor Encoder) technology and SRS (Smooth Ride Suspension) provide stability for smooth obstacle transitions and unmatched driving performance. The 4Front™ features standard speeds of up to 6 mph in the lowered position, and speeds up to 3.2 mph while elevated at 10" with our optional safe seat elevation. Standard LED fender lights allow you to see and be seen.



QUANTUM® 4FRONT™ shown in Electric Blue with optional *Safe Seat Elevation* and Q-Logic 3 Advanced Drive Controls with Bluetooth®



Dedicated to revolutionizing the power chair experience.

QUANTUM®

(US) 866-800-2002

| (CAN) 888-570-1113

| [quantumrehab.com](http://quantumrehab.com)



## Specifications

<b>DRIVE WHEELS</b>	14" flat free
<b>CASTER WHEELS</b>	Rear: 9" flat free
<b>ANTI-TIP WHEELS</b>	Front: 4"
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 6 mph
<b>GROUND CLEARANCE</b>	3" (motor pivot)
<b>TURNING RADIUS</b> <sup>2</sup>	24.5"
<b>OVERALL LENGTH</b> <sup>2</sup>	40"
<b>BASE WIDTH</b>	24"
<b>SEATING SIZES</b>	
<b>TRU-Balance® 3 Static Seat</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU-Balance® 3 Lift &amp; Tilt</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>SEAT-TO-FLOOR HEIGHTS</b> <sup>9</sup>	
<b>TRU-Balance® 3 Static Seat</b> <sup>5</sup>	15.25-17.25"
<b>TRU-Balance® 3 Power Tilt</b> <sup>5</sup>	17.5-19.5"
<b>TRU-Balance® 3 Lift &amp; Tilt</b> <sup>5</sup>	17.75-19.75"
<b>MANUAL RECLINE</b>	Yes
<b>BATTERY SIZE</b> <sup>3</sup>	Group 24 or Group 34
<b>BATTERY WEIGHT</b> <sup>6</sup>	Group 24: 51.5 lbs. Group 34: 38.5 lbs.
<b>AVAILABLE ELECTRONICS</b> <sup>7</sup>	120A Q-Logic 3 EX
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	Two-pole, front-wheel drive
<b>WEIGHT CAPACITY</b> <sup>4</sup>	300 lbs.
<b>BASE WEIGHT</b>	140 lbs.
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 7) All Quantum Rehab® units can be ordered with specialty controls.
- 8) Call for availability
- 9) Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front's suspension system, seat-to-floor heights are approximately 0.25" lower when the unit is loaded.

## Standard Features

- CASE (Caster Angle Sensor Encoder) technology
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Q-Logic 3 EX Advanced Drive Controls
- LED front and rear fender lights
- Quiet ride technology



LED fender lights

## Options

- 10" of power adjustable safe seat elevation at speeds of up to 3.2 mph
- TRU-Comfort 2 Seating
- Trilogy vent tray <sup>8</sup>
- TRU-Balance® 3 Seating
- Compact vent tray <sup>8</sup>

## Colors

SEAT BACKS for units without TRU-Balance® 3 are only available in *Back in Black*.



Pretty in Pink



Lime Light



Yellow Blaze



Deep Purple



Orange Crush



Electric Blue



Red Corvette



Midnight Blue



White Satin



Night Train



Silver Bullet



Grey Street



Back in Black



Military Green

**QUANTUM®**

(US) 866-800-2002

(Canada) 888-570-1113

quantumrehab.com

