



Q6 EDGE® 2.0
shown in *Green Machine* with
12" iLevel® technology and
Q-Logic 3 Advanced Drive
Controls with Bluetooth®

FDA Class II Medical Device *

Q6 edge[®] 2.0

With features like standard 6.25 mph motors and advanced responsiveness, the Q6 Edge® 2.0 is engineered to meet the performance needs of the most active user. The Q6 Edge 2.0 accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel® technology.

QUANTUM[®]

Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic 14" Knobby
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Up to 6.25 mph ¹⁰
GROUND CLEARANCE	2.7" (battery tray)
TURNING RADIUS ²	20.5"
OVERALL LENGTH ²	32.23"
BASE WIDTH	24.4"
SEATING SIZES	
Synergy[®]/TRU-Balance[®] 3 Static Seat	Width: 10-22" Depth: 10-22"
TRU-Balance[®] 3 Power Tilt	Width: 12-22" Depth: 12-22"
TRU-Balance[®] 3 Tilt & iLevel[®]	Width: 12-22" Depth: 12-22"
TRU-Balance[®] Flex Seat	Width: 15-22" Depth: 15-20"
SEAT-TO-FLOOR HEIGHTS	
Synergy[®] Seat ⁵	15.875-20.875"
TRU-Balance[®] 3 Static Seat ⁵	16.5-18.5"
TRU-Balance[®] 3 Power Tilt ⁵	17.5-19.5"
TRU-Balance[®] 3 Tilt & iLevel[®] ⁵	17.5-19.5"
TRU-Balance[®] Flex Seat ⁵	16.5-18.25"
MANUAL RECLINE ⁶	No
BATTERY SIZE ³	NF-22
BATTERY WEIGHT ⁷	38 lbs.
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 EX 120A Q-Logic 3e 75A Q-Logic NE 75A Q-Logic NE+ 75A Q-Logic NES
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole Mid-Wheel 6 [®]
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	138.6 lbs.
MEDICARE CODE COVERAGE	Group 3: Standard (K0848/K0849), Single Power (K0856/K0857), Multiple Power (K0861)

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
2) Without front riggings
3) Battery capacity may change when adding power positioning systems.
4) Weight capacity is dictated by the selection of seating.
5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.
Listed measurements are in the lowest position; does not include elevating system unless specified.
6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
7) Battery weight may vary +/- 3 lbs. based on manufacturer.
8) All Quantum Rehab[®] units can be ordered with specialty controls.
9) Call for availability
10) 6.25 mph speed available only with Q-Logic 3.

Features

- Standard 6.25 mph motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics
- iLevel[®] option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- iLevel is retrofittable on the Q6 Edge[®] 2.0
- Available with WC19 crash-test approved occupied transit kit option on tilt & recline and iLevel

Options

- iLevel[®] option allows 12" of seat elevation while at walking speed (3.5 mph)
- Q-Logic 3 EX Drive Controls
- TRU-Balance[®] 3 Seating
- TRU-Comfort 2 Seating
- Depth-adjustable vent tray⁹
- Compact vent tray⁹
- 4" wide flat-free knobby extended traction tires⁹
- LED fender lights (standard with iLevel[®])
- Accu-Trac Advanced Tracking Technology

Colors



(US) 866-800-2002

(Canada) 888-570-1113

quantumrehab.com

* Quantum Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

