



The compact J6 is the ultimate choice for clients who demand the superb stability of Mid-Wheel 6[®] Design and tight-quarter maneuverability. Highly adaptable, the J6 accepts a wide range of seating and electronics options while its powerful drivetrain delivers high-performance torque.

J6 shown with high-back seat and NE controller

PRIDE[®]
PERFORMANCE MOBILITY
— LIVE YOUR BEST[®] —

Specifications

DRIVE WHEELS	10"
OMNI-CASTERS	Front: 5" solid Rear: 6" solid
SUSPENSION	ATX Suspension (Active-Trac [®] with extra stability)
MAXIMUM SPEED ¹	Up to 4 mph
BRAKES	Intelligent braking (electronic, regenerative disc brakes)
GROUND CLEARANCE	2.25"
TURNING RADIUS ²	21.25"
OVERALL SIZE ²	Length: 34.5" Width: 23.25"
SEATING	Synergy [®] Seat Adjustable Width: 10"-20" Adjustable Depth: 10"-20" High-back Seat (20" x 20" max. dim.) Medium-back Seat (20" x 20" max. dim.) Specialty Seat (20" x 20" max. dim.)
DRIVETRAIN	2-motor, Mid-Wheel 6 [®]
BATTERIES	Two 12 volt, deep cycle U-1
RANGE PER CHARGE ¹	Up to 15 miles
BATTERY CHARGER	5A, off-board
AVAILABLE ELECTRONICS	75A, Q-Logic NE (standard) 70A, VSI 6-key (for use with power elevating seat)
WEIGHT CAPACITY	300 lbs. maximum (user weight)
BASE WEIGHT	104.5 lbs.
SEAT WEIGHT	42.5 lbs.
BATTERY WEIGHT	24.5 lbs. (each)
WARRANTY	Lifetime limited warranty on frame; 2-year warranty on electronics; 18-month warranty on drive motor

Shroud Colors



Jazzy
Red



Jazzy
Blue

FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Dimensions shown are for power base only. Overall measurements will vary based on seating and accessory selections.

Features


- Mid-Wheel 6[®] allows six wheels on the ground for maximum stability
- OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- Side-mounted, easily accessible freewheel levers
- ATX Suspension (Active-Trac[®] with extra stability) incorporates front OMNI-Casters for enhanced performance over more varied terrain
- Easy front access to batteries

Options

- Swing-away joystick mount
- Multi-axis joystick mount
- Elevating leg rests
- Oxygen tank holder
- Neoprene leg rest sleeves
- Gel adductor buttons
- Power elevating seat
- Angle-adjustable footplates
- 70° swing-away leg rests

Pride[®] Lifts & Ramps

- Single-Fold Ramps
- Solid Ramps
- Rubber Threshold Ramps



Single-Fold Ramps

Solid Ramps

These ramps are compatible with the J6.

Beamsville, Ontario
888-570-1113
www.pridemobility.com

Pride[®]
Mobility Products Co.
— LIVE YOUR BEST[®] —