



# Q6 edge®

With features like standard 4-pole motors, Mid-Wheel 6® drive design and ATX Suspension, the Q6 Edge® is engineered to meet the performance needs of the most active user. The Q6 Edge accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value.

[quantumrehab.com](http://quantumrehab.com)

CA-Q6EDGE-33116

**Q**  
QUANTUM®

Beamsville, Ontario

888-570-1113

The Rehab Power Chair  
Solutions Company

# Specifications

<b>DRIVE WHEELS</b>	Standard: 14" Solid Optional: 14" Pneumatic	
<b>CASTER WHEELS</b>	Front: 6" Solid Rear: 6" Solid	
<b>ANTI-TIP WHEELS</b>	N/A	
<b>MAXIMUM SPEED</b> <sup>4</sup>	Up to 5 mph (Optional: 6 mph)	
<b>GROUND CLEARANCE</b>	3" (frame center)	
<b>TURNING RADIUS</b> <sup>5</sup>	20"	
<b>OVERALL LENGTH</b>	35.5" without foot riggings	
<b>BASE WIDTH</b>	24"	
<b>SEATING SIZES</b>	<b>TRU-Balance® 2</b>	<b>TRU-Balance® 3</b>
<b>Synergy®/Static Seat</b>	Width: 10-22" Depth: 10-22"	Width: 12-22" Depth: 12-22"
<b>Power Tilt</b>	Width: 14-22" Depth: 14-22"	Width: 12-22" Depth: 12-22"
<b>Lift &amp; Tilt</b>	N/A	Width: 12-22" Depth: 12-22"
<b>SEAT-TO-FLOOR HEIGHTS</b>	<b>TRU-Balance® 2</b>	<b>TRU-Balance® 3</b>
<b>Synergy®/Static Seat</b> <sup>8</sup>	15.9-17.9"	16.5-18.5"
<b>Power Tilt</b> <sup>8</sup>	17.3-19.3"	17.5-19.5"
<b>Lift &amp; Tilt</b> <sup>8</sup>	N/A	18-20"
<b>MANUAL RECLINE</b> <sup>10</sup>	Yes	
<b>BATTERY SIZE</b> <sup>6</sup>	NF-22	
<b>BATTERY WEIGHT</b> <sup>14</sup>	38 lbs.	
<b>AVAILABLE ELECTRONICS</b> <sup>12</sup>	75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX	
<b>BATTERY CHARGER</b>	8A Off-board	
<b>MOTOR PACKAGES</b>	4-Pole 5 mph <sup>1</sup> or 6 mph <sup>2</sup>	
<b>WEIGHT CAPACITY</b> <sup>7</sup>	300 lbs.	
<b>BASE WEIGHT</b>	138 lbs.	

- 1) Standard
- 2) Optional
- 3) Requires HD Suspension
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems.
- 7) Weight capacity is dictated by the selection of seating.
- 8) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 9) 250 lbs. weight capacity with TRU-Balance® Systems
- 10) Maximum dimensions of the Synergy® Manual Recline vary by power base and seating system. Call for details.
- 11) Requires Group 34 batteries
- 12) All Quantum Rehab® units can be ordered with specialty controls
- 13) Required with high-speed motors
- 14) Battery weight may vary +/- 3 lbs. based on manufacturer.

# Features

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® 3 Power Positioning
- Low impact OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters and semi-independent rear caster beam for enhanced performance over more varied terrain

# Options

- TRU-Balance® 3 Lift & Tilt
- TRU-Comfort 2 Seating System
- Q-Logic 2 EX Drive Controls
- 6 mph motor package
- Accu-Trac Advanced Tracking Technology
- TRU-Balance® 3 (10") Lift Trilogly vent tray
- Compact vent tray
- LED full lighting package
- 4" wide flat-free knobby extended traction tires
- WC19 transit securement package

# Colors



Beamsville, Ontario • 888-570-1113 • [www.quantumrehab.com](http://www.quantumrehab.com) • [www.ilevel.rehab](http://www.ilevel.rehab)

