

i L E V E L R E A D Y

AVAILABLE  
WITH 12" iLEVEL®  
POWER SEAT  
ELEVATION

⋮

OPTIONAL  
Q-LOGIC 3  
ADVANCED  
DRIVE  
CONTROLS

⋮

LED FENDER  
LIGHTS COME  
STANDARD  
WITH iLEVEL®



**Q6 EDGE® 2.0**  
shown in *Lime Light* with  
iLevel® technology and  
Q-Logic 3 drive controls

# Q6 edge® 2.0

With features like standard 6 mph motors and advanced responsiveness, the Q6 Edge® 2.0 is engineered to meet the performance needs of the most active user. The Q6 Edge 2.0 accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel® technology.

QUANTUM®

## Specifications

<b>DRIVE WHEELS</b>	Standard: 14" Solid Optional: 14" Pneumatic 14" Knobby
<b>CASTER WHEELS</b>	Front: 6" Solid (Optional: 8" Solid) Rear: 6" Solid (Optional: 8" Solid)
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 6 mph
<b>GROUND CLEARANCE</b>	2.7" (battery tray)
<b>TURNING RADIUS</b> <sup>2</sup>	20.5" with 6" casters 23.25" with 8" casters
<b>OVERALL LENGTH</b> <sup>2</sup>	35.5" with 6" casters 40.48" with 8" casters
<b>BASE WIDTH</b>	24.4"
<b>SEATING SIZES</b>	
<b>Synergy®/TRU-Balance® 3 Static Seat</b>	<b>Width:</b> 10-22" <b>Depth:</b> 10-22"
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU-Balance® 3 Tilt &amp; iLevel®</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>Synergy®/TRU-Balance® 3 Static Seat</b> <sup>5</sup>	16.5-18.5"
<b>TRU-Balance® 3 Power Tilt</b> <sup>5</sup>	17.5-19.5"
<b>TRU-Balance® 3 Tilt &amp; iLevel®</b> <sup>5</sup>	17.5-19.5"
<b>MANUAL RECLINE</b> <sup>6</sup>	Yes
<b>BATTERY SIZE</b> <sup>3</sup>	NF-22
<b>BATTERY WEIGHT</b> <sup>7</sup>	38 lbs.
<b>AVAILABLE ELECTRONICS</b> <sup>8</sup>	120A Q-Logic 3 EX 75A Q-Logic 2 NE 75A Q-Logic 2 NE+
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-Pole Mid-Wheel 6®
<b>WEIGHT CAPACITY</b> <sup>4</sup>	300 lbs.
<b>BASE WEIGHT</b>	138.6 lbs. with 6" casters 148 lbs. with 8" casters
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability

## Features

- Standard 6 mph motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- iLevel is retrofittable on the Q6 Edge® 2.0

## Options

- iLevel® option allows 12" of seat elevation while at walking speed (3.5 mph)
- Q-Logic 3 EX Drive Controls
- TRU-Balance® 3 Seating
- TRU-Comfort 2 Seating
- Trilogy vent tray<sup>9</sup>
- Compact vent tray<sup>9</sup>
- 4" wide flat-free knobby extended traction tires<sup>9</sup>
- LED fender lights (standard with iLevel®)
- 12" x 12" seat sizes with iLevel (SOF only)
- 8" solid casters
- Accu-Trac Advanced Tracking Technology

## Colors

