

EDGE[®] 3

4.5 MPH iLEVEL[®] READY


NOW AVAILABLE

MAXIMUM SPEED UP TO 7 MPH
with new optional 4-pole motors



FDA Class II
Medical Device*

Upgraded SRS (Smooth Ride Suspension) provides a smooth, comfortable ride and advanced stability. Optional iLevel[®] is available at 4.5 mph with up to 12" of power adjustable seat height. Standard features include front and rear LED fender lights, a USB charger, and drive wheel color accents. TRU Balance[®] 4 Power Positioning with anterior tilt and memory seating.



RASPBERRY BERET

Q
QUANTUM[®]
#1 FOR REHAB POWER



quantumrehab.com

Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic; 14" Knobby
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Standard: Up to 6.25 mph ⁶ Optional: Up to 7 mph ⁹ (4.5 mph with iLevel [®])
GROUND CLEARANCE	2.7" (battery tray)
TURNING RADIUS ²	20.5"
OVERALL LENGTH ²	35.5"
BASE WIDTH	24.4"
SEATING SIZES	
TRU Balance[®] 4 Anterior Tilt	Width: 16-22" Depth: 14-22"
Synergy[®] Seat	Width: 10-22" Depth: 10-22"
TRU Balance[®] 3 Static Seat	Width: 12-22" Depth: 12-22"
TRU Balance[®] 3 Power Tilt TRU Balance[®] 4 Power Tilt	Width: 12-22" Depth: 12-22"
TRU Balance[®] 3 Tilt and iLevel[®]	Width: 12-22" Depth: 12-22"
TRU Balance[®] Flex Seating	Width: 15-22" Depth: 15-20"
SEAT-TO-FLOOR HEIGHTS ¹⁰	
TRU Balance[®] 4	17.5-19.25"
Synergy[®] Seat	15.875-20.875"
TRU Balance[®] 3 Static Seat ⁵	16.5-18.5"
TRU Balance[®] 3 Power Tilt ⁵ TRU Balance[®] 4 Power Tilt ⁵	17.5-19.5"
TRU Balance[®] 3 Tilt & iLevel[®] ⁵	17.5-19.5"
TRU Balance[®] Flex Seating	16.5 - 18.25"
BATTERY SIZE ³	NF-22 (2)
BATTERY WEIGHT ⁷	40 lbs. each
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 75A Q-Logic NES 90A Q-Logic NES 75A Q-Logic NE
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-pole and 4-pole
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	138.6 lbs.
MEDICARE CODE COVERAGE	Group 3: Standard (K0848/K0849), Single Power (K0856/K0857), Multiple Power (K0861)

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) 6.25 mph speed available only with Q-Logic 3.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab[®] units can be ordered with specialty controls.
- 9) Up to 7 mph speed available only with 4-pole motor option and Q-Logic 3 controller.

*The Quantum[®] Edge 3 is an FDA Class II Medical Device designed to aid individuals with mobility impairments.

Features

- Upgraded SRS (Smooth Ride Suspension)
- **Up to 7 mph** with new optional 4-pole motors ⁹
- Up to 4.5 mph with iLevel[®]
- Fender lights and USB charger included
- Made in America of foreign and domestic parts
- Available with optional WC19 crash-test approved occupied transit kit
- Available with optional Accu-Trac Advanced Tracking Technology

Colors



Hydrographics



PROUDLY MADE IN
AMERICA
OF FOREIGN AND DOMESTIC PARTS

QUANTUM[®]

#1 FOR REHAB POWER

(US) 866-800-2002

| (CAN) 888-570-1113

| quantumrehab.com

