

REHAB POWER CHAIR SOLUTIONS



Q6 EDGE® X

shown in *Silver* with TRU-Balance® 3
Power Positioning



Q6 edge®



With features like standard 4-pole motors, Mid-Wheel 6® Drive Design and ATX Suspension, the Q6 Edge® X is engineered to meet the performance needs of the most active user. The Q6 Edge X accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value.

QUANTUM®

Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic	
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid	
ANTI-TIP WHEELS	N/A	
MAXIMUM SPEED ⁴	Up to 6 mph	
GROUND CLEARANCE	3" (frame center)	
TURNING RADIUS ⁵	20"	
OVERALL LENGTH	35.5" without foot riggings	
BASE WIDTH	24"	
SEATING SIZES	TRU-Balance[®] 2	TRU-Balance[®] 3
Synergy[®]/Static Seat	Width: 10-22" Depth: 10-22"	Width: 12-22" Depth: 12-22"
Power Tilt	Width: 14-22" Depth: 14-22"	Width: 12-22" Depth: 12-22"
Lift & Tilt	N/A	Width: 12-22" Depth: 12-22"
SEAT-TO-FLOOR HEIGHTS	TRU-Balance[®] 2	TRU-Balance[®] 3
Synergy[®]/Static Seat ⁸	15.9-17.9"	16.5-18.5"
Power Tilt ⁸	17.3-19.3"	17.5-19.5"
Lift & Tilt ⁸	N/A	18-20"
MANUAL TILT/RECLINE ¹⁰	Yes	
BATTERY SIZE ⁶	NF-22	
BATTERY WEIGHT ¹⁴	38 lbs.	
AVAILABLE ELECTRONICS ¹²	75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX	
BATTERY CHARGER	8A Off-board	
MOTOR PACKAGES	4-Pole 6 mph	
WEIGHT CAPACITY ⁷	300 lbs.	
BASE WEIGHT	138 lbs.	
MEDICARE CODE COVERAGE	Group 3: Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)	

- 1) Standard
- 2) Optional
- 3) Requires HD Suspension
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems.
- 7) Weight capacity is dictated by the selection of seating.
- 8) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 9) 250 lbs. weight capacity with TRU-Balance[®] Systems
- 10) The Synergy[®] Manual Tilt weight capacity is 75 lbs. and maximum dimensions are 14" W x 14" D. Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 11) Requires Group 34 batteries
- 12) All Quantum Rehab[®] units can be ordered with specialty controls
- 13) Required with high-speed motors
- 14) Battery weight may vary +/- 3 lbs. based on manufacturer.

Features

- 6 mph motor package
- Attendant control and mounting brackets
- NE Plus joystick (can control up to two power functions)
- Synergy[®] Seating System
- Two egg switches (can be utilized as an alternate method for access on/off and mode or tilt functions when a tilt is ordered)

Options

- Captain seating
- TRU-Balance[®] 3 Power Positioning Systems
- TRU-Balance 3 Static Seating
- Choice of two complimentary accessories with TRU-Balance 3 seating
- Available with TRU-Balance 3 seating:
 - Complimentary USB charger, remote stop switch, auto style seat belt, and glove box
 - Choice of complimentary push handles or transfer handles
 - Choice of complimentary centermount adjustable foot platform or adjustable angle swing-away footrests

Colors



QUANTUM[®]