

End-to-End Rehab Solutions



	Q6 EDGE® 2.0 <small>AVAILABLE WITH iLEVEL®</small>	4FRONT™ <small>AVAILABLE WITH 10" LIFT</small>	Q6 EDGE®	Q6 EDGE® 2.0X	Q6 EDGE® Z <small>AVAILABLE WITH iLEVEL®</small>	QUANTUM® 610	Q6 EDGE® HD <small>AVAILABLE WITH iLEVEL® (400 LBS. WT. CAPACITY)</small>	QUANTUM® RIVAL®	J6	QUANTUM® 1450
DRIVE CONFIGURATION	Mid-wheel Drive	Front-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Mid-wheel Drive	Rear-wheel Drive	Mid-wheel Drive	Front-wheel Drive
DRIVE WHEELS	14" Solid Optional: 14" Pneumatic; 14" Knobby	14" Flat Free	14" Solid Optional: 14" Pneumatic	14" Solid Optional: 14" Pneumatic; 14" Knobby	14" Run Flat Optional: 14" Pneumatic	10" Solid Optional: 10" Pneumatic	14" Run Flat Optional: 14" Pneumatic	14" Solid Optional: 14" Pneumatic	10" Solid Optional: 10" Pneumatic	14" Solid
CASTER WHEELS	Front: 6" Solid (Optional: 8") Rear: 6" Solid (Optional: 8")	Rear: 9" Flat Free	Front: 6" Solid Rear: 6" Solid	Front: 6" Solid (Optional: 8") Rear: 6" Solid (Optional: 8")	Front: 6" Solid Rear: 6" Solid	Front: 5" Solid Rear: 6" Solid	Front: 6" Solid Rear: 6" Solid	Front: 9" Pneumatic	Front: 5" Solid Rear: 6" Solid	9" Solid
ANTI-TIP WHEELS	N/A	Front: 4"	N/A	N/A	N/A	N/A	N/A	N/A	N/A	4" Solid
MAXIMUM SPEED ⁴	Up to 6 mph	Up to 6 mph	Up to 5 mph (Optional: 6 mph)	Up to 6 mph	Up to 6 mph (Optional: 8 mph) ¹³	Up to 4 mph	Up to 4.7 mph	Up to 6 mph (Optional: 8 mph)	Up to 4 mph	Up to 4.5 mph
GROUND CLEARANCE	2.7" (battery tray)	3" (motor pivot)	3" (frame center)	2.7" (battery tray)	3" (frame base)	2.5" (frame)	3" (frame center)	3.25" (frame base)	2.5" (frame)	1.5" (anti-tip wheels)
TURNING RADIUS ⁵	6" casters: 20.5" 8" casters: 23.25"	24.5"	20"	6" casters: 20.5" 8" casters: 23.25"	19.75"	22"	19.75"	35.75"	22"	26.5"
OVERALL LENGTH ¹¹	6" casters: 35.5" 8" casters: 40.48"	40"	35.5"	6" casters: 35.5" 8" casters: 40.48"	35.6"	34.74"	35.6"	36"	34.74"	38.25"
BASE WIDTH	24.4"	24"	24"	24.4"	26.25"	23.19"	26.25"	25"	23.19"	29.25"
Seating Sizes:	TRU-Balance® 3	TRU-Balance® 3	TRU-Balance® 3	TRU-Balance® 3	TRU-Balance® 3	TRU-Balance® 2	TRU-Balance® 2 HD TRU-Balance® 3 HD	TRU-Balance® 3	TRU-Balance® 2	TRU-Balance® 2 HD Tilt
SYNERGY®	W: 10-22" D: 10-22"	N/A	W: 10-22" D: 10-22"	W: 10-22" D: 10-22"	W: 10-24" D: 10-24"	W: 10-20" D: 10-20"	W: 20-28" D: 20-28" W: 20-32" D: 20-28"	W: 12-22" D: 12-22"	W: 10-20" D: 10-20"	W: 20-32" D: 20-28"
STATIC	W: 10-22" D: 10-22"	W: 12-22" D: 12-22"	W: 10-22" D: 10-22"	W: 10-22" D: 10-22"	W: 12-24" D: 12-24"	W: 10-20" D: 10-20"	W: 20-28" D: 20-28" W: 20-32" D: 20-28"	W: 12-22" D: 12-22"	W: 10-20" D: 10-20"	W: 20-32" D: 20-28"
POWER TILT	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-24" D: 14-24"	W: 14-20" D: 14-20"	W: 20-28" D: 20-28" W: 20-32" D: 20-25"	W: 12-22" ³ D: 12-22" ³	W: 14-20" D: 14-20"	W: 20-32" D: 20-28"
LIFT & TILT	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"	W: 16-22" D: 14-22"	W: 12-24" D: 12-24"	W: 16-20" D: 16-20"	N/A W: 20-24" D: 20-25" ¹⁰	W: 12-22" ³ D: 12-22" ³	N/A	N/A
Seat-to-Floor Heights ¹⁶:	TRU-Balance® 3	TRU-Balance® 3	TRU-Balance® 3	TRU-Balance® 3	TRU-Balance® 3	TRU-Balance® 2	TRU-Balance® 2 HD TRU-Balance® 3 HD	TRU-Balance® 3	TRU-Balance® 2	TRU-Balance® 2 HD Tilt
SYNERGY® / STATIC ⁸	16.5-18.5"	15.25-17.25" <small>NOT AVAILABLE WITH SYNERGY®</small>	16.5-18.5"	16.5-18.5"	16.5-18.5"	16.5-17.5"	16.5-18.5" 17.5-19.5"	16.5-18.1"	16.5-17.5"	16.5-18.5"
POWER TILT ⁸	17.5-19.5"	17.5-19.5"	17.5-19.5"	17.5-19.5"	18.125"-20.125"	16.875"	18.2-20.9" 18.5-20.5"	18-19.6"	16.875"	18.125-20.125"
LIFT & TILT ⁸	17.5-19.5"	17.75-19.75"	18-20"	17.5"-19.5"	19.375"-21.375"	17.75"	N/A 17.75-19.75"	18.5-20.1"	N/A	N/A
MANUAL RECLINE	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes	No
BATTERY SIZE ⁶	NF-22	Group 24: 51.5 lbs. Group 34: 38.5 lbs.	NF-22	NF-22	Group 24	U-1	Group 24	Group 24	U-1	Group 24
BATTERY WEIGHT ¹⁴	38 lbs.	51.5 lbs.	38 lbs.	38 lbs.	52 lbs. each	23.4 lbs.	52 lbs.	52 lbs.	23.4 lbs.	52 lbs.
AVAILABLE ELECTRONICS ¹²	120A Q-Logic 3 EX 75A Q-Logic 2 NE 75A Q-Logic 2 NE+	120A Q-Logic 3 EX	75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX	120A Q-Logic 3 EX	75A Q-Logic 2 NE+ 100A Q-Logic 2 EX	70A PG VR2 70A Q-Logic 2 EX	120A Q-Logic 3 EX 75A Q-Logic 2 NE 75A Q-Logic 2 NE+	75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX	70A Q-Logic 2 NE 70A Q-Logic 2 NE+ 70A Q-Logic 2 EX	90A VR2 6-key 100A Q-Logic 2 EX
BATTERY CHARGER	8A Off-board	8A Off-board	8A Off-board	8A Off-board	8A Off-board	5A Off-board	8A Off-board	8A Off-board	5A Off-board	8A Off-board
MOTOR PACKAGES	2-Pole Mid-Wheel 6®	2-pole, front-wheel drive	4-Pole 5 mph ¹ or 4-Pole 6 mph ²	2-Pole Mid-Wheel 6®	4-Pole Mid-Wheel 6®	2-Pole Glide 4 mph ¹	4-Pole 4.7 mph ¹	4-Pole 6 mph ¹ or 4-Pole 8 mph ²	2-Pole 4 mph	4-Pole 4.5 mph
WEIGHT CAPACITY ⁷	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs. ⁹	450 lbs. ¹⁵	300 lbs.	300 lbs. ⁹	600 lbs.
BASE WEIGHT	138.6 lbs. (6" casters) 148 lbs. (8" casters)	140 lbs.	138 lbs.	138.6 lbs. (6" casters) 148 lbs. (8" casters)	143.5 lbs.	100.2 lbs.	143.5 lbs.	129 lbs.	100.2 lbs.	195 lbs.
MEDICARE CODE COVERAGE - For more information on Medicare code coverage, please consult our Product Planning and Reimbursement Center.										
	Group 3: Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)	Group 3: Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)	Group 3: Standard (K0848/K0849), Single Power (K0856/K0857), Multiple Power (K0861)	Group 3: Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)	Group 3: Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)	Group 2: Standard (K0822/K0823), Single Power (K0835/K0836)	Group 3: Heavy Duty (K0850/K0851), Single Power HD (K0858), Multiple Power HD (K0862)	Group 3: Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)	Group 2: Standard (K0822/K0823), Single Power (K0835)	Group 3: Very Heavy Duty (K0852/K0853), Single Power Very HD (K0860), Multiple Power Very HD (K0863)

1) Standard
2) Optional
3) Requires HD Suspension
4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
5) Without front riggings
6) Battery capacity may change when adding power positioning systems.
7) Weight capacity is dictated by the selection of seating.
8) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Depths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
9) 250 lbs. weight capacity with TRU-Balance® Systems
10) A minimum seating size of 20" x 20" is required with iLevel®
11) Without foot riggings
12) All Quantum units can be ordered with specialty controls
13) iLevel® not available with 8 mph motor package
14) Battery weight may vary +/- 3 lbs. based on manufacturer.
15) Weight capacity is reduced to 400 lbs. with iLevel® power adjustable seat height
16) Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front's suspension system, seat-to-floor heights are approximately 0.25" lower when the unit is loaded.
HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.
The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.

(US) 866-800-2002 • (Canada) 888-570-1113 • quantumrehab.com

