

TRU COMFORT²

CUSHIONS

TRU-Comfort Seating and Positioning products provide adjustable positioning and pressure management and are compatible with both TRU-Balance® Power Positioning Systems and the Synergy® Seat Frame System. A deep contoured, high-density molded foam seat cushion with a viscoelastic total top layer, makes for superb comfort and support.

FEATURES

- Four-way stretch, recovery, fluid-resistant, breathable cover
- Additional 1" of viscoelastic foam in seat cushion for maximum comfort and pressure management
- Pelvic well
- High contour side for femur positioning
- Reduces pressure (ischial tuberosities)



SEAT CUSHION
E2603/E2604 (Skin Protection)



QUANTUM REHAB®

(US) 866-800-2002 • (Canada) 888-570-1113 • www.quantumrehab.com