

THE POWER OF CHOICE

KNOW YOUR OPTIONS

TRU-Balance® 3 Advantages (on the Q6 Edge 2.0)

- ✓ Ability to easily grow or shrink your seat width via overlapping seat pan (4" of built-in adjustment)
- ✓ Easy seat depth and back height adjustment
- ✓ Option for seat width independent of back width
- ✓ Fine adjustments can be made while client is in the seating system
- ✓ Internally or externally rotate cantilever armrests for optimal positioning
- ✓ Seat-to-floor height as low as 17.5"
- ✓ Available in seat sizes as small as 12" wide x 12" deep up to 22" wide x 22" deep
- ✓ Easy center of gravity adjustment
- ✓ Optional attendant control available with enhanced display or SCIM module
- ✓ Tool kit included for quick adjustments



Power Tilt



Power Tilt & Recline



Power Recline



THE POWER OF CHOICE

KNOW YOUR OPTIONS

The iLevel® Advantage (on the Q6 Edge 2.0)

- ✓ Base that can be driven at **3.5 MPH** at 10" of elevation
- ✓ **Only** power adjustable seat that offers 10" of elevation in 16 seconds
- ✓ **Only** system available with WC-19 crash-test approved occupied transit kit option for tilt, recline and iLevel®
- ✓ Assists with transferring to and from wheelchair
- ✓ Client can elevate on the move
- ✓ Option to retrofit iLevel®
- ✓ Available in seat sizes as small as 12"x12"
- ✓ With iLevel, the suspension allows power base to traverse over limited thresholds between 1"-1.5" when elevated

