



DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

**BASE MODEL** **EXPANDABLE CONTROLS**

Description	Price	HCPCS Code
<input type="checkbox"/> Quantum 610 2SP-SS	\$6,171	K0835

**COLORS**

**Colors:**

- Candy Apply Red     Sunburst Orange
- Viper Blue         Pearl Pink
- Onyx Black

**BASE MOUNTED FOOT PLATFORM**

Description	Price
<input type="checkbox"/> Base Mounted Foot Platform Installed	Standard
<input type="checkbox"/> No Base Mounted Foot Platform	No Charge

**DRIVE WHEELS**

Description	Price
<input type="checkbox"/> 10" Flat-Free (pair)	Standard
<input type="checkbox"/> 10" Pneumatic (pair)	No Charge



**BATTERY CHARGER**

Description	Price
<input type="checkbox"/> 5 Amp Off-board Charger	Standard





**LED LIGHT PACKAGE**

Description	Price
<input type="checkbox"/> Through Q-Logic 2 EX w/ Lights	\$995

**NON-EXPANDABLE CONTROLS**



Description	Price
<input type="checkbox"/> PG VR2 4-Key 70 Amp 	Standard
<input type="checkbox"/> PG VR2 6-Key 70 Amp 	No Charge

*Q-Logic 2 EX Joystick or Enhanced Display must be ordered with a Q-Logic 2 EX Controller and Harness For Expandable Controls*

Description	Price
<input type="checkbox"/> 100 Amp Q-Logic 2 EX Controller 	\$695
<input type="checkbox"/> Expandable Controller Harness	\$425
<input type="checkbox"/> Q-Logic 2 EX Joystick 	\$700
<input type="checkbox"/> Q-Logic 2 EX Joystick for Lights 	\$700
<input type="checkbox"/> Stand Alone Joystick w/ External Charger Port 	\$700

**Armor Color Options:**  
 Black Soft Touch (Standard color option)  
**Colors: \$50**  
 Candy Apple Red                       Sunburst Orange  
 Viper Blue                                 Pearl Pink

**PROGRAMMING**

Description	Price
<input type="checkbox"/> Q-Logic Handheld Programmer 	\$545
<input type="checkbox"/> Q-Logic Programming Cable w/ Disk <sup>1</sup>	\$750
<input type="checkbox"/> PG Handheld Programmer 	\$423

1. Works with Q-Logic 2 EX only

**JOYSTICK HANDLES**

Description	Price	HCPCS Code
<input type="checkbox"/> Standard Carrot Handle	Standard	
<input type="checkbox"/> Mushroom Handle	\$85	E2323
<input type="checkbox"/> T-Bar	\$120	E2323
<input type="checkbox"/> Chin Cup	\$85	E2324 <sup>1</sup>
<input type="checkbox"/> Extended I-Handle	\$85	E2323
<input type="checkbox"/> Soft Ball Handle	\$65	E2323
<input type="checkbox"/> Medium Ball Handle	\$65	E2323
<input type="checkbox"/> Bodypoint 3" Goal Post Handle	\$99	E2323
<input type="checkbox"/> Bodypoint 4" Goal Post Handle	\$99	E2323

1. Use HCPCS code E2323 when not used in conjunction with a chin control.

**JOYSTICK MOUNTING BRACKETS**

Description	Price	HCPCS Code
<i>Comes Standard with a Quick Release, Adjustable-Height Joystick Block</i>		
<input type="checkbox"/> Fixed Inline <sup>1</sup>	Standard	
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Swing-away	\$245	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
Offset: <input type="checkbox"/> Inline <input type="checkbox"/> Offset		
<input type="checkbox"/> Retract4 Swing-Away Joystick Arm Kit	\$350	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Multi-axis	\$100	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		

1. Not available with Q-Logic 2 EX Stand Alone Joystick

**TRU-BALANCE® POWER POSITIONING** **BACK OPTIONS (CONT'D)**

Description	Price	HCPCS Code
<input type="checkbox"/> TRU-Balance® Power Tilt (55°) <i>55 degree center of gravity (CG) power tilt seating system. 250lbs Weight Capacity</i>	\$4,995	E1002
<input type="checkbox"/> TRU-Balance® Power Recline (168°) <i>168 degrees of power recline with full drive lock-out set at 20 degrees. 250lbs Weight Capacity</i>	\$5,595	E1004
<input type="checkbox"/> TRU-Balance® Lift (10") <sup>1,2</sup>	\$2,895	E2300

1. 10" of scissor Lift. Must select TB Tilt in conjunction with lift. 1/4 Speed Inhibit, 250 lbs. max. wt., Allowable dimensions 16" x 16" - 20" x 20". Requires Q-Logic 2 EX. Not available with Manual Recline or Power Recline Only.  
2. Not PDAC code verified as a multiple power base.

**SEATING SPECIFICATIONS**

**Seat Width:**

- 14"     15"     16"     17"     18"  
 19"     20"

**Seat Width Packages:**

- 14"W-18"W     16"W-20"W

**Seat Depth:**

- 14"     15"     16"     17"     18"  
 19"     20"

**Seat Depth Package:**

- 14"D-18"D     15"D-20"D

**Back Cane Height:**

- 14"     15"     16"<sup>1</sup>     17"<sup>1</sup>     18"  
 19"     20"     21"     22"     23"  
 24"     25"

1. Selection of a Synergy® Back & Flip Up Armrests disqualified 16" or 17"H back canes.

**BACK OPTIONS**

Description	Price	HCPCS Code
-------------	-------	------------

**Tilt Back:**

- Omit Back    No Charge  
 Synergy® Back with 1.5" Thick Back Cushion    No Charge  
Overall Back     18" H     19" H     20" H     21" H  
Height:     22" H     23" H     24" H     25" H  
Choose Lateral Extensions<sup>1</sup>:  
 6" for 13"-16"H Back    \$25  
 10" for 17"-18"H Back    \$25  
 12" For 19"-22"H Back    \$25  
 16" For 23"-25"H Back    \$25

- TRU-Comfort™ Back 4-way Stretch Material Cover  
Width:     16"W     18" W     20"W    \$675    E2620

- Overall Back     20" H     21" H     22" H     23" H  
Height<sup>2</sup>:     24" H     25" H     27" H     28" H

- Width:     17"W     19"W    \$675    E2620

- Overall Back     22" H     23" H     24" H     25" H  
Height<sup>2</sup>:     27" H     28" H

- Width:     22"W<sup>3</sup>     24"W<sup>3</sup>    \$775    E2621

- Overall Back     24" H     25" H     27" H     28" H  
Height<sup>2</sup>:

1. Lateral extensions are required to mount lateral thoracic supports.  
2. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".  
3. Available with offset (meso) brackets.

Description	Price	HCPCS Code
-------------	-------	------------

**Tilt Back Options:**

- Standard Back Canes    Standard  
 Omit Push Handles    \$125  
 Flip Forward Back Canes<sup>4</sup>    \$150  
 Omit Push Handles    \$125  
 Offset (Meso)<sup>5</sup>    \$195  
 2" (1" Per Side)     4" (2" Per Side)  
 Inset (Endo)<sup>5</sup>    \$195  
 2" (1" Per Side)     4" (2" Per Side)  
 Attendant Manual Recline<sup>6</sup>    \$750    E1226

**Power Recline Back:**

- Synergy® Reclining Back    Standard  
with 1.5" Thick Back Cushion  
 TRU-Comfort™ Back 4-way Stretch Material Cover  
Width:     16"W     17"W<sup>7</sup>     18"W    \$595/\$795    E2620  
 19"W<sup>7</sup>     20"W  
 22"W    \$695/\$795    E2621  
Overall Back  
Height<sup>8</sup>:     20" H     21" H     22" H     23" H  
 24" H     25" H     27" H     28" H  
 After-Market Back Interface<sup>9</sup>    \$200

**Recline Back Options:**

- Offset (Meso)<sup>10</sup>    \$195  
 2" (1" Per Side)     4" (2" Per Side)  
 Inset (Endo)<sup>10</sup>    \$195  
 2" (1" Per Side)     4" (2" Per Side)

**Recline Color Shroud Option:**

- Black (Standard color option)

**Metallic Colors:N/C**

- Candy Apple Red  
 Viper Blue  
 Sunburst Orange  
 Pearl Pink

4. Not available with Attendant Manual Recline, Vent Tray, O2 Holder, cane mounted or cantilever armrests. May interfere with Lateral Supports and Waterfall Armrests.  
5. Dual Post Flip-back Armrests interfere with offset canes. Not available with Manual Recline.  
6. Not available with Flip Back Arms, Cantilever Arms, Shape Back, Tru Comfort back, AFP, rear accessories or 10" Lift. Requires 22"H canes. Pivot point 2.5" above seat pan.  
7. SOF is required with no extended lead time.  
8. Overall Back height measurement is measured from the seat pan to the top of the back. All back heights except 25" and 28" are SOF with no extended lead time.  
9. Overall Seat depth may need to be adjusted to accommodate certain backs. Contact your Quantum Specialist for details. Does not include brackets to mount after market back.  
9. Requires Cantilever Arms. Back width can not be less than 14" or greater than 24".

**SYNERGY® ARMRESTS**

Description	Price
-------------	-------

**Tilt Armrests:**

- Flip Back, Cane Mounted Armrests<sup>1</sup>    No Charge  
 2-Post, Flip-back, Height Adjustable, Removable<sup>2</sup>    No Charge  
 Single Post Height Adjustable Removable<sup>3</sup>    No Charge  
 Quad Release Lever  
 Removable, Single Post, Height Adjustable, Heavy Duty    No Charge  
Height:     8"H - 12"H     10"H - 14"H  
 Pediatric, Removable<sup>4</sup>    No Charge  
 Quad Release Lever  
 Height Adjustable Cantilever Armrests<sup>5</sup>    \$195  
*Angle Adjustable: -30° through +30°*

**Power Recline Armrests:**

- Reclining, Flip Back    Standard  
Height:     8"H - 9"H     10"H - 11"H  
 Height Adjustable Cantilever Armrests<sup>1</sup>    \$195  
*Angle Adjustable: -30° through +30°*

1. Interferes with Swing-Away laterals. Not available w/ manual recline or flip forward back canes.  
2. Height adjustment from 9"-13" in 1" increments. Interferes with swing-away laterals. Not available with a depth less than 16".  
3. Height adjustable from 11" - 17" in 1/2" increments.  
4. Height adjustment from 7.75" - 9.75" in 1/2" increments. Not available with full length armrests.  
5. Height adjustment from 11"-17" depends on cane height. Not available with Manual Recline, Flip Forward Back Canes, Cane Clamp Pushbar, or Shape Back.

**SYNERGY® ARMPADS**

Description	Price
<input type="checkbox"/> Straight, Full Length (2"x14") <input type="checkbox"/> Left <input type="checkbox"/> Right	Standard
<input type="checkbox"/> Straight, Desk Length (2"x10") <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Waterfall, Full Length (2"x14") <sup>1</sup> <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Waterfall, Desk Length (2"x10") <sup>1</sup> <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Falcon Gel, Full Length, (3.5"x14")* <input type="checkbox"/> Left <input type="checkbox"/> Right	\$92.50 ea.
<input type="checkbox"/> Full Length Gel Armpad Mounting Plate (Required)*	\$56.69 ea.
<input type="checkbox"/> Falcon Gel, Desk Length, (3.5"x10")* <input type="checkbox"/> Left <input type="checkbox"/> Right	\$70.00 ea.
<input type="checkbox"/> Desk Length Gel Armpad Mounting Plate (Required)*	\$53.52 ea.

1. Not available w/ Cane Mounted, Flip Back Armrests combined with a Synergy Back.

**POWER POSITIONING ELECTRONICS**

Description	Price	HCPCS Code
<b>Single Actuator Electronics:</b>		
<input type="checkbox"/> Tilt Only Through Toggle <input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount	No Charge	
<input type="checkbox"/> Through VR2 6-Key	No Charge	
<input type="checkbox"/> Through Single Switch (Tilt Only & Q-Logic 2 EX Only)	\$500	
<input type="checkbox"/> Through Q-Logic 2 EX Controller <input type="checkbox"/> AM1 For Single Actuator	\$1,895	E2310
<input type="checkbox"/> AAM For Future Actuator Expansion	\$2,895	
<b>Multiple Actuator Electronics:</b> <i>One TRU Balance Power Positioning Function with AFP</i>		
<input type="checkbox"/> Through Q-Logic 2 EX Controller	\$2,895	E2311
<input type="checkbox"/> Q-Logic 2 EX Through 4-Way Toggle <input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount	\$2,895	
<input type="checkbox"/> Q-Logic 2 EX Through 4-Way Button Toggle <input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount	\$2,895	
<input type="checkbox"/> Q-Logic 2 EX Through Quad Pushbutton <input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount	\$2,895	

**HEADRESTS**

Description	Price	HCPCS Code
<input type="checkbox"/> No Headrest Required		
<b>Stealth:</b>		
<input type="checkbox"/> Small (9.75"W x 3"H) (STL-OTSM250)	\$64.20	E0955
<input type="checkbox"/> Large (11"W x 6"H) (STL-OTLG150)	\$80.25	E0955
<input type="checkbox"/> Pediatric Comfort Plus (6"W) (CPN350)	\$151.94	E0955
<input type="checkbox"/> Small Comfort Plus (10"W) (CP250)	\$151.94	E0955
<input type="checkbox"/> Large Comfort Plus (14"W) (CP150)	\$151.94	E0955
<input type="checkbox"/> Removable Headrest Bracket <sup>2</sup> (HMO475P)	\$160.50	E1028
<b>Aftermarket Bracket:</b>		
<input type="checkbox"/> Stealth's World's Best Removable Headrest Hardware (TWB485)	\$220.50	E1028

**NON-POWERED LEGRESTS**


Description	Price	HCPCS Code
<input type="checkbox"/> <b>60-90° Swing-away Legrests<sup>1</sup></b> (Adjustable in 5° increments)	Standard	
<input type="checkbox"/> Standard Extruded Footplates	No Charge	
Standard Extruded Extensions: <input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
Standard Extruded Footplates: <input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle-Adjustable Footplates <sup>2</sup>	\$85 ea.	
Angle Adjustable Extensions: <input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
Angle-Adjustable Footplates: <input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8" <input type="checkbox"/> 6" x 8" <sup>3</sup> <input type="checkbox"/> 6" x 11" <sup>3</sup> <input type="checkbox"/> 8" x 11" <sup>4</sup>		
<input type="checkbox"/> <b>High Mount, Clamp-On Footrest*</b>		
<input type="checkbox"/> Angle Adjustable		
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<input type="checkbox"/> Multi-Axis		
<input type="checkbox"/> Small, 4" x 7.75"	\$456	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	
<input type="checkbox"/> <b>Heavy Duty Drop-In Legrests</b>	<b>\$420</b>	
<input type="checkbox"/> Angle Adjustable Footplates <sup>2</sup>	\$85 ea.	
Angle Adjustable Extension Length: <input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
Angle Adjustable Footplate Size: <input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8" <input type="checkbox"/> 6" x 8" <sup>3</sup> <input type="checkbox"/> 6" x 11" <sup>3</sup> <input type="checkbox"/> 8" x 11" <sup>4</sup>		
<input type="checkbox"/> <b>Elevating Legrests<sup>5</sup></b>	<b>\$325</b>	E0990
<input type="checkbox"/> Standard Footplates		
Standard Extruded Extension Length: <input type="checkbox"/> 11"-14" <input type="checkbox"/> 14"-17"		
Standard Extruded Footplate Size: <input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle Adjustable Footplates <sup>2</sup>	\$85 ea.	
Angle Adjustable Extension Length: <input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
Angle Adjustable Footplate Size: <input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8" <input type="checkbox"/> 6" x 8" <sup>3</sup> <input type="checkbox"/> 6" x 11" <sup>3</sup> <input type="checkbox"/> 8" x 11" <sup>4</sup>		
<input type="checkbox"/> <b>90° Pediatric Swing-Away Legrests<sup>5</sup></b>	<b>\$162</b>	
<input type="checkbox"/> Angle Adjustable Footplates <sup>2</sup>		
Angle Adjustable Extension Length: <input type="checkbox"/> 7"-10"		
Angle Adjustable Footplate Size: <input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6"		

1. Possible caster wheel interference when adjusted to less than 70 degrees.  
 2. Includes Heel Loops  
 3. Causes interference on seat widths less than 18".  
 4. Causes interference on seat widths less than 20".  
 5. Extensions are measured from the seatpan to footplate.

Items marked with a (\*) should be chosen in the Positioning Components screen in iQuantum

**ADJUSTABLE CENTER MOUNT FOOT PLATFORM**


Description	Price	HCPCS Code
<input type="checkbox"/> Adjustable Center Mount Foot Platform	\$250	
<i>Footplates:</i>		
<input type="checkbox"/> 10"W 11"D Split Footplates		
<input type="checkbox"/> 10"W 8"D Split Footplates		
<input type="checkbox"/> 10"W 6"D Split Footplates		
<i>Extensions:<sup>1</sup></i>		
<input type="checkbox"/> 10"-16" Ext.		
<i>Calf Supports<sup>2:</sup></i>		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small (6"H x 6"W)	\$150	
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)	\$150	
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) <sup>3</sup>	\$200	
<i>Optional Modular Footplate Extension Components/Kits<sup>2,4</sup></i>		
<input type="checkbox"/> 4" Width Extension Kit	\$50	
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/> 4" Side Guard <sup>5</sup>	\$50	
<input type="checkbox"/> 6" Side Guard <sup>6</sup>	\$50	
<i>Additional Optional Accessories<sup>2,4:</sup></i>		
<input type="checkbox"/> Footplate Connector	\$15	
<input type="checkbox"/> Wheel Bumper	\$12.50 ea.	



1. Extension length is measured from top of the seat pan to footplate.  
 2. Choose in Rigging Accessory Section in iQuantum.  
 3. Only available on 20"W or greater.  
 4. Assembly required. May require seating configuration adjustments to avoid front caster interference.  
 5. Only available with 10"W x 6"D split footplates.  
 6. Only available with 10"W x 8"D and 10"W x 11"D split footplates.

**CENTER MOUNT POWER ARTICULATING FOOT PLATFORM**

Description	Price	HCPCS Code
<input type="checkbox"/> Center Mount Power Articulating Foot Platform	\$2995	See note <sup>1</sup>
<i>Footplates:</i>		
<input type="checkbox"/> 10"W 11"D Split Footplates		
<input type="checkbox"/> 10"W 8"D Split Footplates		
<input type="checkbox"/> 10"W 6"D Split Footplates		
<i>Extensions:<sup>2:</sup></i>		
<input type="checkbox"/> 10"-16" Ext.		
<i>Calf Supports<sup>3:</sup></i>		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small (6"H x 6"W)		
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)		
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) <sup>4</sup>	\$200	
<i>Optional Modular Footplate Extension Components/Kits<sup>3,5</sup></i>		
<input type="checkbox"/> 4" Width Extension Kit	\$50	
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/> 4" Side Guard <sup>6</sup>	\$50	
<input type="checkbox"/> 6" Side Guard <sup>7</sup>	\$50	
<i>Additional Optional Accessories<sup>2:</sup></i>		
<input type="checkbox"/> Footplate Connector	\$15	
<input type="checkbox"/> Wheel Bumper	\$12.50 ea.	



1. Please verify the appropriate billing code to be used for this item with all third-party payors.  
 2. Extension length is measured from top of the seat pan to footplate.  
 3. Choose in Rigging Accessory Section in iQuantum.  
 4. Only available on 20"W or greater.  
 5. Assembly required. May require seating configuration adjustments to avoid front caster interference.  
 6. Only available with 10"W x 6"D split footplates.  
 7. Only available with 10"W x 8"D and 10"W x 11"D split footplates.

**FRONT RIGGING ACCESSORIES**

Description	Price	HCPCS Code
<input type="checkbox"/> Wheel Bumper	\$12.50 ea	
<input type="checkbox"/> Heel Loops (pair) <sup>1</sup>	\$50	E0951
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) <sup>2</sup>	\$80	K0108
<input type="checkbox"/> ELR/ALR Adductor Buttons (pair)	\$125	K0108
<input type="checkbox"/> Residual Limb Support	\$300	E1020
<input type="checkbox"/> Left <input type="checkbox"/> Right		

1. Standard with Angle Adjustable or Multi-Axis footplates.  
 2. Only available with Swing-Away legrests.

**SEAT CUSHIONS**

Description	Price	HCPCS Code
<input type="checkbox"/> No Cushion Required		
<i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
<b>Simplicity Cushion</b> General Use		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$98	E2601
<b>Solution@ 1 Cushion</b> Skin Protection:		
<input type="checkbox"/> 14" x 14" - 20" x 20"	\$430	E2603
<b>Solution@ Cushion</b> Skin Protection & Positioning		
<input type="checkbox"/> 10" x 10" - 20" x 20"	\$430	E2607
<b>Spectrum@ Gel Cushion</b> Positioning		
<input type="checkbox"/> 12" x 10" - 20" x 20"	\$530	E2605
<b>Spectrum@ Visco Cushion</b> Skin Protection & Positioning		
<input type="checkbox"/> 10" x 10" - 20" x 20"	\$530	E2607
<b>TRU-Comfort™ 2 4-Way Stretch</b> Skin Protection		
<input type="checkbox"/> 16" x 14" - 20" x 20"	\$400	E2603

**LATERALS**

Description	Price	HCPCS Code
<input type="checkbox"/> <b>Stealth Thoracic Laterals (pair)<sup>1*</sup></b>		
<i>Bracket:</i>		
<input type="checkbox"/> 125° Swing-Away Bracket	\$600	E1028
<input type="checkbox"/> Swing-Away Rail Mounts	\$472.50	E1028
<i>Pad:</i>		
<input type="checkbox"/> Curved Pad, 3"x5", 4-Way Stretch Cover	\$210	E0956
<input type="checkbox"/> Curved Pad, 4"x6", 4-Way Stretch Cover	\$210	E0956
<input type="checkbox"/> <b>Stealth Pelvic/Thigh Laterals<sup>*</sup></b>		
<i>Side:</i>		
<input type="checkbox"/> Left Side <input type="checkbox"/> Right Side		
<i>Bracket:</i>		
<input type="checkbox"/> Removable Bracket (ea.)	\$262.50	E1028
<input type="checkbox"/> Removable Bracket w/ 1" Offset (ea.)	\$400	E1028
<input type="checkbox"/> Flip Down Bracket (ea.)	\$236.25	E1028
<input type="checkbox"/> Flip Down, Full Surface Contact Bracket (ea.)	\$236.25	E1028
<i>Pad:</i>		
<input type="checkbox"/> Short Pad, 4"x4" (ea.)	\$75	E0956
<input type="checkbox"/> Medium Pad, 4"x8" (ea.)	\$75	E0956
<input type="checkbox"/> Long Pad, 4"x12" (ea.)	\$75	E0956
<i>Adapter:</i>		
<input type="checkbox"/> 2" Depth/Height Adapter Plate (ea.) <sup>2</sup>	\$26.25	K0108
<input type="checkbox"/> 2 1/4" Depth/Height Adapter Plate (ea.) <sup>3</sup>	\$26.25	K0108

1. Swing-away laterals interfere with flip-up armrests.  
 2. Height adapter plate is needed w/ Tru Comfort seat cushion. Not available on the flip down bracket.  
 3. Only available on the flip down bracket.



LAPBELTS	
Description	Price
<input type="checkbox"/> 38" Lap Belt	Standard
<input type="checkbox"/> 50" Lap Belt	Standard
<input type="checkbox"/> 60" Lap Belt	Standard
<input type="checkbox"/> 60" Padded Lap Belt	\$115
<input type="checkbox"/> 70" Lap Belt	Standard
<input type="checkbox"/> 80" Lap Belt	Standard
<input type="checkbox"/> 90" Lap Belt	Standard

I-FIT PELVIC BELTS	
Description	Price
<input type="checkbox"/> i-Fit 60" Single Pull, Push Button, Medium (7.5") Padded 1.5" Belt	\$86
<input type="checkbox"/> i-Fit 60" Single Pull, Push Button, Large (9.5") Padded 2" Belt	\$86
<input type="checkbox"/> i-Fit 60" Single Pull, Push Button, Extra Large (11.5") Padded 2" Belt	\$86

ACCESSORIES		
Description	Price	HCPCS Code
<input type="checkbox"/> Oxygen Holder <sup>1</sup>	\$150	E2208
<input type="checkbox"/> Recline Accessory Mounting Bracket <sup>2</sup>	\$50	
<input type="checkbox"/> Weather Cover, Medium, Black	\$168	
<input type="checkbox"/> Cup Holder	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Tablet Holder	\$315	
<input type="checkbox"/> Deluxe Cell Phone/Media Holder	\$140	

1. Not available with a Power Elevating Seat or a Lift/Tilt. Oxygen holder on Tilt with Tru Comfort back requires SOF. N/A on Power Recline.  
 2. 1" adaptor tube allows cane mounted accessories to mount to power recline back or siderail.

BATTERY OPTIONS		
Description	Price	HCPCS Code
<input type="checkbox"/> U1 32 Amp, AGM <i>Installation:</i>	\$133 ea.	E2365
<input type="checkbox"/> Batteries Installed <input type="checkbox"/> Do Not Install		

ADDITIONAL ACCESSORIES			
Description	Price	HCPCS Code	
<input type="checkbox"/> <b>Stealth Foot Boxes (pair)</b>	\$498.75	K0108	
<input type="checkbox"/> Small			
<input type="checkbox"/> Medium			
<input type="checkbox"/> Large			
<input type="checkbox"/> <b>Stealth Padded Shoe Holders (pair)</b>	\$266	K0108	
<input type="checkbox"/> Small			
<input type="checkbox"/> Medium			
<input type="checkbox"/> X-Large			
<input type="checkbox"/> <b>Stealth Elbow Support</b>	\$183.75 ea.	E1028	
<i>Side:</i>			
<input type="checkbox"/> Bracket - Left			
<input type="checkbox"/> Bracket - Right			
<i>Pad:</i>			
<input type="checkbox"/> 3" x 3" (ARMS-ES-3x3)	\$52.50 ea.	K0108	
<input type="checkbox"/> 4" x 4" (ARMS-ES-4x4)	\$52.50 ea.	K0108	
<input type="checkbox"/> <b>Stealth i-Fit Dynamic Chest Straps</b>	\$142	E0960	
<input type="checkbox"/> Upper Classic, Medium (IFPS1103)			
<input type="checkbox"/> Upper Classic, Large (IFPS1104)			
<input type="checkbox"/> Upper Contour, Medium (IFPS2103)			
<input type="checkbox"/> Upper Contour, Large (IFPS2104)			
<input type="checkbox"/> <b>Stealth i-Fit Chest Strap Belt Anchor</b>	\$90		
<input type="checkbox"/> Medium (IFPS4113)			
<input type="checkbox"/> Large (IFPS4114)			

ADDITIONAL HEADREST ACCESSORIES			
Description	Price	HCPCS Code	
<input type="checkbox"/> <b>Stealth Facial Laterals<sup>1</sup></b>			
<input type="checkbox"/> Right Side, 4.5" x 2" Pad (CPS9-R-HW & 900-Q)	\$347.75	E1028	
<input type="checkbox"/> Left Side, 4.5" x 2" Pad (CPS9-L-HW & 900-Q)	\$347.75	E1028	
<input type="checkbox"/> Right and Left 4.5" x 2" Pad (CPS9-B-HW & 900-Q)	\$453.75	E1028	
<input type="checkbox"/> Egg Switch Mount (ESM)	\$52.50 ea.		
<input type="checkbox"/> Facial Pad, 6" x 2" (Q920-Q)	\$53.50 ea.	K0108	
<input type="checkbox"/> Figure Eight Bracket with 13" Rod and Egg Switch Mount (F8678-13-Q)	\$220	K0108	
<input type="checkbox"/> Swing-away Joystick Mount, 18" Rod (SUS18-Q)	\$250	K0108	

1. Only available with Stealth Comfort Series Headrests.