

# Quantum® 1402

All Inclusive

Order Form / Price List  
Revised: December 2005



## Pride Mobility Products Corporation®

182 Susquehanna Ave., Exeter, PA 18643  
Sales: Phone: (866) 800-2002 • Fax: (866) 707-3422

This product is available for purchase online at [www.mypride.com](http://www.mypride.com)

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:			P.O. #:
Address:	Address:			Discount:
Phone Contact:	Marked For:			Terms:
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.				Freight:
				Taken By:

All pricing is MSRP.

### BASE

- Quantum 1402 K0014 \$8695  
*Top speed of 4.5 MPH,  
400 lb. weight capacity*
  - High-Speed Motors NC  
*Programmable speed up to 6 MPH,  
325 lb. weight capacity.*
- Comes standard with the Synergy Seat, Std. seating configuration of 18"W x 16"D x 18"H includes 2" foam cushion and sling back*

### COLOR

- Champagne NC
- Candy Apple Red NC
- Viper Blue NC
- Onyx Black NC
- Sunburst Orange NC

### CONTROLLER

- PG Electronics**
- 3-Key Remote Plus 70 amp <sup>13,19</sup> Std.  
*(Select Box for Profile)*  
 1 Profile  Multiple Profile
  - 7-Key Access. Remote Plus 70 amp <sup>19,13,32</sup> **E2320** \$283
  - Remote Plus 70 amp with Speed Dial <sup>15,13</sup> **E2320** \$197
  - Remote Plus 70 amp with Toggles <sup>15,13</sup> **E2320** \$230
  - Remote Plus 70 amp with Switch Jacks <sup>15,13</sup> **E2320** \$185
  - Microdrive, 70 amp <sup>5,15</sup> **E2320** \$545

### SPECIALTY CONTROLS

- See Specialty Controls Order Form <sup>4</sup>

**RECOMMENDED INFORMATION**

Height: \_\_\_\_\_

Weight: \_\_\_\_\_

Leg Length: \_\_\_\_\_

Lower Leg Length: \_\_\_\_\_

### JOYSTICK HANDLE OPTIONS

- PG Controls**
- Carrot Handle Std.
  - Mushroom Handle **E2323** \$71
  - T-Bar **E2323** \$107
  - Chin Cup **E2323** \$71
  - Extended I-Handle **E2323** \$71
  - Soft Ball Handle **E2323** \$50
  - Medium Ball Handle **E2323** \$50

### CONTROLLER MOUNT

- Fixed, In-Line Joystick Mount Std.
    - Left  Right
  - Quick Disconnect, Joystick Mount <sup>34</sup> NC
    - Left  Right
  - Swing-away Joystick Mount **E1028** \$196
    - Left  Right
  - Multi-axis Joystick Mount **K0108** \$100
    - Left  Right
- (Fixed, does not swing-away)*

### PROGRAMMER

- Remote Plus Programmer \$423

### DRIVE WHEELS

- 14" Pneumatic (pair) Std.
- Flat-Free Inserts (pair) **K0093** \$213

### BATTERIES (2 required)

- Group 24, 70 amp, Gel **E2363** \$315
- Group 34, 60 amp, Gel **K0108** \$300  
Required w/ Super Low Tilt
- NF-22, 50-amp Batteries **E2361**
  - Maintenance-free AGM \$225
  - Maintenance-free Gel Cell \$265Required with power tilt / recline elevating system
- Ship unit with batteries installed \$200
- Ship unit with batteries, but do not install NC

### BATTERY CHARGER

- 8 amp Off-board Charger Std.
- 5 amp On & 8 amp Off-board Charger \$220  
8 amp Off-board charger required with super low seating.

### ACCESSORIES

- Rear Basket <sup>2,24</sup> \$103
  - Cane/Crutch Holder <sup>2,24</sup> \$98
  - Oxygen Holder <sup>2,24</sup> **K0104** \$103
  - Oxygen Holder <sup>2,24,29</sup> **K0104** \$150  
*(TRU-Balance Tilt & Synergy Seat only)*
  - Stroller handle \$85  
*(Synergy Manual Tilt only)*
  - Accessory Mounting Bracket \$50  
*(1" adaptor tube for use w/ TB Power Recline)*
  - Walker Holder <sup>2,24</sup> \$92
  - Flag Holder \$17
  - Dust Cover \$192
    - Gray  Burgundy
  - Saddle bag (contoured seating only) \$28
  - Cup holder <sup>33</sup> \$17
    - Left  Right
- N/A on desk length arm with S/A mount or toggle box.
- Lap Belt **E0978** \$38
  - Pediatric 38" <sup>27</sup>  50"  60"
  - 70"  80"  90"

### CUSHIONS<sup>11</sup>

- Simplicity Cushion
  - 16 x 16 - 20 x 18 **E2601** \$98
- Solution Cushion
  - 10 x 10 - 20 x 18 **E2607** \$380
  - 20 x 20 - 20 x 22 **E2607** \$430
  - 22 x 16 - 24 x 24 **E2608** \$430
- Spectrum Cushion
  - Gel
    - 10 x 10 - 20 x 18 **E2605** \$480
    - 20 x 20 - 20 x 24 **E2605** \$530
    - 22 x 16 - 24 x 24 **E2606** \$530
  - Visco
    - 10 x 10 - 20 x 18 **E2607** \$480
    - 20 x 20 - 20 x 24 **E2607** \$530
    - 22 x 16 - 24 x 24 **E2608** \$530
- Structure Cushion
  - Gel  Visco
    - 10 x 10 - 20 x 18 **E2607** \$560
    - 20 x 20 - 20 x 22 **E2607** \$610
    - 22 x 16 - 24 x 24 **E2608** \$610

### POSITIONING COMPONENTS

- See Positioning Components Order Form

# Quantum® 1402

Order Form / Price List  
Revised: December 2005

## FRONT RIGGINGS

- Swing-away footrests Std.  
(\$162 w/ contoured or Synergy seating)
- Elevating Leg Rests E0990 \$325  
Adj. Calf Pads w/ ELRS (1 hole std.)  
 1 hole  3 holes\*  6 holes\*
- \*Potential for interference w/ shrouds
- HD Drop-Ins, 70<sup>07,10</sup> K0108 \$420
- HD Drop-Ins, 80<sup>07,10</sup> K0108 \$420  
(HD Drop-ins N/A w/ contoured or Specialty seating)
- Swing-away, Pediatric K0040 \$162  
 4"x6"  5"x6" Angle-Adj. (pair)
- Drop-ins, Pediatric K0040 \$162  
 4"x6"  5"x6" Angle-Adj. (pair)
- HD Foot Platform<sup>25</sup> NC
- No Front Rigging NC

Power legrests are only available on TRU-Balance and Synergy Seating Systems. Only the combined motion can be selected on the Synergy Seat.

- Power-elevating, high-pivot<sup>21,25</sup> E1010 \$1695  
Legrests, (calf pads included)
- Independent  Combined<sup>21</sup>  Synchronized<sup>36</sup>
- Power-elevating, Articulating<sup>21,25,35</sup> E1010 \$2395  
Articulating Range:  4" L Std.  6" L NC  
Set at:  70°  80°
- Independent  Combined<sup>21</sup>  Synchronized<sup>36</sup>
- Calf Pads Left & Right (6 x 5.35) Std.
- Calf Pads Left & Right (8.00 x 5.86) NC
- Footplates (pair)**
- Reinforced Footplate Std.  
 7.25" x 6"  8.5" x 6"
- 5"x8" Angle-Adj. Footplate<sup>37</sup> K0040 \$150
- 6"x8" Angle-Adj. Footplate<sup>30</sup> K0040 \$150
- 6"x11" Angle-Adj. Footplate<sup>30</sup> K0040 \$150
- 8"x11" Angle-Adj. Footplate<sup>31</sup> K0040 \$150
- Wheel Bumpers (Ang-adj. only) \$25

### Extensions:<sup>12</sup>

- 11"-14" Ext.<sup>22</sup> NC
- 14"-17" Ext. NC
- 17"-20" Ext.<sup>28</sup> NC
- High mount Clamp-on; Ang. Adj. K0040 \$326  
 Small  Large
- High mount Clamp-on; Multi-axis K0040 \$412  
 Small  Large

- High mount footplatform K0108 \$450  
 9.5" deep  11.5" deep
- Lower Extensions**
- 10"-14" L
- 14"-17" L
- 17"-20" L

### Rigging Options

- Heel Loops (pair)<sup>14</sup> E0951 \$20
- Gel Neoprene Footrest K0108  
Sleeve (pair)<sup>16</sup> \$65
- ELR/ ALR Adductor Buttons (pair) K0108 \$100
- POSAlinc Knee Adductor (pair) K0108  
for S/A legrests \$464
- Residual limb support left E1020 \$300
- Residual limb support right E1020 \$300

\*SEAT DIMENSION CODING FOR SPECIALTY & SYNERGY SEATING 20" AND LARGER IS K0108

## Synergy Seating Systems<sup>23</sup>

- Synergy Seating™ System Std.  
(Price with Base)
- Synergy Manual Recline Pkg.<sup>1,3,6,17</sup> \$700  
Maximum weight 300 lbs. E1226
- Synergy Manual Tilt Pkg.<sup>6</sup> \$550  
(Max Dimensions 18"W x 18"D)  
Maximum weight 160 lbs. K0108
- Synergy Manual Systems come standard with mercury drive lock-out and are in addition to the MSRP of the Synergy Seat.**
- Synergy Quick Release<sup>8,20</sup> \$200  
(Max Dimensions 20"W x 18"D)

**Synergy Quick Release is in addition to the MSRP of the Synergy Seat.**

## WIDTH

### Set Seat Width At:

- |                                   |                                   |
|-----------------------------------|-----------------------------------|
| <input type="checkbox"/> 10" \$20 | <input type="checkbox"/> 11" \$20 |
| <input type="checkbox"/> 12" \$20 | <input type="checkbox"/> 13" \$20 |
| <input type="checkbox"/> 14" \$20 | <input type="checkbox"/> 15" \$20 |
| <input type="checkbox"/> 16" \$20 | <input type="checkbox"/> 17" \$20 |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" \$20 |
| <input type="checkbox"/> 20" \$20 | <input type="checkbox"/> 21" \$20 |
| <input type="checkbox"/> 22" \$20 |                                   |

### Seat Width Package

- |                                    |                                    |
|------------------------------------|------------------------------------|
| <input type="checkbox"/> 10"W-14"W | <input type="checkbox"/> 14"W-18"W |
| <input type="checkbox"/> 16"W-20"W | <input type="checkbox"/> 20"W-24"W |

### Seat Width Package Manual Tilt

- |                                    |                                    |
|------------------------------------|------------------------------------|
| <input type="checkbox"/> 10"W-12"W | <input type="checkbox"/> 12"W-14"W |
| <input type="checkbox"/> 14"W-16"W | <input type="checkbox"/> 16"W-18"W |

## DEPTH

### Set Seat Depth At:

- |                                   |                                   |
|-----------------------------------|-----------------------------------|
| <input type="checkbox"/> 10" \$20 | <input type="checkbox"/> 11" \$20 |
| <input type="checkbox"/> 12" \$20 | <input type="checkbox"/> 13" \$20 |
| <input type="checkbox"/> 14" \$20 | <input type="checkbox"/> 15" \$20 |
| <input type="checkbox"/> 16" Std. | <input type="checkbox"/> 17" \$20 |
| <input type="checkbox"/> 18" \$20 | <input type="checkbox"/> 19" \$20 |
| <input type="checkbox"/> 20" \$20 |                                   |

### Seat Depth Package

- |                                    |                                    |
|------------------------------------|------------------------------------|
| <input type="checkbox"/> 10"D-14"D | <input type="checkbox"/> 14"D-18"D |
| <input type="checkbox"/> 16"D-20"D |                                    |

### Seat Depth Package Manual Tilt

- |                                    |                                    |
|------------------------------------|------------------------------------|
| <input type="checkbox"/> 10"D-14"D | <input type="checkbox"/> 14"D-16"D |
| <input type="checkbox"/> 16"D-18"D |                                    |

## BACK CANES

Synergy Manual Recline comes standard with 22" non-height adjustable back canes.

### Back Cane Height

- |                                   |                                 |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC   | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" NC   | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC   | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC   | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC   | <input type="checkbox"/> 25" NC |

## Footnotes for page 1 & 2:

1. Cannot use flip-back, cane mounted arms.
2. Includes mounting bracket and hardware.
3. Includes removable, adjustable, quick-height armrests.
4. May require special order, call for lead-time, availability, and pricing.
5. Requires multi-axis joystick mount; No charge for mount.
6. Solid back not included.
7. Standard w/ heel loops and angle-adjustable footplates, check size below.
8. Cannot be used with Synergy manual tilt or recline, units with power elevating seat or single seat post.
9. For use when controlling power elevating seat or TRU-Balance Power Tilt.
10. May cause interference w/ front anti-tips w/ certain seat widths and lower leg lengths. Call Quantum Customer Service.
11. See Synergy cushion order form for detailed pricing, sizes, and coding.
12. Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured seating or rehab seating w/ cushions are selected. The overall length may vary approximately 3" depending on choice of seat and /or cushion.
13. Needed for controlling power positioning seating systems.
14. Std. when angle-adjustable footplates are selected.
15. Available w/ five individually programmable drive profiles.
16. Available w/ swing-away footrests, 70°-80° heavy-duty drop-ins.
17. N/A with Shape back.
18. Reflects price when purchased w/ this unit.
19. Std. w/ one programmable drive profile.
20. St. rear mount accessories are N/A.
21. Requires the addition of an ALM.
22. Special order on ALR's.
23. Manual systems do not include a headrest.
24. N/A on TRU-Balance Power Recline System or fully reclining seating.
25. Power legrests are only available on TRU-Balance and Synergy Seating Systems. Only the combined & independent motion can be selected on the Synergy Seat.
26. No charge when substituting for S/A legrests.
27. N/A w/ limited recline seating.
28. N/A with Super-low systems.
29. N/A w/ flip-back arms.
30. Only avail. on 16" W or wider seats.
31. Only avail. on 20" W or wider seats.
32. Required when selecting full light package.
33. N/A on desk length arm with S/A mount or toggle box.
34. N/A w/ waterfall armpads; only available on Quantum Seating Systems.
35. Std. w/ wheel bumpers.
36. Only avail. with TB Recline.
37. Only avail. on 14" W or wider seats.

# Quantum® 1402

Order Form / Price List  
Revised: December 2005

## Synergy Seating Systems Cont.

### BACK OPTIONS

- See Synergy Shape Back Order Form
- Solid Curved Back up to 20" W \$281
- Solid Curved Back 21"-24" W \$338
- Back Height:**
- 18"H     19"H     20"H     21"H
- 22"H     23"H     24"H     25"H
- Sling Back<sup>1</sup> \$98  
*Back height will be the same as back can height. Includes gap for cushion.*
- Push Bar selection must match back width; center tube standard when selecting flip-up arms.**

### ARMRESTS

- Flip-back, cane mounted<sup>2,4,8</sup> Std.
- 2-Post, flip-back, height<sup>3,4,8</sup> \$150ea  
adjustable, removable armrests
- Removable, adjustable<sup>5</sup> E0973 \$75ea  
quick height  
*(3" of height adjustment in 1/2" increments)*
- Removable, height adjustable E0973 \$75ea  
heavy duty armrests
- 8"H-12"H                       10"H-14"H
- Pediatric removable \$75ea  
adjustable, quick height.  
*(7 3/4" - 10 3/4" height adjustment in 1/2" increments from seat pan to top of armrest)*

### ARMPADS

- Full length, Straight, 2" X 14" Std.
- Left                               Right
- Desk Length, Straight, 2" X 10" Std.
- Left                               Right
- Full Length, Waterfall, 2" X 14"<sup>6</sup> NC
- Left                               Right
- Desk Length, Waterfall, 2" X 10"<sup>6,7</sup> NC
- Left                               Right

### HEAD REST<sup>9</sup>

- Stealth** <sup>23,24</sup> E0955  
Large (11"wide x 6" high) \$210
  - Stealth** <sup>23,24</sup> E0955  
Small (9 3/4" wide x 3" high) \$210
  - Stealth, swing-away (Flip-Down)**  
Head Rest Bracket E1028 \$150
- Interfaces w/ any headrest that uses 1/2" square down-tube.
- Medical Applications** <sup>23,24</sup> E0955  
Large (11"wide x 6" high) \$275
  - Medical Applications** <sup>23,24</sup> E0955  
Small (9 3/4" wide x 3" high) \$275
  - Fixed, Offset Mounting Bracket NC

### TRU-BALANCE TILT E1002

- TRU-Balance Power Tilt \$4995
- TRU-Balance Super Low Power Tilt<sup>22</sup> \$4995  
*55 degree center of gravity (CG) power tilt seating system. 300 lbs. weight capacity. Requires Group 34 batteries and off-board charger*
- HD package-400 lbs.weight capacity NC  
will raise seat-to-floor 1/2"
- TRU-Balance \$1895  
Power-Elevating Seating Option<sup>11</sup>  
**7" of Lift & 1/4 Speed Inhibit- 300 lbs. max. wt.**  
*Requires NF-22 batteries*

### SYSTEM SPECIFICATIONS

#### Tower Selection<sup>12</sup>

- Adjustable-angle Std.  
*Will raise minimum STF by 1".*
- Fixed NC

#### Set Seat Width At:

- |                                   |                                 |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC   | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" NC   | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC   | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC   |                                 |

#### Seat Width Package

- 14"W-18"W                       16"W-20"W
- 20"W-22"W

#### Set Seat Depth At:

- |                                   |                                 |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC   | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" Std. | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" NC   | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC   |                                 |

#### Seat Depth Package

- 14"D-18"D                       16"D-20"D

*Push Bar selection must match back width; center tube standard when selecting flip-up arms.*

#### Set Back Cane Height At:

- |                                   |                                 |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC   | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" Std. | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" NC   | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC   | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC   | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC   | <input type="checkbox"/> 25" NC |

#### Seat Overall Back Height<sup>13</sup>

- |                                   |                                 |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC   | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" NC.  | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC   | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC   | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC   | <input type="checkbox"/> 25" NC |

#### Back Type

- See Synergy Shape Back Order Form
- Solid Curved Back Std.  
Full width, fully T-nutted, 1" foam with naugahyde sides and back and a rubatex face.
- Attendant manual recline<sup>14,15</sup> \$750  
180 degrees of recline with 2 1/2" raised pivot point. Controlled through dual gas struts. 300 lb. weight capacity.
- Omit Back -(\$100)

### ARMRESTS

- Flip-back, cane mounted<sup>8,16,25</sup> Std.  
*Mounted on back canes with full width spreader bar.*
- 2-Post, flip-back, hgt. adj, removable<sup>8,25</sup> NC
- Removable, adj. quick height<sup>5</sup> NC
- Removable, single-post, height adjustable, heavy duty NC
- 8"-12" Adjustment
- 10"-14" Adjustment
- Pediatric, removable, adjustable, quick height NC

#### Armrest Pads

- Full Length, straight, 2"x14" NC
- Left                               Right
- Desk Length, straight, 2"x10" NC
- Left                               Right
- Full Length, Waterfall, 2"x14"<sup>17</sup> NC
- Left                               Right
- Desk Length, Waterfall, 2"x10"<sup>17,18</sup> NC
- Left                               Right

### POWER POSITIONING ELECTRONICS

#### PG Electronics

- Through Toggle Box<sup>18,21</sup> NC
- Left                               Right (mounts to side of armrest)
- Through Single Switch Tilt Module<sup>21</sup> E2310 \$500

#### Single Motor Controller<sup>21</sup>:

- Through Remote Plus Controller<sup>21</sup> E2310\* \$1175
- Through ALM w/ 2 Jacks<sup>21,25</sup> E2310\* \$1175

#### Multiple Motor Controller<sup>17</sup>:

- Through Remote Plus Controller<sup>17</sup> E2311\* \$2375
- Through Toggle<sup>17</sup> (Select one) E2311\* \$2375
- 4-way Toggle
- 4-way Button Toggle
- Quad Push-Button Control

*Toggle box mounts to side of armrest.*

- Left                               Right

- Attendant Control Pod NC  
*(avail. only w/ power elevating seat)*

*\*E2310 coding may vary due to configuration*

### HEAD REST

- Stealth** <sup>23,24</sup> E0955  
Large (11"wide x 6" high) \$210
- Stealth** <sup>23,24</sup> E0955  
Small (9 3/4" wide x 3" high) \$210
- Stealth, swing-away (Flip-Down)**  
Head Rest Bracket E1028 \$150

Interfaces w/ any headrest that uses 1/2" square down-tube.

#### Medical Applications <sup>23,24</sup>

- Medical Applications** <sup>23,24</sup> E0955 \$275  
Large (11"wide x 6" high)
- Medical Applications** <sup>23,24</sup> E0955 \$275  
Small (9 3/4" wide x 3" high)
- Omit Head Rest - To be supplied by provider
- Fixed, Offset Mounting Bracket NC
- Bracket for after market headrest \$180

# Quantum® 1402

Order Form / Price List  
Revised: December 2005

## TRU-BALANCE RECLINE

TRU-Balance Power Recline--300lbs. \$5595

**Full Drive lock-out at 20 degrees** E1004

TRU-Balance \$1895

Power-Elevating Seating Option<sup>11,19</sup>

**7" of Lift & 1/4 Speed Inhibit- 300 lbs. max. wt.**

## SYSTEM SPECIFICATIONS

### Set Width At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	

### Seat Width Package

<input type="checkbox"/> 14"W-18"W	<input type="checkbox"/> 16"W-20"W
<input type="checkbox"/> 20"W-22"W	

### Set Depth At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	

### Seat Depth Package

<input type="checkbox"/> 14"D-18"D	<input type="checkbox"/> 16"D-20"D
------------------------------------	------------------------------------

### Set Back Height At:<sup>13</sup>

<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" Std.	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC
<input type="checkbox"/> 26" NC	<input type="checkbox"/> 27" NC
<input type="checkbox"/> 28" NC	

### Back Type

TRU-Balance Reclining Back Std.

## ARMRESTS

Reclining, Flip-back<sup>8</sup> Std.

8"-9" H  10"-11" H

### Armrest Pads

Full Length, straight, 2"x14" NC

Left  Right

Desk Length, straight, 2"x10" NC

Left  Right

Full Length, Waterfall, 2"x14"<sup>17</sup> NC

Left  Right

Desk Length, Waterfall, 2"x10"<sup>17</sup> NC

Left  Right

## POWER POSITIONING ELECTRONICS

All options require motor control module.

PG Electronics- Multiple Motor controller<sup>20</sup>

### Single Function:

Operation through joystick E2310 \$1175

Operation through toggle E2310 NC

(Select one)

4-way Toggle

4-way Button Toggle

Quad Push-Button Control

Toggle box mounts to side of armrest.

Left  Right

## Power Positioning Electronics Continued

PG Electronics- Multiple Motor controller

### Multiple Function:

Operation through joystick E2311 \$2375

Operation through toggle E2311 \$2375

(Select one)

4-way Toggle

4-way Button Toggle

Quad Push-Button Control

Toggle box mounts to side of armrest.

Left  Right

## HEAD REST

Reclining Head Rest (6.5" x 7")<sup>23,24</sup> E0955 \$180

Omit Head Rest - To be supplied by provider

Bracket for after market headrest \$180

## TRU-BALANCE TILT & RECLINE

TRU-Balance Power E1007  
Tilt & Recline \$9595

TRU-Balance Super Low Power  
Tilt & Recline<sup>10,22</sup> E1007 \$9595

**Full Drive lock-out at 20 degrees**

**55 degree center of gravity (CG) power tilt seating system. 300 lbs. weight capacity**

TRU-Balance \$1895

Power-Elevating Seating Option<sup>10,11,19</sup>

**7" of Lift & 1/4 Speed Inhibit- 300 lbs. max. wt.**

## SYSTEM SPECIFICATIONS

### Tower Selection<sup>12</sup>

Adjustable-angle Std.

Will raise minimum STF by 1".

Fixed NC

### Set Width At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	

### Seat Width Package

<input type="checkbox"/> 14"W-18"W	<input type="checkbox"/> 16"W-20"W
<input type="checkbox"/> 20"W-22"W	

### Set Depth At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	

### Seat Depth Package

<input type="checkbox"/> 14"D-18"D	<input type="checkbox"/> 16"D-20"D
------------------------------------	------------------------------------

### Set Back Height At<sup>13</sup>:

<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" Std.	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC
<input type="checkbox"/> 26" NC	<input type="checkbox"/> 27" NC
<input type="checkbox"/> 28" NC	

### Back Type

TRU-Balance Reclining Back Std.

## ARMRESTS

Reclining, Flip-back<sup>8</sup> Std.

8"-9" H  10"-11" H

### Armrest Pads

Full Length, straight, 2"x14" NC

Left  Right

Desk Length, straight, 2"x10" NC

Left  Right

Full Length, Waterfall, 2"x14"<sup>17</sup> NC

Left  Right

Desk Length, Waterfall, 2"x10" NC

Left  Right

## POWER POSITIONING ELECTRONICS

PG Electronics--multiple motor controller

Operation through joystick E2311 \$2375

Operation through toggle E2311 \$2375

(Select one)

4-way Toggle

4-way Button Toggle

Quad Push-Button Control

Toggle box mounts to side of armrest.

Left  Right

## HEAD REST

Reclining Head Rest (6.5" x 7")<sup>23,24</sup> E0955 \$180

Omit Head Rest - To be supplied by provider

Bracket for after market headrest \$180

### Footnotes for page 3 & 4:

- Not available on on Synergy manual recline, tilt & power positioning..
- N/A on Synergy Manual Recline.
- May cause interference with certain dimension configurations. Call Quantum Customer Service.
- Call for configuration avail. w/ Shape back.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17" from seat pan to top of armrest pad.
- Only available on removable, adjustable, quick height arms or 2-post arms.
- Swing-away joystick mount is N/A.
- Stealth S/A laterals interfere w/ flip back arms.
- N/A with sling back.
- Requires Group 34 batteries.
- Seat-to-floor height will increase.
- Tower selection is N/A when ordering super-low systems.
- Measured from seat pan to top of back pan with a gap to accommodate the cushion.
- N/A w/ flip-back, cane mounted arms.
- Includes price for removable, adjustable, quick-height arms and 2-post arms.
- N/A w/ Attendant manual recline.
- Required w/ power elevating seat.
- N/A with waterfall armpads.
- Requires NF-22 batteries.
- Recline only--N/A w/ power legrest or power elevating seat.
- N/A w/ power elevating seat.
- The selection of a foot platform on a super low system requires a special order.
- Comes with swing-away / removable hardware; HCPCS code is E1028.
- These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
- Consult specialty controls order form for a wide variety of switch type and sizes-- 2 required.

**Quantum® 1402**

Order Form / Price List  
Revised: December 2005

\*SEAT DIMENSION CODING FOR CONTOURED SEATING 20" AND LARGER IS K0108

**CONTOURED SEATING ONLY**

**Medium-Back (-\$364)**

**Medium-Back<sup>1,2,4</sup> w/ Solid Seat Pan (-\$364)**

- Brown Vinyl
- Beige Vinyl
- Gray Vinyl
- Black Vinyl
- Gray Fabric
- Brown Fabric

**Medium-Back Seating Dimensions**

(Dimensions 18"W x 18"D)

<b>Width:</b>		<b>Depth:</b>	
<input type="checkbox"/> 16" NC		<input type="checkbox"/> 16" NC	
<input type="checkbox"/> 18" Std.		<input type="checkbox"/> 18" Std.	

**UPGRADED CONTOURED SEATING**

**High-Back with Headrest (-\$364)**

**High-Back with Headrest<sup>1,2,4</sup> w/ Solid Seat Pan (-\$364)**

- Black Vinyl
- Gray Vinyl
- Beige Vinyl
- Charcoal Fabric
- Brown Fabric

**High-Back Seating Dimensions**

(Dimensions 18"W x 18"D)

<b>Width:</b>		<b>Depth:</b>	
<input type="checkbox"/> 16" NC		<input type="checkbox"/> 16" NC	
<input type="checkbox"/> 18" Std.		<input type="checkbox"/> 18" Std.	
<input type="checkbox"/> 20"* \$325		<input type="checkbox"/> 20"* \$454	
<input type="checkbox"/> 22"* \$325		<input type="checkbox"/> 22"* \$454	

**ARMRESTS**

**For Medium and High-Back Seats**

- Flip-up, removable Std.
- Flip-up, height-adjustable removable<sup>3</sup> E0973 \$170

**Armpads**

- Desk Length, 10" NC
  - Left
  - Right
- Full Length, 14" NC
  - Left
  - Right
- Black Vinyl Armpads Std.

**POWER ELEVATING SEAT<sup>5,6,7,8</sup>**

**400 lbs. weight capacity**

- Power Elevating Seat \$975
  - Through Toggle
- Power Elevating Seat Through Remote Controller<sup>7</sup> \$1375
- Lap Belt Std.

**MANUAL RECLINER<sup>8,9,10</sup>**

**Manual Recline High Back (-\$54)**

**Manual Recline High Back<sup>1,2,4</sup> w/ Solid Seat Pan (-\$54)**

**(180 degrees of Recline)**

(Dimensions 18"W x 18"D)

<input type="checkbox"/> Gray Vinyl	<input type="checkbox"/> Beige Vinyl
<b>Width:</b>	<b>Depth:</b>
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 16" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 18" Std.
<input type="checkbox"/> 20"* \$325	<input type="checkbox"/> 20"* \$454

When ordering a lap belt with this seat, an additional 10" is added to accommodate

**Footnotes for page 5:**

1. Includes black vinyl seat pan, back available in all listed colors.
2. 18" x 18" std. Special order available. Call Customer Service.
3. N/A with manual recliner.
4. Special measurements are not available.
5. N/A with any power positioning systems.
6. High-torque motors have a half-speed inhibitor. High-speed motors have a quarter-speed inhibitor.
7. Includes appropriate actuator module.
8. N/A with manual recliner.
9. Rear mount accessories are not applicable.

Lowest seat-to-floor height

Synergy Seat:	17.125"
TB Tilt:	18.625"
SL TB Tilt:	16.75"
TB Recline:	17.125"
TB Tilt/Recline:	18.625"
High-back Seat:	21"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Increased widths may increase seat-to-floor height. Contact customer service for details. Listed measurements are in the lowest position. It does not include elevating systems.

**Quantum Custom Delivery Service**

Includes: liftgate service and delivery of your product(s) inside your facility. Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.