

Jazzy 1420

Order Form

500lb wt cap












Pride Mobility Products Corporation
 182 Susquehanna Ave., Exeter, PA 18643
 Phone: (866)800-2002 • Fax: (866) 707-3422

Discount: _____
 Terms: _____
 Freight: _____
 Taken By: _____





Account Number: _____
 Provider Name: _____
 Contact: _____
 Phone: _____
 Fax: _____
 Date: _____
 PO#: _____
 Marked For: _____
 Ship To Address: _____
 City: _____ State: _____ Zip: _____

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

BASE MODEL			EXPANDABLE CONTROLS				
Description	Price	HCPCS Code	Description	Price	HCPCS Code		
<input type="checkbox"/> Jazzy 1420 2HD-C	\$8,700	K0825	<i>All Joysticks or Omni Plus must be used in conjunction with PG Controller and Harness for Expandable Controller</i>				
<input type="checkbox"/> Jazzy 1420 2HD-SS	\$8,700	K0824	<input type="checkbox"/> PG Controller	\$700	E2377		
<input type="checkbox"/> Jazzy 1420 2VHD-C	\$8,700	K0827	<input type="checkbox"/> Harness for Expandable Controller	\$350	E2313		
<input type="checkbox"/> Jazzy 1420 2VHD-SS	\$8,700	K0826	<input type="checkbox"/> Remote Plus 3-Key <input type="checkbox"/> 1 Profile <input type="checkbox"/> Multiple Profile	\$700			
COLOR			<input type="checkbox"/> Remote Plus 3-Key ¹ With Switch Jacks	\$700			
Description	Price	<input type="checkbox"/> Champagne	No Charge	<input type="checkbox"/> Remote Plus 3-Key ¹ With Speed Dial	\$700		
<input type="checkbox"/> Candy Apple Red	No Charge	<input type="checkbox"/> Viper Blue	No Charge	<input type="checkbox"/> Remote Plus 3-Key ¹ With Toggles	\$700		
<input type="checkbox"/> Onyx Black	No Charge	<input type="checkbox"/> Sunburst Orange	No Charge	<input type="checkbox"/> Remote Plus 3-Key, With Environmental Controls <input type="checkbox"/> 1 Profile <input type="checkbox"/> Multiple Profile	\$700		
DRIVE WHEEL			<input type="checkbox"/> Remote Plus 3-Key, With Toggles And Environmental Controls ¹	\$700			
Description	Price	<input type="checkbox"/> 16" Pneumatic	No Charge	<input type="checkbox"/> Remote Plus 3-Key, With Switch Jacks And Environmental Controls ¹	\$700		
<input type="checkbox"/> 14" Flat Free	\$213	BATTERY CHARGER			<input type="checkbox"/> Microdrive Joystick ¹	\$700	
Description	Price	<input type="checkbox"/> 8 Amp Off-board Charger	Standard	NON-EXPANDABLE CONTROLS			
Description	Price	<input type="checkbox"/> Base Mounted Foot Platform Installed	Standard	Description	Price		
<input type="checkbox"/> No Base Mounted Foot Platform	No Charge	<input type="checkbox"/> No Base Mounted Foot Platform	No Charge	<input type="checkbox"/> P&G 4-Key VR2, 60 Amp	Standard		
POWER ELEVATING SEAT			<input type="checkbox"/> Power Elevating Seat ¹	\$975	<input type="checkbox"/> P&G 6-Key VR2, 60 Amp	\$200	
<input type="checkbox"/> Through 6-Key VR2	No Charge	<input type="checkbox"/> Through Toggle	No Charge				
<input type="checkbox"/> Left <input type="checkbox"/> Right		<input type="checkbox"/> Power Elevating Seat ¹	\$1375				
<input type="checkbox"/> Through Remote +		<input type="checkbox"/> Through Remote +					
<p>1. 300 lb wt capacity. Includes 1/4 Speed Inhibit. Raises Seat-To-Floor Height. Not available with a manual recline. Max seat size 20"W x 20"D.</p>			<p>1. Standard with 5 individual drive profiles.</p>				

PROGRAMMER

Description	Price
<input type="checkbox"/> Q-Logic Handheld Programmer 	\$545
<input type="checkbox"/> Q-Logic Programming Cable w/ Disk ¹	\$750
<input type="checkbox"/> P&G Programmer 	\$423

1. Works with Q-Logic EX and NE+ 6-Key only.

JOYSTICK MOUNTS

Description	Price	HCPCS Code
<i>Comes Standard with a Quick Release, Adjustable-Height Joystick Block</i>		
<input type="checkbox"/> Fixed Inline Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	Standard	
<input type="checkbox"/> Swing-away Side: <input type="checkbox"/> Left <input type="checkbox"/> Right Style: <input type="checkbox"/> Inline <input type="checkbox"/> Offset Length: <input type="checkbox"/> 4" Ext. <input type="checkbox"/> 7" Ext. ¹	\$210	E1028
<input type="checkbox"/> Multi-axis Fixed Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	\$100	K0108

1. In iQuantum choose in addition to Swing-Away 4" mount.

JOYSTICK HANDLES

Description	Price	HCPCS Code
<input type="checkbox"/> Standard Carrot Handle	Standard	
<input type="checkbox"/> Mushroom Handle	\$71	E2323
<input type="checkbox"/> T-Bar	\$107	E2323
<input type="checkbox"/> Chin Cup	\$71	E2324 ¹
<input type="checkbox"/> Extended I-Handle	\$71	E2323
<input type="checkbox"/> Soft Ball Handle	\$50	E2323
<input type="checkbox"/> Medium Ball Handle	\$50	E2323
<input type="checkbox"/> Bodypoint 3" Goal Post Handle	\$66	E2323
<input type="checkbox"/> Bodypoint 4" Goal Post Handle	\$66	E2323

1. Use HCPCS code E2323 when not used in conjunction with a chin control.

CONTOURED SEATING

Description	Price
High Back: ¹	
<input type="checkbox"/> Limited Recline High-Back with Headrest & Solid Seat Pan ¹	\$275
Seat Width:	
<input type="checkbox"/> 18" <input type="checkbox"/> 20"	
Seat Depth:	
<input type="checkbox"/> 18" <input type="checkbox"/> 20"	
Upholstery:	
<input type="checkbox"/> Black Vinyl	
Depth Adjustable High Back:	
<input type="checkbox"/> Limited Recline Depth Adjustable High Back w/ Headrest ¹	Standard
Color: <input type="checkbox"/> Gray Vinyl <input type="checkbox"/> Black Vinyl	
Size: <input type="checkbox"/> 18"W x 18-20"D <input type="checkbox"/> 20"W x 18-20"D <input type="checkbox"/> 22"W x 20-22"D	
Armrests:	
<input type="checkbox"/> Flip-Up, Removable	No Charge
<input type="checkbox"/> Flip-Up, Height-Adjustable, Removable	\$85 ea. E0973
Arm pads:	
Left: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length	
Right: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length	

1. Seat comes automatically with a 60" belt. Requires Full Length Arms.

CONTOURED SEATING (CONT'D)

Special Order Seat

_____ "W x _____ "D

Seat size can not exceed 24"W x 24"D

\$450 MSRP

MANUAL RECLINER

Description	Price	HCPCS Code
<input type="checkbox"/> Manual Recliner High Back with Headrest	\$600	E1226
<input type="checkbox"/> Manual Recliner High Back with Headrest and Solid Seat Pan ¹	\$600	E1226

Upholstery:

Black Vinyl

Seat Width:

16" 18" 20"

Seat Depth:

16" 18" 20"

Arm pads:

Left: Full Length Desk Length

Right: Full Length Desk Length

1. Includes black vinyl seatpan. Back available in listed colors.

SPECIALTY SEATING

Description	Price	HCPCS Code
<input type="checkbox"/> Specialty Seating Kit	No Charge	
Seat Width:		
<input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"		
Seat Depth:		
<input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20"		
Back Cane Height:		
<input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/> 22"		
Back Rest Height:		
<input type="checkbox"/> 18" <input type="checkbox"/> 20" <input type="checkbox"/> 22"		
Armrests:		
<input type="checkbox"/> 2-Post, Fixed, Flip Back ¹	No Charge	
<input type="checkbox"/> Flip-Back, Cane Mounted ¹	\$95 ea.	
<input type="checkbox"/> Removable, Adjustable, Quick Height ²	\$275 ea.	E0973
<input type="checkbox"/> Removable, Single-Post, Fixed Height Adjustable Heavy Duty	\$150 ea.	E0973
<input type="checkbox"/> 8-12" Adjustment <input type="checkbox"/> 10-14" Adjustment		
<input type="checkbox"/> Pediatric, Removable ³	\$275	E0973
Color: <input type="checkbox"/> Match Shroud <input type="checkbox"/> Black		
<input type="checkbox"/> Quad Release Lever		
Arm pads:		
Left: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length		
Right: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length		
Backs:		
<input type="checkbox"/> Synergy [®] Shape Back Positioning Back	Various	E2613/E2614
<i>See Synergy Shape Back Order Form For Detailed Pricing, Sizes and Codes</i>		
<input type="checkbox"/> Omit Back	No Charge	
<input type="checkbox"/> Sling Back	\$98	
<input type="checkbox"/> Synergy [®] Back with 1.5" Thick Back Cushion	\$275	
Overall Back: <input type="checkbox"/> 18" H <input type="checkbox"/> 19" H <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H		
Height: <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
Choose Lateral Extensions:		
<input type="checkbox"/> 6" for 13"-16"H Back	\$25	
<input type="checkbox"/> 10" for 17"-18"H Back	\$25	
<input type="checkbox"/> 12" For 19"-22"H Back	\$25	

1. Swing-away Laterals Interfere with flip-back arms. Call for configuration with Shape back.
 2. Adjustable from 11" to 17"H in 1/2" increments.
 3. Height adjustable from 7.75"-9.75" in 1/2" increments. Full length arm pads not available.

SYNERGY® SEATING SYSTEMS

BACK OPTIONS (CONT'D)

Description	Price	HCPCS Code
<input type="checkbox"/> Synergy Seating Frame	\$434	
<input type="checkbox"/> Standard Stem Mounts		
<input type="checkbox"/> 2" Extended Vertical Stem Mounts ¹		
<input type="checkbox"/> Synergy Manual Recline Package ²	\$700	E1226
<input type="checkbox"/> Synergy Manual Tilt Package ³	\$550	K0108
<i>Max. Seat Size 14"W x 14"D, 75 lbs Weight Capacity</i>		

1. Increases seat to floor height by 2". Lowest seat-to-floor height not available.
 2. Not available with Flip Back Arms, Cantilever Arms or Shape Back. Requires 22"H canes. Pivot point 2.5" above seat pan.
 3. In i-Quantum this option is chosen under "Width Package Type"

SEATING SPECIFICATIONS

Seat Width:

<input type="checkbox"/> 10"	<input type="checkbox"/> 11"	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"
<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
<input type="checkbox"/> 20"				

Seat Width Package:

10"W-14"W 14"W-18"W 16"W-20"W

Seat Depth:

<input type="checkbox"/> 10"	<input type="checkbox"/> 11"	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"
<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
<input type="checkbox"/> 20"				

Seat Depth Package:

10"D-14"D 14"D-18"D 15"D-20"D

Back Canes Height:

<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"
<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"
<input type="checkbox"/> 24"	<input type="checkbox"/> 25"			

BACK OPTIONS

Description	Price	HCPCS Code
Synergy Backs:		
<input type="checkbox"/> Synergy® Shape Back <i>Positioning Back</i>	Various	E2613/E2614
<i>See Synergy Shape Back Order Form For Detailed Pricing, Sizes and Codes</i>		
<input type="checkbox"/> Omit Back	No Charge	
<input type="checkbox"/> Sling Back	\$98	
<input type="checkbox"/> Synergy® Back <i>with 1.5" Thick Back Cushion</i>	\$275	
Overall Back <input type="checkbox"/> 18" H <input type="checkbox"/> 19" H <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H		
Height: <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
Choose Lateral Extensions:		
<input type="checkbox"/> 6" for 13"-16"H Back	\$25	
<input type="checkbox"/> 10" for 17"-18"H Back	\$25	
<input type="checkbox"/> 12" For 19"-22"H Back	\$25	
<input type="checkbox"/> 16" For 23"-25"H Back	\$25	
New Design! <input type="checkbox"/> TRU-Comfort™ Back <i>4-way Stretch Material Cover</i>		
Width: <input type="checkbox"/> 16"W <input type="checkbox"/> 18" W <input type="checkbox"/> 20"W	\$595	E2620
Overall Back <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H		
Height ¹ : <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Width: <input type="checkbox"/> 17"W <input type="checkbox"/> 19"W	\$595	E2620
Overall Back <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
Height ² : <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Width: <input type="checkbox"/> 22"W ² <input type="checkbox"/> 24"W ⁴	\$695	E2621
Overall Back <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Height ² :		

1. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".
 2. Available with offset (meso) brackets.

Description	Price
Tru Comfort Color Shroud Option:	
<input type="checkbox"/> Black <i>(Standard color option)</i>	
Metallic Colors:\$100	
<input type="checkbox"/> Champagne	
<input type="checkbox"/> Candy Apple Red	
<input type="checkbox"/> Viper Blue	
<input type="checkbox"/> Sunburst Orange	
<input type="checkbox"/> Pearl Pink	
<input type="checkbox"/> Silver	
<input type="checkbox"/> Pewter Gray	
<input type="checkbox"/> Emerald Green	
<input type="checkbox"/> Graphite Blue	
<input type="checkbox"/> Sapphire Blue	
Solid Colors:\$100	
<input type="checkbox"/> Mariner Blue	
<input type="checkbox"/> Cardinal Red	
<input type="checkbox"/> Spring Green	
<input type="checkbox"/> Solar Yellow	
<input type="checkbox"/> Violet	
<input type="checkbox"/> Tangerine	
<input type="checkbox"/> Kelly Green	
Tilt Back Options:	
<input type="checkbox"/> Standard Back Canes	Standard
<input type="checkbox"/> Omit Push Handles	\$125
<input type="checkbox"/> Offset (Meso) ³	\$195
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)	
<input type="checkbox"/> Inset (Endo) ³	\$195
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)	

3. Dual Post Flip-back Armrests interfere with offset canes.

SYNERGY® ARMRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> Flip Back, Cane Mounted Armrests ¹	No Charge	
<input type="checkbox"/> 2-Post, Flip-back, Height Adjustable, Removable ²	\$150 ea.	E0973
<input type="checkbox"/> Single Post Height Adjustable Removable ³	\$75 ea.	E0973
<input type="checkbox"/> Quad Release Lever		
<input type="checkbox"/> Removable, Single Post, Height Adjustable, Heavy Duty	\$75 ea.	E0973
Height: <input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H - 14"H		
<input type="checkbox"/> Pediatric, Removable ⁴	\$75 ea.	E0973
Color: <input type="checkbox"/> Black <input type="checkbox"/> Match Shroud		
<input type="checkbox"/> Quad Release Lever		
<input type="checkbox"/> Height Adjustable Cantilever Armrests ⁵	\$195	K0020
<i>Angle Adjustable: -30° through +30°</i>		

1. Interferes with Swing-Away laterals.
 2. Height adjustment from 9"-13" in 1" increments. Interferes with swing-away laterals. Not available with a depth less than 15".
 3. Height adjustable from 11" - 17" in 1/2" increments.
 4. Height adjustment from 7.75" - 9.75" in 1/2" increments. Not available with full length armpads.
 5. Height adjustment from 11"-17" depends on cane height. Not available with Flip Forward Back Canes, Cane Clamp Pushbar, or Shape Back. Call for details when selecting after market arm troughs.

ARMPADS

Description	Price
<input type="checkbox"/> Straight, Full Length (2"x14")	Standard
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Straight, Desk Length (2"x10")	No Charge
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Waterfall, Full Length (2"x14")	No Charge
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Waterfall, Desk Length (2"x10")	No Charge
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Falcon Gel, Full Length, (3.5"x14")*	\$160 pr.
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Full Length Gel Armpad Mounting Plate <i>(Required)*</i>	\$56.69 ea.
<input type="checkbox"/> Falcon Gel, Desk Length, (3.5"x10")*	\$125 pr.
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Desk Length Gel Armpad Mounting Plate <i>(Required)*</i>	\$53.52 ea.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

LED LIGHT PACKAGE

Description	Price
<input type="checkbox"/> Through Remote Plus	\$995

LAP BELTS

Description	Price
<input type="checkbox"/> 38" Lap Belt	Standard
<input type="checkbox"/> 50" Lap Belt	Standard
<input type="checkbox"/> 60" Lap Belt	Standard
<input type="checkbox"/> 70" Lap Belt	Standard
<input type="checkbox"/> 80" Lap Belt	Standard
<input type="checkbox"/> 90" Lap Belt	Standard

HEADRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> No Headrest Required		
TRU-Comfort™ Headrest¹:		
<input type="checkbox"/> 8" W Leatherette	\$205	E0955
<input type="checkbox"/> 8"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> 10" W Leatherette	\$205	E0955
<input type="checkbox"/> 10"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> Removable Headrest Bracket	\$210	E1028
Stealth:		
<input type="checkbox"/> Small (9.75"W x 3"H) (STL-OTSM250)	\$205	E0955
<input type="checkbox"/> Large (11"W x 6"H) (STL-OTLG150)	\$205	E0955
<input type="checkbox"/> Pediatric Comfort Plus (6"W) (CPN350)	\$142	E0955
<input type="checkbox"/> Small Comfort Plus (10"W) (CP250)	\$142	E0955
<input type="checkbox"/> Large Comfort Plus (14"W) (CP150)	\$142	E0955
<input type="checkbox"/> Removable Headrest Bracket ² (HMO475P)	\$210	E1028
<input type="checkbox"/> Swing-away (Flip-down) Headrest Bracket (FDM350)	\$210	
Used in Addition to Removable Headrest Bracket		
Aftermarket Bracket:		
<input type="checkbox"/> Stealth's World's Best Removable Headrest Hardware (TWB485)	\$210	E1028

1. Measurement is taken from the inside of the headrest.

NON-POWERED LEGRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> 70° Swing-away Legrests¹	\$162	
<input type="checkbox"/> Straight <input type="checkbox"/> Tapered		
<input type="checkbox"/> Standard Extruded Footplates	No Charge	
<i>Standard Extruded Extensions:</i>		
<input type="checkbox"/> 11"-13" <input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
<i>Standard Extruded Footplates:</i>		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle-Adjustable Footplates ²	\$150	
<i>Angle Adjustable Extensions:</i>		
<input type="checkbox"/> 11"-13" <input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
<i>Angle-Adjustable Footplates:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ⁴ <input type="checkbox"/> 8" x 11" ⁴		
<input type="checkbox"/> Multi-Axis Footplates	\$495	
<i>Multi-Axis Extensions:</i>		
<input type="checkbox"/> 11"-16.5" <input type="checkbox"/> 15.5"-19.5"		
<i>Multi-Axis Footplates:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ⁴ <input type="checkbox"/> 6" x 11" ⁴ <input type="checkbox"/> 8" x 11" ⁵		
<input type="checkbox"/> High Mount, Clamp-On Footrest*		
<input type="checkbox"/> Angle Adjustable		K0040
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<input type="checkbox"/> Multi-Axis		K0040
<input type="checkbox"/> Small, 4" x 7.75"	\$456	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	

1. Legrest Extensions are measured from seatpan to footplate.
 2. Includes Heel Loops
 3. Causes interference on seat widths less than 18".
 4. Causes interference on seat widths less than 20".

NON-POWERED LEGRESTS (CONT'D)

Description	Price	HCPCS Code
<input type="checkbox"/> Heavy Duty Drop-In Legrests^{5,6}	\$420	
<input type="checkbox"/> Angle Adjustable Footplates ²	\$150	
<i>Angle Adjustable Extension Length:</i>		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5"-18"		
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ⁴ <input type="checkbox"/> 6" x 11" ⁴ <input type="checkbox"/> 8" x 11" ⁵		
<input type="checkbox"/> Elevating Legrests⁶	\$325	K0195/E0990
<input type="checkbox"/> Standard Footplates		
<i>Standard Extruded Extension Length:</i>		
<input type="checkbox"/> 11"-14" <input type="checkbox"/> 14"-17"		
<i>Standard Extruded Footplate Size:</i>		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle Adjustable Footplates ⁷	\$150	
<i>Angle Adjustable Extension Length:</i>		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5"-18"		
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ⁸ <input type="checkbox"/> 6" x 11" ⁸ <input type="checkbox"/> 8" x 11" ⁹		
<input type="checkbox"/> Calf Pad Adjustment		
<input type="checkbox"/> 1 Hole <input type="checkbox"/> 3 Hole <input type="checkbox"/> 6 Hole		
<input type="checkbox"/> 90° Pediatric Swing-Away Legrests⁶	\$162	
<input type="checkbox"/> Angle Adjustable Footplate ⁷		
<i>Angle Adjustable Extension Length:</i>		
<input type="checkbox"/> 7"-10"		
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6"		

5. Only available on Synergy Seating.
 6. Legrest Extensions are measured from seatpan to footplate.
 7. Includes Heel Loops.
 8. Causes interference on seat widths less than 18".
 9. Causes interference on seat widths less than 20".

SYNERGY® CUSHIONS

Description	Price	HCPCS Code
<input type="checkbox"/> No Cushion Required		
<i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
Simplicity Cushion General Use		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$98	E2601
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$98	E2602
Solution® 1 Cushion Skin Protection:		
<input type="checkbox"/> 14" x 14" - 20" x 18"	\$380	E2603
<input type="checkbox"/> 20" x 20" - 20" x 24"	\$430	E2603
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2604
Solution® Cushion Skin Protection & Positioning		
<input type="checkbox"/> 10" x 10" - 20" x 18"	\$380	E2607
<input type="checkbox"/> 20" x 20" - 20" x 24"	\$430	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2608
Spectrum® Gel Cushion Positioning		
<input type="checkbox"/> 12" x 10" - 20" x 19"	\$480	E2605
<input type="checkbox"/> 20" x 20" - 20" x 22"	\$530	E2605
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2606
Spectrum® Visco Cushion Skin Protection & Positioning		
<input type="checkbox"/> 10" x 10" - 20" x 18"	\$480	E2607
<input type="checkbox"/> 20" x 20" - 20" x 22"	\$530	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2608
Spectrum® Air Cushion Adjustable Skin Protection		
<input type="checkbox"/> 14" x 12" - 20" x 20"	\$480	E2622
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2623
Spectrum® Air Contour Cushion Adjustable Skin Protection & Positioning		
<input type="checkbox"/> 14" x 12" - 20" x 20"	\$499	E2624
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$549	E2625
TRU-Comfort™ Seat Cushion General Use		
<i>Includes Allure leather-like cover</i>		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$380	E2601
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2602
TRU-Comfort™ Plus Seat Cushion Skin Protection		
<i>Includes 4-way stretch material cover</i>		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2604

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

FRONT RIGGING ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Wheel Bumpers	\$25	
<input type="checkbox"/> Heel Loops (pair) ¹	\$38	E0951
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) ²	\$65	K0108
<input type="checkbox"/> ELR/ALR Adductor Buttons (pair)	\$125	K0108
<input type="checkbox"/> POSAlinc Knee Adductor for Swing-Away Legrests (pair)*	\$256	K0108
<input type="checkbox"/> POSAlinc Rail Mounted Knee Adductor (pair) *	\$287	K0108
<input type="checkbox"/> Residual Limb Support	\$300	E1020
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Bodypoint Ankle Huggers (pair)*	\$134	K0108
<input type="checkbox"/> Small <input type="checkbox"/> Medium		

1. Standard with Angle Adjustable or Multi-Axis footplates.
2. Only available with Swing-Away legrests.

LATERALS (SYNERGY ONLY)

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Thoracic Laterals (pair)^{1*}		
Bracket:		
<input type="checkbox"/> Swing-Away Rail Mounted	\$450	E1028
<i>For Use on Synergy Back and Tru-Comfort</i>		
Pad:		
<input type="checkbox"/> Flat Pad, 3"x5", Leather-Like Cover	\$200	E0956
<input type="checkbox"/> Flat Pad, 4"x6", Leather-Like Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 2"x3", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 3"x5", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 4"x6", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Stealth Pelvic/Thigh Laterals*		
Side:		
<input type="checkbox"/> Left Side <input type="checkbox"/> Right Side		
Bracket:		
<input type="checkbox"/> Removable Bracket (ea.)	\$250	E1028
<input type="checkbox"/> Removable Bracket w/ 1" Offset (ea.)	\$325	E1028
<input type="checkbox"/> Flip Down Bracket (ea.)	\$225	E1028
Pad:		
<input type="checkbox"/> Short Pad, 4"x4" (ea.)	\$90	E0956
<input type="checkbox"/> Medium Pad, 4"x8" (ea.)	\$90	E0956
<input type="checkbox"/> Long Pad, 4"x12" (ea.)	\$90	E0956
Adapter:		
<input type="checkbox"/> 2" Depth/Height Adapter Plate (ea.) ²	\$25	K0108

1. Swing-away Laterals Interfere with flip-up armrests.
2. Height adapter plate is needed w/ Tru Comfort seat cushion. Not available on the flip down bracket.

MISCELLANEOUS POSITIONING COMPONENTS (SYN ONLY)

Description	Price	HCPCS Code
<input type="checkbox"/> Therafin Flip-down Knee Abductor Bracket*	\$172	E1028
Pad:		
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large	\$135	E0957
<input type="checkbox"/> Stealth Foot Boxes (pair)*	\$475	K0108
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large		
<input type="checkbox"/> Stealth Padded Shoe Holders (pair)*	\$166	K0108
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Large		
<input type="checkbox"/> Otto Bock Arm Channels (ea.)*	\$104	E2209
Side:		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
Size:		
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large		
<input type="checkbox"/> Otto Bock Flat Hand Pad(ea.)*	\$52	
Side:		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
Size:		
<input type="checkbox"/> Medium <input type="checkbox"/> Large		

Please see positioning components order form for additional types, sizes and details.

ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Rear Basket ¹	\$103	
<input type="checkbox"/> Cane/Crutch Holder ¹	\$98	E2207
<input type="checkbox"/> Weather Cover, Medium, Black	\$158	
<input type="checkbox"/> Cup Holder	\$17 ²	
Side:		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Saddle Bag	\$28	
<input type="checkbox"/> Flag Holder ¹	\$17	
<input type="checkbox"/> Quantum Backpack ³	\$60	
<input type="checkbox"/> Oxygen Holder Contoured Seating ¹	\$103	E2208
<input type="checkbox"/> Oxygen Holder Synergy Seat ⁴	\$150	E2208
<input type="checkbox"/> Walker Holder ¹	\$92	K0108
<input type="checkbox"/> Tablet Holder ⁵	\$315	
<input type="checkbox"/> Deluxe Cell Phone/Media Holder ⁵	\$140	
<input type="checkbox"/> Cell Phone/Media Holder ⁵	No Charge	

1. Includes mounting hardware. Not available on Manual Recliner.
2. No charge on Synergy Seating.
3. Backpack may interfere w/ other rear accessories. Not available w/ contoured seating.
4. Includes mounting hardware.
5. Only available with Synergy or Specialty seating.

SEAT TO FLOOR HEIGHTS:

Seating System	Low	Middle	Tall
Contoured Seating	22.5"	23.25"	24"
Specialty Seat	17.875"	18.625"	19.375"
Synergy [®] Seat	18.125"	19.125"	20.125"

Seat to floor height are +/- 1/8" with 18"W x 18"D with Flat-Free tires (add 3/8" for pneumatic tires). Increased widths may increase lowest seat to floor heights.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum