

Quantum® 1122 Series
 Order Form / Price List
 Revised: November 2007



Pride Mobility Products Corporation®
 182 Susquehanna Ave., Exeter, PA 18643
 Sales: Phone: (866) 800-2002 • Fax: (866) 707-3422

This product is available for purchase online at www.mypride.com

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:		P.O. #:	Discount:
Address:	Address:			Terms:
Phone Contact:	Marked For:			Freight:
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.				All pricing is MSRP. Taken By:

☐ QUANTUM 1122 3SP-SS (single power) REFERENCE GUIDE #1-3

Quantum 1122 3SP-SS \$6,256
 Group 3 / Single Power
K0856
 300 lbs. weight capacity

Seat Options:
 TB Tilt E1002
 TB Recline E1004

-Flip-back, cane mounted, full length arms
 -70 amp, 6-Key PG VSI controller
 -14" Solid tires
 -8 amp off-board charger
 -60" lapbelt
 -Swing-away legrests
 -Reinforced footplates
 -Tilt through 6-key VSI
 Omit Solid Tires (-\$213)
 Omit 6-Key VSI (-\$200)
 Omit contoured seat (-\$275)

Alternate Seating Options:
 Power Elevating Seat

Select seating , front rigging options in appropriate section & base color.

☐ QUANTUM 1122 3SP-C (single power) REFERENCE GUIDE #2-3

Quantum 1122 3SP-C \$6,256
 Group 3 / Single Power
K0857
 300 lbs. weight capacity

*The following single power options are available at a later date:
 1. Power Tilt
 2. Power Recline

-18" W x 18" D Ltd. Recline Gray vinyl high-seat
 -Full length armrests
 -70 amp 6-key PG VSI
 -14" Solid tires
 -8 amp off-board charger
 -60" lapbelt
 -Foot Platform
 Omit Solid Tires (-\$213)
 Omit 6-Key VSI (-\$200)
 Omit contoured seat (-\$275)

Alternate Seating Options:
 Power Elevating Seat

Select seating , front rigging options in appropriate section & base color.

☐ QUANTUM 1122 4MP-SS (multiple power) REFERENCE GUIDE #3-4

Quantum 1122 4MP-SS \$6,056
 Group 4 / Multiple power
K0884
 300 lbs. weight capacity

-Reclining, full length arms
 70 amp Remote Plus controller E2377 \$700
 -14" Solid tires
 -8 amp off-board charger
 -60" lapbelt
 -Swing-away legrests
 -Reinforced footplates
 Omit Solid Tires (-\$213)
 Omit contoured seat (-\$275)

Seat Options:
 TB Tilt E1002
 TB Recline E1004
 TB Tilt & Recline E1007

Alternate Seating Options:
 Power Elevating Seat

Select seating , front rigging options in appropriate section & base color.

When configuring power chair use reference guides to determine accessories included in the MSRP pricing package

Model number with suffix must be used for all chairs billed through Medicare.

SPECIFIC CONFIGURATION

- Quantum 1122 \$5,843
- Quantum 1122 3SP-SS \$6,256
- Quantum 1122 3SP-C \$6,256
- Quantum 1122 4MP-SS \$6,056
- Omit contoured seat²⁰ (-\$275)

COLOR

- Champagne NC
- Candy Apple Red NC
- Viper Blue NC
- Onyx Black NC
- Sunburst Orange NC

Quantum® 1122 Series

Order Form / Price List
Revised: November 2007

BATTERIES(each)

- Group 24, 70 amp, Gel **E2363** \$352ea
- Group 34, 60 amp, Gel **K0108** \$330ea
Required w/ Super Low Tilt
- NF-22, 50-amp Batteries **E2361**
- Maintenance-free AGM \$255ea
- Maintenance-free Gel Cell \$295ea

Required with power tilt / recline elevating system

BATTERY CHARGER

- 8 amp Off-board Charger Std.
- 5 amp On & 8 amp Off-board Charger \$220

Off-board charger required with super low tilt

NON-EXPANDABLE CONTROLS

- VSI Integral 70 amp^{1,25} NC
- 6-Key VSI Integral 70 amp^{1,3} \$200

EXPANDABLE CONTROLS

All joysticks must be used in conjunction
w/ a PG controller

CONTROLLER OPTIONS

- PG Controller **E2377** \$695

All joysticks or enhanced display must
be used in conjunction a PG controller

JOYSTICK OPTIONS**PG Electronics**

- 3-Key Remote Plus 70 amp^{1,3,22} **E2399** \$700
(Select Box for Profile)
- 1 Profile Multiple Profile

- 7-Key Access. Remote Plus
70 amp^{1,3,4,22} **E2399** \$700

- Remote Plus 70 amp with
Speed Dial^{3,5,22} **E2399** \$700

- Remote Plus 70 amp
with Toggles^{3,5,22} **E2399** \$700

- Remote Plus 70 amp with
Switch Jacks^{3,5,22} **E2399** \$700

- Microdrive, 70 amp^{3,5,6,22} **E2399** \$700

- Omni Plus²⁶ \$1721.63

(Requires a secondary control to
operate the power chair)

Joystick includes electronics, cables, &
junction boxes.

CONTROLLER MOUNT

- Fixed, In-Line Joystick Mount Std.
 Left Right
- Quick Disconnect, Joystick Mount^{23,29} -NC
 Left Right
- Swing-away Joystick Mount²⁴ **E1028** \$196
 4" L² 7" L
 Left Right
- Multi-axis Joystick Mount **K0108** \$100
 Left Right
(Fixed, does not swing-away)
- Vertical Adj. Clamp⁸ **K0108** \$95
Provides vertical joystick adjustment
w/o adjusting armrest height.

PROGRAMMER

- PG Programmer \$423

JOYSTICK HANDLE OPTIONS**PG Controls**

- Std. Knob Std.
- Mushroom Handle **E2323** \$71
- T-Bar **E2323** \$107
- Chin Cup **E2323** \$71
- Extended I-Handle **E2323** \$71
- Soft Ball Handle **E2323** \$50
- Medium Ball Handle **E2323** \$50
- Goal Post Handle **E2323** \$66
 3" W 4" W

DRIVE WHEELS

- 14" Pneumatic (pair)²⁵ NC
- 14" Flat-Free (pair) \$213

Quantum® 1122 Series

Order Form / Price List

Revised: November 2007

Legrest measurements have changed. Lengths are now measured according to industry standard from the seat pan to the back of the footplate

NEW SWING-AWAY LEGRESTS

<input type="checkbox"/> 70° Swing-away footrests Style #8 ²¹	\$162
<input type="checkbox"/> 80° Swing-away footrests Style #8 ²¹	\$162
<input type="checkbox"/> 70° Tapered S/A footrests Style #8 ²¹	\$162
(N/C with power positioning systems)	
<input type="checkbox"/> Quad Release Lever	NC

S/A Standard Extensions:	NC
<input type="checkbox"/> 11"-13" Ext.	NC
<input type="checkbox"/> 13"-16.5" Ext.	Std.
<input type="checkbox"/> 15.5"-19.5" Ext.	NC
<input type="checkbox"/> 18.5"-22" Ext.	NC
(Lengths over 17" may cause interference)	

<input type="checkbox"/> Aluminum Footplate	Std.
<input type="checkbox"/> 7.25" x 6" ¹⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 8.5 x 6" ⁹ (Tapered legs require 18"W or greater seat)	

Angle Adjustable Ext.:	\$150
<input type="checkbox"/> 12" - 14.5" Ext.	
<input type="checkbox"/> 14.5" - 18" Ext.	

(Lengths over 17" may cause interference)

Angle Adjustable Footplates (pair)

<input type="checkbox"/> 4"x6" ³⁰ (Tapered legs require 14"W or greater seat)	
<input type="checkbox"/> 5"x6" ¹⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 5"x8" ¹⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 6"x8" ⁹ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 6"x11" ⁹ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 8"x11" ¹¹ (Tapered legs require 22"W or greater seat)	
<input type="checkbox"/> Wheel Bumpers	\$25

Multi-Axis Ext. (ball & socket):	\$495
<input type="checkbox"/> 11"-16.5" Ext.	
<input type="checkbox"/> 15.5"-19.5" Ext.	

(Lengths over 17" may cause interference)

Multi-Axis Footplates (pair)

<input type="checkbox"/> 4"x6" ³⁰ (Tapered legs require 14"W or greater seat)	
<input type="checkbox"/> 5"x6" ¹⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 5"x8" ¹⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 6"x8" ⁹ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 6"x11" ⁹ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 8"x11" ¹¹ (Tapered legs require 22"W or greater seat)	
<input type="checkbox"/> Wheel Bumpers	\$25

FRONT RIGGINGS

<input type="checkbox"/> Elevating Leg Rests E0990	\$325
Adj. Calf Pads w/ ELRS (1 hole std.)	
<input type="checkbox"/> 1 hole	<input type="checkbox"/> 3 holes* <input type="checkbox"/> 6 holes*

***Potential for interference w/ shrouds**

<input type="checkbox"/> HD Drop-Ins, 70° ^{08,15,20}	\$420
<input type="checkbox"/> HD Drop-Ins, 80° ^{08,15,20}	\$420

N/A w/ contoured or Specialty seating

Angle Extensions:	NC
<input type="checkbox"/> 12" - 14.5" Ext.	
<input type="checkbox"/> 14.5" - 18" Ext.	
<input type="checkbox"/> 17.5" - 21" Ext.	

(Lengths over 17" may cause interference)

FRONT RIGGINGS CONT.

<input type="checkbox"/> Swing-away, Pediatric	\$162
<input type="checkbox"/> 4"x6" <input type="checkbox"/> 5"x6" Angle-Adj. (pair)	
<input type="checkbox"/> Drop-ins, Pediatric	\$162
<input type="checkbox"/> 4"x6" <input type="checkbox"/> 5"x6" Angle-Adj. (pair)	
<input type="checkbox"/> Foot Platform, 9.5" deep ²⁸ (Substitute for S/A)	NC
<input type="checkbox"/> Foot Platform, 11.5" deep ²⁸ (Substitute for S/A)	NC
<input type="checkbox"/> Power-elevating, high-pivot ^{16,17} Legrests, (calf pads included) ¹⁹	E1010 \$1695
<input type="checkbox"/> Independent <input type="checkbox"/> Combined ¹⁶ <input type="checkbox"/> Synchronized ¹⁸	
<input type="checkbox"/> Power-elevating, Articulating ^{7,16,17}	E1010 \$2395
Articulating Range: <input type="checkbox"/> 4" L Std. <input type="checkbox"/> 6" L	NC
Set at: <input type="checkbox"/> 70° <input type="checkbox"/> 80°	
<input type="checkbox"/> Independent <input type="checkbox"/> Combined ¹⁶ <input type="checkbox"/> Synchronized ¹⁸	
<input type="checkbox"/> Swing-away Power-elevating, Articulating ^{7,16,17,27}	E1010 \$2395
Articulating Range: <input type="checkbox"/> 4" L	Std.
Set at: <input type="checkbox"/> 70° <input type="checkbox"/> 80°	
<input type="checkbox"/> Quad Release Lever	NC
<input type="checkbox"/> Independent <input type="checkbox"/> Combined ¹⁶ <input type="checkbox"/> Synchronized ¹⁸	
<input type="checkbox"/> Calf Pads Left & Right (6 x 5.35)	Std.
<input type="checkbox"/> Calf Pads Left & Right (8.00 x 5.86)	NC
<input type="checkbox"/> Power Articulating Foot Platform K0108	\$2,895
<input type="checkbox"/> Synchronized <input type="checkbox"/> Non-Synchronized	
w/ split ang. adj. foot platform	
-8" of Articulation	
- Elevation range 90°-170°	
Adjustment Packages	NC
<input type="checkbox"/> 12"-14" <input type="checkbox"/> 14"-16"	
Select Foot platform for above: <input type="checkbox"/> Adult <input type="checkbox"/> Pediatric	
<input type="checkbox"/> 16"-18" (depth adj) <input type="checkbox"/> 18"-20" (depth adj)	
(Lengths over 17" may cause interference)	
<input type="checkbox"/> Connection Bracket	\$49
<input type="checkbox"/> Heel Cups	\$39
Footplates (pair)	
<input type="checkbox"/> Reinforced Footplate	Std.
<input type="checkbox"/> 7.25" x 6" ¹⁰ <input type="checkbox"/> 8.5 x 6" ⁹	
<input type="checkbox"/> 5"x8" Angle-Adj. Footplate ¹⁰	\$150
<input type="checkbox"/> 6"x8" Angle-Adj. Footplate ⁹	\$150
<input type="checkbox"/> 6"x11" Angle-Adj. Footplate ⁹	\$150
<input type="checkbox"/> 8"x11" Angle-Adj. Footplate ¹¹	\$150
<input type="checkbox"/> Wheel Bumpers (Ang-adj. only)	\$25
Extensions: ¹²	
<input type="checkbox"/> 11"-14" Ext. ²⁷ (\$125 w/ power legs)	NC
<input type="checkbox"/> 14"-17" Ext	Std.
<input type="checkbox"/> 17"-20" Ext	NC
<input type="checkbox"/> High mount Clamp-on; Ang. Adj.	\$326
<input type="checkbox"/> Small <input type="checkbox"/> Large	
<input type="checkbox"/> High mount Clamp-on; Multi-axis	\$412
<input type="checkbox"/> Small <input type="checkbox"/> Large	
<input type="checkbox"/> High mount footplatform ¹⁹ K0108	\$450
<input type="checkbox"/> 9.5" deep <input type="checkbox"/> 11.5" deep	
Lower Extensions	
<input type="checkbox"/> 10"-14" L	
<input type="checkbox"/> 14"-17" L	
<input type="checkbox"/> 17"-20" L	
Rigging Options	
<input type="checkbox"/> Heel Loops (pair) ¹³ E0951	\$20
<input type="checkbox"/> Gel Neoprene Footrest K0108 Sleeve (pair) ¹⁴	\$65
<input type="checkbox"/> ELR / ALR Adductor Buttons (pair) K0108	\$100
<input type="checkbox"/> POSAlinc Knee Adductor (pair) for S/A legrests K0108	\$464
<input type="checkbox"/> Residual limb support left E1020	\$300
<input type="checkbox"/> Residual limb support right E1020	\$300

Footnotes for pages 1 & 2

- Std. w/ one programmable drive profile.
- N/A with waterfall armpads and dual post armrests.
- Needed for controlling power positioning systems.
- Needed for lighting packages.
- Avail. w/ five individual drive profiles.
- Requires multi-axis joystick mount; N/C for mount.
- Std. w/ heel loops & ang-adj. footplates, check size below.
- N/A on Specialty, contoured seating or waterfall armpads.
- Only avail. on 16" W or greater seats.
- Only avail. on 14" W or greater seats.
- Only avail. on 20" W or greater seats.
- Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured or rehab seating w/ cushions is selected. The overall length may vary approx. 3" depending on choice of seat and /or cushion.
- Std. w/ ang-adj. footplates.
- Avail. w/ S/A footrests & 70,80 degree HD drop-ins.
- May cause interference w/ front anti-tips w/ certain seat widths and lower leg lengths. Call Customer Service.
- Requires the addition of an ALM.
- Power legrests are only available on TRU-Balance and Synergy Seating Systems. Only the combined & independent motion can be selected on the Synergy Seat.
- Only avail. on TB Recline.
- Std. calf pad is not depth adjustable.
- Applicable when ordering no seat, specialty, Synergy seat or any power positioning system.
- Please note swing away footrests are standard when ordered with Power Positioning.
- Separately reimbursable through Medicare for multiple power units.
- N/A w/ waterfall armpads; only available on Quantum Seating Systems and only one option maybe mounted.
- Can only mount one armrest accessory in combination w/ S/A joystick.
- If ordered w/ reference guides or configuration packages discounts may apply.
- Please consult specialty controls order form for a variety of secondary controls and mounting options.
- 11" - 14" Ext. are N/A w/ power S/A legrests.
- Fixed bracket, does not elevate.
- N/A w/ 2-post, desk length arms w/ a toggle.
- Only avail. on 12" W or greater seats.

Quantum® 1122 Series

Order Form / Price List
Revised: November 2007

STANDARD SEATING

Medium-Back **NC**

Medium-Back
w/ Solid Seat Pan^{1,2} **NC**

<input type="checkbox"/> Brown Vinyl	<input type="checkbox"/> Beige Vinyl
<input type="checkbox"/> Gray Vinyl	<input type="checkbox"/> Black Vinyl
<input type="checkbox"/> Gray Fabric	<input type="checkbox"/> Brown Fabric

Medium-Back Seating Dimensions

(Dimensions 18"W x 18"D)

Width:	Depth:
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 16" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 18" Std.

UPGRADED SEATING

Ltd. Recline High-Back w/ Headrest **NC**

Ltd. Recline High-Back with Headrest^{1,2}
w/ Solid Seat Pan **NC**

<input type="checkbox"/> Black Vinyl
<input type="checkbox"/> Gray Vinyl
<input type="checkbox"/> Beige Vinyl
<input type="checkbox"/> Gray Fabric
<input type="checkbox"/> Brown Fabric

High-Back Seating Dimensions

(Dimensions 18"W x 18"D)

Width:	Depth:
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 16" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 18" Std.
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 20" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 22" NC

22" W AND ABOVE WILL NOT BE LTD RECLINE.

ARMRESTS

For Medium and High-Back Seats

<input type="checkbox"/> Flip-up, removable	Std.
<input type="checkbox"/> Flip-up, height-adjustable removable ⁴ E0973	\$170

Armpads

<input type="checkbox"/> Desk Length, 10"	NC
<input type="checkbox"/> Left NC	<input type="checkbox"/> Right NC
<input type="checkbox"/> Full Length, 14"	NC
<input type="checkbox"/> Left Std.	<input type="checkbox"/> Right Std.

POWER ELEVATING SEAT^{17,18,19,20}

<input type="checkbox"/> Power Elevating Seat	\$975
<input type="checkbox"/> Through Toggle	NC
<input type="checkbox"/> Through 6-Key VSI ¹⁷	NC
<input type="checkbox"/> Remote Plus	\$1375

MANUAL RECLINER^{2,5}

Manual Recline High Back **E1226** \$600

Manual Recline High Back¹ **E1226**
w/ Solid Seat Pan \$600

(180 degrees of Recline)

(Dimensions 18"W x 18"D)

<input type="checkbox"/> Gray Vinyl	<input type="checkbox"/> Beige Vinyl
Width:	Depth:
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 16" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 18" Std.
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 20" NC

When ordering a lap belt with this seat, an additional 10" is added to accommodate usable length

CUSHIONS⁶

<input type="checkbox"/> Simplicity Cushion		
16 x 16 - 20 x 20	E2601	\$98
<input type="checkbox"/> Solution 1 Cushion		
14 x 14 - 20 x 18	E2603	\$380
20 x 20 - 20 x 24	E2603	\$430
22 x 16 - 24 x 24	E2604	\$430
<input type="checkbox"/> Solution Cushion		
10 x 10 - 20 x 18	E2607	\$380
20 x 20 - 20 x 24	E2607	\$430
22 x 16 - 24 x 24	E2608	\$430
<input type="checkbox"/> Spectrum Cushion		
<input type="checkbox"/> Gel		
10 x 10 - 20 x 18	E2605	\$480
20 x 20 - 20 x 24	E2605	\$530
22 x 16 - 24 x 24	E2606	\$530
<input type="checkbox"/> Visco		
10 x 10 - 20 x 18	E2607	\$480
20 x 20 - 20 x 24	E2607	\$530
22 x 16 - 24 x 24	E2608	\$530
<input type="checkbox"/> Structure 2 Cushion		
14 x 16 - 20 x 18	K0734	\$480
20 x 20 - 20 x 24	K0734	\$580
22 x 18 - 24 x 24	K0735	\$580

SPECIALTY SEATING

Specialty Seat Kit (limited adjustability) **NC**

SYSTEM SPECIFICATIONS

Specialty Seating Kit Width

<input type="checkbox"/> 16" NC
<input type="checkbox"/> 18" Std.
<input type="checkbox"/> 20" NC

Specialty Seating Kit Depth

<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	

Specialty Seating Kit Back Cane Height

<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 20" NC
<input type="checkbox"/> 22" NC	

Back Selection

See Synergy Shape Back Order Form^{11,12,13}

Solid Plate Back up to 20"W²¹ \$275

See Positioning Components Order Form for additional backs

Back Height Solid Pan¹⁴

<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

Sling Back

<input type="checkbox"/> 16" wide	\$98
<input type="checkbox"/> 18" wide	\$98
<input type="checkbox"/> 20" wide	\$98

Sling Back Height

(Hgt must match cane hgt.)

<input type="checkbox"/> 18" height	NC
<input type="checkbox"/> 20" height	NC
<input type="checkbox"/> 22" height	NC

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

ARMRESTS

<input type="checkbox"/> 2-Post, fixed, flip-back ^{13,16}	Std.
<input type="checkbox"/> Flip-back, cane mounted ^{13,16}	\$95ea
<input type="checkbox"/> Removable, adjustable, quick height ^{15,22} E0973	\$275ea
11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments	
<input type="checkbox"/> Removable, height adjustable ²²	\$150ea
heavy duty armrests E0973	
<input type="checkbox"/> 8"H-12"H	<input type="checkbox"/> 10"H-14"H

Armrest Pads

<input type="checkbox"/> Full Length, straight, 2"x14"	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Desk Length, straight, 2"x10"	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Full Length, Waterfall, 2"x14" ⁹	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Desk Length, Waterfall, 2"x10" ^{9,10}	NC

HEADRESTS³

<input type="checkbox"/> Stealth ^{7,8} E0955	
Large (11" wide x 6" high) \$210	
<input type="checkbox"/> Stealth ^{7,8} E0955	
Small (9 3/4" wide x 3" high) \$210	
<input type="checkbox"/> Stealth, swing-away (Flip-Down)	
Head Rest Bracket E1028	\$150

Interfaces w/ any headrest that uses 1/2" square down-tube.

<input type="checkbox"/> Medical Applications ^{7,8} E0955	
Large (11" wide x 6" high) \$275	

<input type="checkbox"/> Medical Applications ^{7,8} E0955	
Small (9 3/4" wide x 3" high) \$275	

Fixed, Offset Mounting Bracket **NC**

Footnotes for pages 3:

- Includes black vinyl seat pan, back available in all listed colors.
- 18" x 18" std. Special measurements are not available.
- Only available w/ solid backs.
- N/A on manual recliner.
- Rear mount accessories are N/A.
- See Synergy Cushion order forms for detailed pricing, sizes & codes.
- These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
- Comes with swing-away / removable hardware; HCPCS code is E1028.
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- S/A joystick mount N/A.
- 14" & 22" width are N/A.
- Micro-drive, Shape Back & high-mount foot platform are N/A in the same configuration.
- Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- S/A laterals interfere w/ flip-back arms.
- Includes appropriate actuator module.
- PES is N/A w/ manual recliner.
- STF will increase approx 1.5"
- Includes 1/4 speed inhibit.
- 20" D is N/A.
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.

Quantum® 1122 Series

Order Form / Price List

Revised: November 2007

Synergy Seating Systems¹¹

Synergy Seating™ System \$434
(Price with Base)

Synergy Manual Recline Pkg.^{12,13,14,15,23} \$700
Maximum weight 300 lbs. E1226

Synergy Manual Tilt Pkg.^{13,16} \$550
(Max Dimensions 18"W x 18"D)
Maximum weight 160 lbs. K0108

Synergy Manual Systems come standard with mercury drive lock-out and are in addition to the MSRP of the Synergy Seat.

Synergy Quick Release²¹ \$200
(Max Dimensions 20"W x 18"D)

Synergy Quick Release is in addition to the MSRP of the Synergy Seat.

WIDTH

Set Seat Width At:

<input type="checkbox"/> 10" NC	<input type="checkbox"/> 11" NC
<input type="checkbox"/> 12" NC	<input type="checkbox"/> 13" NC
<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	

Seat Width Package

<input type="checkbox"/> 10"W-14"W	<input type="checkbox"/> 14"W-18"W
<input type="checkbox"/> 16"W-20"W	<input type="checkbox"/> 20"W-24"W

Seat Width Package Manual Tilt

<input type="checkbox"/> 10"W-12"W	<input type="checkbox"/> 12"W-14"W
<input type="checkbox"/> 14"W-16"W	<input type="checkbox"/> 16"W-18"W

DEPTH

Set Seat Depth At:

<input type="checkbox"/> 10" NC	<input type="checkbox"/> 11" NC
<input type="checkbox"/> 12" NC	<input type="checkbox"/> 13" NC
<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	

Seat Depth Package

<input type="checkbox"/> 10"D-14"D	<input type="checkbox"/> 14"D-18"D
<input type="checkbox"/> 16"D-20"D	<input type="checkbox"/> 20"D-24"D

Seat Depth Package Manual Tilt

<input type="checkbox"/> 10"D-14"D	<input type="checkbox"/> 14"D-16"D
<input type="checkbox"/> 16"D-18"D	

BACK CANES

Synergy Manual Recline comes standard with 22" non-height adjustable back canes.

Back Cane Height

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

BACK OPTIONS

See Synergy Shape Back Order Form^{1,2,3}

Solid Plate Back up to 20" W²³ \$275

Solid Plate Back 21"-24" W²³ \$325

Back Height:⁴

18"H 19"H 20"H 21"H

22"H 23"H 24"H 25"H

Sling Back²⁰ \$98

*Back height will be the same as back can height.
Includes gap for cushion.*

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

TRU-Comfort Back¹²

Back Height: 22" 25"

Back up to 20" W E2620 \$595

22" W E2621 \$695

TRU-Comfort Cushion

16 x 16 - 20 x 18 E2601 \$380

20 x 20 - 20 x 22 E2601 \$380

22 x 16 - 22 x 24 E2602 \$480

See TRU-Comfort Back & Cushion order forms for specific sizes

ARMRESTS

Flip-back, cane mounted^{3,5,18} Std.

2-Post, flip-back, height^{3,5,19} E0973 \$150ea
adjustable, removable armrests

(9"-13" of height adjustment in 1" increments)

Single Post Height adj.^{6,22} E0973 \$75ea

11"-17" of height adj. 3" on top in 1/2"

increments & 3" on bottom in 1" increments

Quad Release Lever NC

Removable, height adjustable \$75ea

heavy duty armrests²² E0973

8"H-12"H 10"H-14"H

Pediatric removable²² \$75ea

adjustable, quick height. E0973

Match Shroud Color Black

Quad Release Lever NC

(7 3/4" - 9 3/4" height adjustment in 1/2

increments from seat pan to top of armrest)

ARMPADS

Full length, Straight, 2" X 14" Std.

Left Right

Desk Length, Straight, 2" X 10" Std.

Left Right

Full Length, Waterfall, 2" X 14"⁷ NC

Left Right

Desk Length, Waterfall, 2" X 10"^{7,17} NC

Left Right

HEAD REST⁸

Stealth^{9,10} E0955

Large (11" wide x 6" high) \$210

Stealth^{9,10} E0955

Small (9 3/4" wide x 3" high) \$210

Stealth, swing-away (Flip-Down)

Head Rest Bracket E1028 \$150

Interfaces w/ any headrest that uses 1/2" square down-tube.

Medical Applications^{9,10} E0955

Large (11" wide x 6" high) \$275

Medical Applications^{9,10} E0955

Small (9 3/4" wide x 3" high) \$275

Fixed, Offset Mounting Bracket NC

TRU-Comfort^{9,10} E0955

8" Wide w/ Cobra hardware \$417

10" Wide w/ Cobra hardware \$417

Bracket for after market headrest \$180

Footnotes for pages 4:

- 14" & 22" width are N/A.
- Micro-drive, Shape Back & high-mount foot platform are N/A in the same configuration.
- Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- S/A laterals interfere w/ flip-back arms.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- Available w/ solid backs only.
- Comes with swing-away / removable hardware; HCPCS code is E1028.
- These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
- Manual systems do not include a headrest.
- N/A w/ flip-back, cane mounted arms.
- Solid back not included.
- Includes removable, adjustable, quick-height armrests.
- N/A w/ Shape Back.
- Manual Tilt w/ Shape back & oxygen tank holder requires the use of a pushbar.
- N/A w/ S/A joystick mount.
- N/A w/ Synergy manual recline.
- May cause interference with certain dimension configurations. Call Quantum Customer Service.
- N/A w/ Synergy manual systems.
- N/A w/ manual systems or power elevating seat.
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- Solid plate back is N/A w/manual recline.

Quantum® 1122 Series

Order Form / Price List
Revised: November 2007

RECOMMENDED INFORMATION

Height: _____
Weight: _____
Leg Length: _____
Lower Leg Length: _____

TRU-BALANCE 2 TILT

TRU-Balance 2 Power Tilt E1002 \$4995
55 degree center of gravity (CG) power tilt seating system. 300lbs weight capacity.

TRU-Balance \$1895
Power-Elevating Seating Option¹
7" of Lift & 1/4 Speed Inhibit- 300 lbs. max. wt.
Remote Plus electronics only

TRU-BALANCE 2 RECLINE

TRU-Balance Power Recline E1004 \$5595
Full Drive lock-out at 20 degrees
300 lbs. weight capacity

TRU-Balance \$1895
Power-Elevating Seating Option¹
7" of Lift & 1/4 Speed Inhibit- 300 lbs. max. wt.

TRU-BALANCE 2 TILT & RECLINE

TRU-Balance 2 Power Tilt & Recline E1007 \$9595

Full Drive lock-out at 20 degrees
300 lbs. weight capacity
55 degree center of gravity (CG) power tilt seating system.

TRU-Balance \$1895
Power-Elevating Seating Option¹
7" of Lift & 1/4 Speed Inhibit- 300 lbs. max. wt.

SYSTEM SPECIFICATIONS

Tower Selection (Tilt or Tilt/Recline only)

- Adjustable-angle Std.
Will raise minimum STF by 7/8".
 Fixed NC

Set Seat Width At:

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" NC | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC | |

Seat Width Package

- 14"W-18"W 16"W-20"W
 20"W-24"W

Set Seat Depth At:

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" Std. | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" NC | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC | |

Seat Depth Package

- 14"D-18"D 16"D-20"D
 20"D-24"D

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

Set Back Cane Height At:

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" Std. | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" NC | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC | <input type="checkbox"/> 25" NC |

Seat Overall Back Height³

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" NC. | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC | <input type="checkbox"/> 25" NC |

BACKS

See Synergy Shape Back Order Form⁶

- Solid Plate Back⁴ Std.
 Attendant manual recline^{4,11,24,25,26} \$750
180 degrees of recline with 2 1/2" raised pivot point. Controlled through dual gas struts. 300 lb. weight capacity.

(For use on TB Tilt Systems only)

- TRU-Balance Reclining Back Std.
(For use on Recline Systems only)

TRU-Comfort Back⁴

- Back Height: 22" 25"
 Back up to 20" W E2620 \$595
 22" W E2621 \$695
 TRU-Comfort Cushion
16 x 16 - 20 x 18 E2601 \$380
20 x 20 - 20 x 22 E2601 \$380
22 x 16 - 22 x 24 E2602 \$480
 Omit Back -(\$100)

See TRU-Comfort Back & Cushion order forms for specific sizes

ARMRESTS

- Flip-back, cane mounted^{5,9,14} Std.
Mounted on back canes with full width spreader bar.
 2-Post, flip-back, hgt. adj, removable^{5,9} NC
(9"-13" of height adjustment in 1" increments)
 Single Post Height adj.^{15,16} NC
11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments
 Quad Release Lever NC
 Removable, single-post,¹⁶ height adjustable, heavy duty NC
 8"-12" Adjustment
 10"-14" Adjustment
 Pediatric, removable,¹⁶ adjustable, quick height NC
(7 3/4" - 9 3/4" height adjustment in 1/2 increments from seat pan to top of armrest)
 Quad Release Lever NC
 Match Shroud Kid Color Black
 Reclining, Flip-back⁹ Std.
 8"-9" H 10"-11" H

(For use on Recline Systems only)

Armrest Pads

- Full Length, straight, 2"x14" NC
 Left Right
 Desk Length, straight, 2"x10" NC
 Left Right
 Full Length, Waterfall, 2"x14"¹⁷ NC
 Left Right
 Desk Length, Waterfall, 2"x10"^{2,17} NC
 Left Right

HEAD REST

- Stealth^{18,19} E0955
Large (11" wide x 6" high) \$210
 Stealth^{18,19} E0955
Small (9 3/4" wide x 3" high) \$210
 Stealth, swing-away (Flip-Down)
Head Rest Bracket E1028 \$150
Interfaces w/ any headrest that uses 1/2" square down-tube.
 Medical Applications^{18,19}
Large (11" wide x 6" high) E0955 \$275
 Medical Applications^{18,19}
Small (9 3/4" wide x 3" high) E0955 \$275
 Omit Head Rest - To be supplied by provider
 Fixed, Offset Mounting Bracket NC
 Reclining Head Rest (6.5" x 7")^{18,19} E0955 \$180
(For use on Recline Systems only)
 TRU-Comfort^{18,19} E0955
 8" Wide w/ Cobra hardware \$417
 10" Wide w/ Cobra hardware \$417
 Bracket for after market headrest \$180

Quantum® 1122 Series

Order Form / Price List
Revised: November 2007

POWER POSITIONING ELECTRONICS**Select Electronics to be used**

- PG Remote Plus
 PG 4-Key VSI
 PG 6-Key VSI

SINGLE POWER ELECTRONICS**Power Tilt or Power Recline**

- Through Toggle Box²⁰ (Tilt only) NC
 Left Right (mounts to side of armrest)
- Through 6-Key VSI (Tilt only) NC
- Through Single Switch (1/8" Jack) Tilt Module^{21,22} (Tilt only) \$500
- Through 2-Directional Control Box w/ dual Switch Jacks (1/8" Jack)^{21,22} \$1175
- Through Remote Plus Controller E2310 \$1175

Multiple Motor/Actuator Electronics*:

* At least one TB2 seat function plus power riggings

- Operation Remote Plus Controller E2311 \$2375
- Operation through toggle¹⁴ \$2375
(Select one)
 4-way Toggle
 4-way Button Toggle
 Quad Push-Button Control
Toggle box mounts to side of armrest.
 Left Right

MULTIPLE POWER ELECTRONICS^{23,29}**Power Tilt & Recline**

- Operation Remote Plus Controller E2311 \$2375
- Operation through toggle¹⁴ \$2375
(Select one)
 4-way Toggle
 4-way Button Toggle
 Quad Push-Button Control
Toggle box mounts to side of armrest.
 Left Right

POSITIONING COMPONENTS**AES Laterals**

- Thoracic Curved-S^{18,19} \$532pr
 Thoracic Curved-M^{18,19} \$532pr
 Thoracic Curved-L^{18,19} \$532pr
 Thoracic Tapered-S^{18,19} \$532pr
 Thoracic Tapered-M^{18,19} \$532pr
 Thoracic Tapered-L^{18,19} \$532pr
 Hip & Thigh Removable Short^{18,19} \$474pr
 Hip & Thigh Removable Long^{18,19} \$474pr
 Thoracic S/A Curved Rail-S^{18,19,27} \$532pr
 Thoracic S/A Curved Rail-M^{18,19,27} \$532pr
 Thoracic S/A Curved Rail-L^{18,19,27} \$532pr

Stealth

- Lateral Thoracic Curved-S/A-Pediatric \$650pr
Flush to 1" Offset, Back cane mount^{11,18,19,28}
- Lateral Thoracic Curved-S/A-Medium \$650pr
Flush to 1" Offset, Back cane mount^{11,18,19,28}
- Lateral Thoracic Curved-S/A-Large \$650pr
Flush to 1" Offset, Back cane mount^{11,18,19,28}

ACCESSORIES

- Rear Basket^{6,7,12} \$103
 Cane/Crutch Holder^{6,7,12} \$98
 Oxygen Holder^{6,7,12,13} E2208 \$103
 Oxygen Holder^{6,7,10,11,13} E2208 \$150
(TRU-Balance 2 Tilt & Synergy Seat only)
 Stroller handle \$85
(Synergy Manual Tilt only)
 Accessory Mounting Bracket \$50
(1" adaptor tube for use w/ TB Power Recline)
 Walker Holder^{2,6,7,12} \$92
 Flag Holder (SOF on rehab seating) \$17
 Weather Cover, Large, Black \$202
 Saddle bag(contoured seating only) \$28
 Cup holder \$17
 Left Right

N/A on desk length arm with S/A mount or toggle box.

- Lap Belt Std.
 Pediatric 38" N/A w/ limited recline seating
 50"⁸ 60" 70" 80" 90"

Quantum Custom Delivery Service

Includes: liftgate service and delivery of your product(s) inside your facility
Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.

Footnotes for pages 5 & 6:

- STF will increase.
- S/A joystick mount N/A.
- Measured from seat pan to top of back pan w/ 3-4" gap on tilt & 7-8" gap on recline to accommodate the cushion.
- N/A w/ flip-back, cane mounted arms.
- Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
- Includes mounting bracket and hardware.
- N/A on reclining seats.
- N/A on manual recline high-back.
- S/A laterals interfere w/ flip-back arms.
- Manual Tilt w/ Shape Back and oxygen holder requires the use of a pushbar.
- S/A laterals interfere w/ flip back arms.
- N/A w/ tilt systems.
- N/A w/ flip-back arms.
- N/A w/ Attendant manual recline.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- Comes with swing-away / removable hardware; HCPCS code is E1028.
- These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
- N/A w/ waterfall armpads.
- Consult specialty control order form for a wide variety of switch types and sizes-- 2 required.
- Does not include switch.
- Multiple motor controller options are only available w/ a PES.
- N/A w/ Shape Back.
- Includes price for removable, adjustable, quick-height arms & 2-post arms.
- Solid plate back is N/A w/manual recline.
- For use w/ recline systems only.
- For use with Synergy Seat and TB Tilt only.
- Tilt and elevate combo will not have an attendant control

Lowest seat-to-floor height

Specialty Seat:	16.875"
Synergy Seat:	17.125"
TB 2 Tilt:	17.5"
TB 2 Recline:	17.125"
TB 2 Tilt/Recline:	17.5"
High-back Seat:	21"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers (please see Tilt information above), 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires).

Increased widths may increase seat-to-floor height. Contact customer service for details. Listed measurements are in the lowest position. It does not include elevating systems.