


**Quantum® 6000 Series
Group 3**
Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.
Pride Mobility Products Corporation®
182 Susquehanna Ave., Exeter, PA 18643
Sales: Phone: (866) 800-2002 • Fax: (866) 707-3422



This product is available for purchase online at www.mypride.com

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:		P.O. #:	Discount:
Address:	Address:			Terms:
Phone Contact:	Marked For:			Freight:
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.				Taken By:

All pricing is MSRP.

☐ QUANTUM 6000 3SP-SS (single power) REFERENCE GUIDE #1-3

Quantum 6000 3SP-SS \$8,421 -Flip-back, cane mounted, full length arms
Group 3 / Single Power -100 amp, Remote Plus controller
K0856 -14" Solid tires
300 lbs. weight capacity -8 amp off-board charger
-60" lapbelt
Seat Options:
-Swing-away legrests
-Reinforced footplates
 TB Tilt E1002 Omit Solid Tires (-\$213)
 TB Recline E1004 Omit contoured seat (-\$275)
Select seating, front rigging options in appropriate section & base color.

☐ QUANTUM 6000 3SP-C (single power) REFERENCE GUIDE #2-3

Quantum 6000 3SP-C \$8,421 -18" W x 18" D Ltd. Recline Gray vinyl high seat
Group 2 / Single Power -Full length armrests
K0857 -100 amp, Remote Plus controller
300 lbs. weight capacity -14" Solid tires
*For use when a single -8 amp off-board charger
power option is added at a later date. -60" lapbelt
-Foot Platform
 Omit Solid Tires (-\$213)
 Omit contoured seat (-\$275)
Alternate Seating Options:
 Power Elevating Seat
Select seating, front rigging options in appropriate section & base color.

☐ QUANTUM 6000 3MP-SS (multiple power) REFERENCE GUIDE #3-3

Quantum 6000 3MP-SS \$8,421 -Flip-back, cane mounted, full length arms
Group 3 / Multiple Power -100 amp, Remote Plus controller
K0861 -14" Solid tires
300 lbs. weight capacity -8 amp off-board charger
-60" lapbelt
Seat Options: -Swing-away legrests
-Reinforced footplates
 TB Tilt E1002 Omit Solid Tires (-\$213)
 TB Recline E1004 Omit contoured seat (-\$275)
 TB Tilt & Recline E1007
Select seating, front rigging options in appropriate section & base color.

When configuring power chair use reference guides to determine accessories included in the MSRP pricing package

Model number with suffix must be used for all chairs billed through Medicare.

SPECIFIC CONFIGURATION

- Quantum 6000 \$8,208
- Quantum 6000 3SP-SS \$8,421
- Quantum 6000 3SP-C \$8,421
- Quantum 6000 3MP-SS \$8,421
- Omit contoured seat²³ (-\$275)

- High-Speed Motor Package \$1,574

*Programmable speed up to 8 MPH, Hammer motors
300 lbs. weight capacity 1/4 speed inhibit with
all power positioning seating systems.
Power elevating seat N/A.
Range will decrease w/ High-Speed Package.*

FRAME COLOR

These colors are on a painted metal frame.

- Metallic Red NC
- Metallic Blue NC
- Sunburst Orange NC
- Metallic Black NC
- Metallic Silver NC

*Due to the nature of the painted frame design,
color substitutions may not be made in the field.*

BATTERIES (Installed)

- Group 24, 70 amp, Gel E2363 \$315
- NF-22 Battery Required w/ Vent Pkg.²¹
(Select below)
- NF-22, 50-amp Batteries E2361
 - Maintenance-free AGM \$225
 - Maintenance-free Gel Cell \$265
- Required with power elevating seat**
- Group 34, 60 amp, Gel K0108 \$300
Required w/ Super Low Tilt

Quantum® 6000 Series Group 3

Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

BATTERY CHARGER

8 amp Off-board Charger Std.

CONTROLLER

PG Electronics

3-Key Remote Plus 100 amp^{3,5} Std.
 7-KeyAccess. Remote Plus
 100 amp^{3,4,5} \$275
 Remote Plus 100 amp with
 Speed Dial^{3,5} \$275
 Remote Plus 100 amp
 with Toggles^{3,5} \$275
 Remote Plus 100 amp with
 Switch Jacks^{3,5} \$275
 Microdrive, 100 amp^{3,5,6} \$275
 Omni Plus²⁸ \$1721.63

Q-Logic Electronics

80 amp Q-Logic Hand Control^{3,5} \$370
 80 amp Q-Logic
 Stand Alone Control^{3,7} \$370
 Q-logic Enhanced Display²⁸ \$1720

PROGRAMMER

Remote Plus Programmer \$423
 Q-Logic Programmer \$423

CONTROLLER MOUNT

Fixed, In-Line Joystick Mount Std.
 Left Right
 Quick Disconnect, Joystick Mount²⁵ NC
 Left Right
 Swing-away Joystick Mount²⁶ E1028 \$196
 4" L²⁹ 7" L
 Left Right
 Multi-axis Joystick Mount K0108 \$100
 Left Right
 (Fixed, does not swing-away)
 Vertical Adj. Clamp⁸ K0108 \$95
 Provides vertical joystick adjustment
 w/o adjusting armrest height.

JOYSTICK HANDLE OPTIONS

PG & Q-Logic Controls

Std. Knob Std.
 Mushroom Handle E2323 \$71
 T-Bar E2323 \$107
 Chin Cup E2323 \$71
 Extended I-Handle E2323 \$71
 Soft Ball Handle E2323 \$50
 Medium Ball Handle E2323 \$50
 Goal Post Handle E2323 \$66
 3" W 4" W

DRIVE WHEELS

14" Pneumatic (pair)²⁷ NC
 14" Flat-Free (pair) \$213

FRONT RIGGINGS

Swing-away footrests²⁴ Std.
 (\$162 w/ captains seating)
 Elevating Leg Rests E0990 \$325
 Adj. Calf Pads w/ ELRS (1 hole std.)
 1 hole 3 holes* 6 holes*
 *Potential for interference w/ shrouds
 HD Drop-Ins, 70°^{8,15,20} K0108 \$420
 (HD Drop-ins N/A w/ contoured or Specialty seating)
 Swing-away, Pediatric K0040 \$162
 4"x6" 5"x6" Angle-Adj. (pair)
 Drop-ins, Pediatric K0040 \$162
 4"x6" 5"x6" Angle-Adj. (pair)
 Foot Platform (Substitute for S/A) NC
 No Front Rigging NC

Power-elevating, high-pivot^{16,17,22} E1010 \$1695
 Legrests, (calf pads included)¹⁹
 Independent Combined¹⁶ Synchronized¹⁸
 Power-elevating, Articulating^{16,17,20,22} E1010 \$2395
 Articulating Range: 4" L Std. 6" L NC
 Set at: 70°
 Independent Combined¹⁶ Synchronized¹⁸
 Swing-away Power-elevating, E1010
 Articulating^{16,17,20,22} \$2395
 Articulating Range: 4" L Std.
 Set at: 70°
 Quad Release Lever NC
 Independent Combined¹⁶ Synchronized¹⁸
 Calf Pads Left & Right (6 x 5.35) Std.
 Calf Pads Left & Right (8.00 x 5.86) NC

Footplates (pair)

Reinforced Footplate Std.
 7.25" x 6" 8.5 x 6"
 5"x8" Angle-Adj. Footplate¹⁰ K0040 \$150
 6"x8" Angle-Adj. Footplate⁹ K0040 \$150
 6"x11" Angle-Adj. Footplate⁹ K0040 \$150
 8"x11" Angle-Adj. Footplate¹¹ K0040 \$150
 Wheel Bumpers (Ang-adj. only) \$25

Extensions:⁴²

11"-14" Ext.² NC
 14"-17" Ext Std.
 High mount Clamp-on; Ang. Adj. K0040 \$326
 Small Large
 High mount Clamp-on; Multi-axis K0040 \$412
 Small Large
 High mount footplatform⁹ K0108 \$450
 Mid-wheel 6 foot platform
 Lower Extensions
 10"-14" L
 14"-17" L

Rigging Options

Heel Loops (pair)¹³ E0951 \$20
 Gel Neoprene Footrest K0108
 Sleeve (pair)¹⁴ \$65
 ELR / ALR Adductor Buttons (pair) K0108 \$100
 POSAlinc Knee Adductor (pair)
 for S/A legrests K0108 \$464
 Residual limb support left E1020 \$300
 Residual limb support right E1020 \$300

Footnotes for pages 1 & 2

- Std. w/ one programmable drive profile.
- Special order on on ALR's.
- Needed for controlling power positioning systems.
- Needed for lighting packages.
- Avail. w/ five individual drive profiles.
- Requires multi-axis joystick mount; N/C for mount.
- Avail. w/ four profiles.
- N/A on Specialty, contoured seating or waterfall armpads.
- Only avail. on 16" W or greater seats.
- Only avail. on 14" W or greater seats.
- Only avail. on 20" W or greater seats.
- Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured or rehab seating w/ cushions is selected. The overall length may vary approx. 3" depending on choice of seat and /or cushion.
- Std. w/ ang-adj. footplates.
- Avail. w/ S/A footrests & 70 degree HD drop-ins.
- May cause interference w/ front anti-tips w/ certain seat widths and lower leg lengths. Call Customer Service.
- Requires Remote Plus.
- Power legrests are only available on TRU-Balance and Synergy Seating Systems. Only the combined & independent motion can be selected on the Synergy Seat.
- Only avail. on TB Recline.
- Std. calf pad is not depth adjustable.
- Std. w/ heel loops & ang-adj. footplates, check size below.
- Vent configuration consists of (2) Group 24 batteries and (1) NF-22.
- Power legrests are N/A w/ Q-Logic controls.
- Applicable when ordering no seat, specialty, Synergy seat or any power positioning system
- Please note swing away footrests are standard when ordered with Power Positioning.
- N/A w/ waterfall armpads; only available on Quantum Seating Systems and only one option maybe mounted.
- Can only mount one armrest accessory in combination w/ S/A joystick.
- If ordered w/ reference guides or configuration packages discounts may apply.
- Please consult specialty controls order form for a variety of secondary controls and mounting options.
- N/A with waterfall armpads and dual post armrests.

Quantum® 6000 Series Group 3

Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

STANDARD SEATING

Medium-Back NC

Medium-Back
w/ Solid Seat Pan^{1,2} NC

- Brown Vinyl
- Beige Vinyl
- Gray Vinyl
- Black Vinyl
- Gray Fabric
- Brown Fabric

Medium-Back Seating Dimensions

(Dimensions 18"W x 18"D)

Width:		Depth:	
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 16"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 18"	Std.

UPGRADED SEATING

Ltd. Recline High-Back w/ Headrest NC

Ltd. Recline High-Back with Headrest^{1,2}
w/ Solid Seat Pan NC

- Black Vinyl
- Gray Vinyl
- Beige Vinyl
- Gray Fabric
- Brown Fabric

High-Back Seating Dimensions

(Dimensions 18"W x 18"D)

Width:		Depth:	
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 16"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 18"	Std.
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 20"	NC
<input type="checkbox"/> 22"	NC		

22" W AND ABOVE WILL NOT BE LTD RECLINE.

ARMRESTS

For Medium and High-Back Seats

- Flip-up, removable Std.
- Flip-up, height-adjustable \$170
- removable⁴ E0973

Armpads

- Desk Length, 10" NC
- Left NC Right NC
- Full Length, 14" NC
- Left Std. Right Std.

MANUAL RECLINER^{2,5}

Manual Recline High Back E1226 \$600

Manual Recline High Back¹
w/ Solid Seat Pan E1226 \$600

(180 degrees of Recline)

(Dimensions 18"W x 18"D)

- Gray Vinyl Beige Vinyl
- | Width: | | Depth: | |
|------------------------------|------|------------------------------|------|
| <input type="checkbox"/> 16" | NC | <input type="checkbox"/> 16" | NC |
| <input type="checkbox"/> 18" | Std. | <input type="checkbox"/> 18" | Std. |
| <input type="checkbox"/> 20" | NC | <input type="checkbox"/> 20" | NC |

When ordering a lap belt with this seat, an additional 10" is added to accommodate

CUSHIONS⁶

<input type="checkbox"/> Simplicity Cushion		
16 x 16 - 20 x 20	E2601	\$98
<input type="checkbox"/> Solution Cushion		
10 x 10 - 20 x 18	E2607	\$380
20 x 20 - 20 x 20	E2607	\$430
22 x 16 - 22 x 20	E2608	\$430
<input type="checkbox"/> Spectrum Cushion		
<input type="checkbox"/> Gel		
10 x 10 - 20 x 18	E2605	\$480
20 x 20 - 20 x 20	E2605	\$530
22 x 16 - 22 x 20	E2606	\$530
<input type="checkbox"/> Visco		
10 x 10 - 20 x 18	E2607	\$480
20 x 20 - 20 x 20	E2607	\$530
22 x 16 - 22 x 20	E2608	\$530
<input type="checkbox"/> Structure Cushion		
<input type="checkbox"/> Gel <input type="checkbox"/> Visco		
10 x 10 - 20 x 18	E2607	\$560
20 x 20 - 20 x 20	E2607	\$610
22 x 16 - 22 x 20	E2608	\$610
<input type="checkbox"/> TRU-Comfort Cushion		
16 x 16 - 20 x 18		\$380
20 x 20 - 20 x 22		\$380
22 x 16 - 22 x 24		\$480

SPECIALTY SEATING

Specialty Seat Kit (limited adjustability) NC

SYSTEM SPECIFICATIONS

Specialty Seating Kit Width

- 16" NC
- 18" Std.
- 20" NC

Specialty Seating Kit Depth

- 16" Std. 17" NC
- 18" NC 19" NC
- 20" NC

Specialty Seating Kit Back Cane Height

- 18" Std. 20" NC
- 22" NC

Back Selection

- See Synergy Shape Back Order Form^{11,12,13}
- Solid Plate Back up to 20"W¹⁷ \$275
- See Positioning Components Order Form for additional backs

Back Height Solid Pan¹⁴

- 18" Std. 19" NC
- 20" NC 21" NC
- 22" NC 23" NC
- 24" NC 25" NC

Sling Back

- 16" wide \$98
- 18" wide \$98
- 20" wide \$98

Sling Back Height

(Hgt must match cane hgt.)

- 18" height NC
- 20" height NC
- 22" height NC

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

ARMRESTS

- 2-Post, flip-back^{13,16} Std.
- Flip-back, cane mounted^{13,16} \$95ea
- Removable, adjustable, quick height^{15,18} E0973 \$275ea
- Removable, height adjustable¹⁸ \$150ea
- heavy duty armrests E0973
- 8"H-12"H 10"H-14"H
- Pediatric, removable,¹⁸ adjustable, quick height E0973 \$250ea

Armrest Pads

- Full Length, straight, 2"x14" NC
- Left Right
- Desk Length, straight, 2"x10" NC
- Left Right
- Full Length, Waterfall, 2"x14"⁹ NC
- Left Right
- Desk Length, Waterfall, 2"x10"^{9,10} NC
- Left Right

HEADRESTS³

- Stealth^{7,8}** E0955
 - Large (11" wide x 6" high) \$210
 - Stealth^{7,8}** E0955
 - Small (9 3/4" wide x 3" high) \$210
 - Stealth**, swing-away (Flip-Down)
 - Head Rest Bracket E1028 \$150
- Interfaces w/ any headrest that uses 1/2" square down-tube.
- Medical Applications^{7,8}** E0955
 - Large (11" wide x 6" high) \$275
 - Medical Applications^{7,8}** E0955
 - Small (9 3/4" wide x 3" high) \$275
 - Fixed, Offset Mounting Bracket NC

Footnotes for pages 3:

1. Includes black vinyl seat pan, back available in all listed colors.
2. 18" x 18" std. Special measurements are not available.
3. Only available w/ solid backs.
4. N/A on manual recliner.
5. Rear mount accessories are N/A.
6. See Synergy Cushion order forms for detailed pricing, sizes & codes.
7. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
8. Comes with swing-away / removable hardware; HCPCS code is E1028.
9. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
10. S/A joystick mount N/A.
11. 14" & 22" width are N/A.
12. Micro-drive, Shape Back & high-mount foot platform are N/A in the same configuration.
13. Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
14. Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
15. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
16. S/A laterals interfere w/ flip-back arms.
17. 20"D is N/A.
18. Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.

Quantum® 6000 Series Group 3

Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

Synergy Seating Systems¹¹

Synergy Seating™ System \$434
(Price with Base)

Synergy Manual Recline Pkg.^{12,13,14,15,22} \$700
Maximum weight 300 lbs. E1226

Synergy Manual Tilt Pkg.^{13,16} \$550
(Max Dimensions 18"W x 18"D) K0108
Maximum weight 160 lbs.

Synergy Manual Systems come standard with mercury drive lock-out and are in addition to the MSRP of the Synergy Seat.

WIDTH

Set Seat Width At:

<input type="checkbox"/> 10"	NC	<input type="checkbox"/> 11"	NC
<input type="checkbox"/> 12"	NC	<input type="checkbox"/> 13"	NC
<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC		

Seat Width Package

<input type="checkbox"/> 10"W-14"W	<input type="checkbox"/> 14"W-18"W
<input type="checkbox"/> 16"W-20"W	<input type="checkbox"/> 20"W-22"W

Seat Width Package Manual Tilt

<input type="checkbox"/> 10"W-12"W	<input type="checkbox"/> 12"W-14"W
<input type="checkbox"/> 14"W-16"W	<input type="checkbox"/> 16"W-18"W

DEPTH

Set Seat Depth At:

<input type="checkbox"/> 10"	NC	<input type="checkbox"/> 11"	NC
<input type="checkbox"/> 12"	NC	<input type="checkbox"/> 13"	NC
<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	Std.	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	NC	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC		

Seat Depth Package

<input type="checkbox"/> 10"D-14"D	<input type="checkbox"/> 14"D-18"D
<input type="checkbox"/> 16"D-20"D	

Seat Depth Package Manual Tilt

<input type="checkbox"/> 10"D-14"D	<input type="checkbox"/> 14"D-16"D
<input type="checkbox"/> 16"D-18"D	

BACK CANES

Synergy Manual Recline comes standard with 22" non-height adjustable back canes.

Back Cane Height

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC	<input type="checkbox"/> 25"	NC

BACK OPTIONS

See Synergy Shape Back Order Form^{1,2,3}

Solid Plate Back up to 20" W²² \$275

Solid Plate Back 21-24" W²² \$325

Back Height:⁴

18"H 19"H 20"H 21"H

22"H 23"H 24"H 25"H

Sling Back²⁰ \$98

Back height will be the same as back can height.
Includes gap for cushion.

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

TRU-Comfort Back¹²

Back Height: 22" 25"

Back up to 20" W \$595

22" W \$695

ARMRESTS

Flip-back, cane mounted^{3,5,18} Std.

2-Post, flip-back, height^{3,5,19} E0973 \$150ea
adjustable, removable armrests

Removable, adjustable^{6,21} E0973 \$75ea
quick height

Quad Release Lever NC
(3" of height adjustment in 3/8" increments)

Removable, height adjustable²¹ E0973 \$75ea
heavy duty armrests

8"H-12"H 10"H-14"H
 Pediatric removable²¹ E0973 \$75ea
adjustable, quick height.

Match Shroud Color Black

Quad Release Lever NC
(7 3/4" - 10 3/4" height adjustment in 1/2" increments from seat pan to top of armrest)

ARMPADS

Full length, Straight, 2" X 14" Std.

Left Right

Desk Length, Straight, 2" X 10" Std.

Left Right

Full Length, Waterfall, 2" X 14"⁷ NC

Left Right

Desk Length, Waterfall, 2" X 10"^{7,17} NC

Left Right

HEAD REST⁸

Stealth^{9,10} E0955
Large (11" wide x 6" high) \$210

Stealth^{9,10} E0955
Small (9 3/4" wide x 3" high) \$210

Stealth, swing-away (Flip-Down)
Head Rest Bracket E1028 \$150

Interfaces w/ any headrest that uses 1/2" square down-tube.

Medical Applications^{9,10} E0955
Large (11" wide x 6" high) \$275

Medical Applications^{9,10} E0955
Small (9 3/4" wide x 3" high) \$275

Fixed, Offset Mounting Bracket NC

TRU-Comfort^{9,10} E0955
 8" Wide w/ Cobra hardware \$417

10" Wide w/ Cobra hardware \$417

Bracket for after market headrest \$180

Footnotes for pages 4:

- 14" & 22" width are N/A.
- Micro-drive, Shape Back & high-mount foot platform are N/A in the same configuration.
- Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- S/A laterals interfere w/ flip-back arms.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- Available w/ solid backs only.
- Comes with swing-away / removable hardware; HCPCS code is E1028.
- These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
- Manual systems do not include a headrest.
- N/A w/ flip-back, cane mounted arms.
- Solid back not included.
- Includes removable, adjustable, quick-height armrests.
- N/A w/ Shape Back.
- Manual Tilt w/ Shape back & oxygen tank holder requires the use of a pushbar.
- N/A w/ S/A joystick mount.
- N/A w/ Synergy manual recline.
- May cause interference with certain dimension configurations. Call Quantum Customer Service.
- N/A w/ Synergy manual systems.
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- Solid plate back is N/A w/ manual recline.

Quantum® 6000 Series Group 3

Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

RECOMMENDED INFORMATION

Height: _____
Weight: _____
Leg Length: _____
Lower Leg Length: _____

TRU-BALANCE 2 TILT

TRU-Balance 2 Power Tilt E1002 \$4995
55 degree center of gravity (CG) power tilt seating system. 300lbs weight capacity.

TRU-Balance Super Low Tilt \$4995
55 degree center of gravity (CG) power tilt seating system. 300lbs weight capacity.
Utilizes TB 1 Tilt, Requires Group 34 batteries
Fixed Towers STF: 16.625"
Ang-Adj. Towers STF: 17.375"

TRU-BALANCE 2 RECLINE

TRU-Balance Power Recline E1004 \$5595
**Full Drive lock-out at 20 degrees
300 lbs. weight capacity**

TRU-Balance \$1895
Power-Elevating Seating Option¹
**7" of Lift & 1/4 Speed Inhibit- 300 lbs. max. wt.
Remote Plus operation only.**

TRU-BALANCE 2 TILT & RECLINE

TRU-Balance 2 Power Tilt & Recline E1007 \$9595
**Full Drive lock-out at 20 degrees
300 lbs. weight capacity
55 degree center of gravity
(CG) power tilt seating system.**

TRU-Balance Super Low Tilt & Recline \$9595
55 degree center of gravity (CG) power tilt seating system. 300lbs weight capacity.
Utilizes TB 1 Tilt, Requires Group 34 batteries
Fixed Towers STF: 16.625"
Ang-Adj. Towers STF: 17.375"

TRU-BALANCE PES

TRU-Balance \$1895
Power-Elevating Seating Option¹
**7" of Lift & 1/4 Speed Inhibit- 300 lbs. max. wt.
Remote Plus operation only.**

SYSTEM SPECIFICATIONS

Tower Selection(Tilt or Tilt/Recline only)

- Adjustable-angle Std.
Will raise minimum STF by 7/8".
 Fixed NC

Set Seat Width At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	

Seat Width Package

- 14"W-18"W 16"W-20"W
 20"W-22"W

Set Seat Depth At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	

Seat Depth Package

- 14"D-18"D 16"D-20"D

*Push Bar selection must match back width;
center tube standard when selecting flip-up arms.*

Set Back Cane Height At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

Seat Overall Back Height³

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

BACKS

- See Synergy Shape Back Order Form⁶
 Solid Plate Back⁴ Std.
 Attendant manual recline^{4,24,25,26,27} \$750
180 degrees of recline with 2 1/2" raised pivot point. Controlled through dual gas struts. 300 lb. weight capacity.
(For use on TB Tilt Systems only)
 TRU-Balance Reclining Back Std.
(For use on Recline Systems only)
 TRU-Comfort Back⁴
Back Height: 22" 25"
 Back up to 20" W \$595
 22" W \$695
 Omit Back -(\$100)

ARMRESTS

- Flip-back, cane mounted^{5,9,14} Std.
Mounted on back canes with full width spreader bar.
 2-Post, flip-back, hgt. adj, removable^{9,15} NC
 Removable, adj. quick height^{15,16} NC
 Quad Release Lever NC
(3" of height adjustment in 3/8" increments)
 Removable, single-post,¹⁶
height adjustable, heavy duty NC
 8"-12" Adjustment
 10"-14" Adjustment
 Pediatric, removable,¹⁶ NC
adjustable, quick height
 Quad Release Lever NC
 Match Shroud Kid Color Black
 Reclining, Flip-back⁹ Std.
 8"-9" H 10"-11" H

(For use on Recline Systems only)

Armrest Pads

- Full Length, straight, 2"x14" NC
 Left Right
 Desk Length, straight, 2"x10" NC
 Left Right
 Full Length, Waterfall, 2"x14"¹⁷ NC
 Left Right
 Desk Length, Waterfall, 2"x10"^{17,19} NC
 Left Right

HEAD REST

- Stealth^{18,19} E0955
Large (11" wide x 6" high) \$210
 Stealth^{18,19} E0955
Small (9 3/4" wide x 3" high) \$210
 Stealth, swing-away (Flip-Down)
Head Rest Bracket E1028 \$150
Interfaces w/ any headrest that uses 1/2" square down-tube.
 Medical Applications^{18,19}
Large (11" wide x 6" high) E0955 \$275
 Medical Applications^{18,19}
Small (9 3/4" wide x 3" high) E0955 \$275
 Omit Head Rest - To be supplied by provider
 Fixed, Offset Mounting Bracket NC
 Reclining Head Rest (6.5" x 7")^{18,19} E0955 \$180
(For use on Recline Systems only)
 TRU-Comfort^{18,19} E0955
 8" Wide w/ Cobra hardware \$417
 10" Wide w/ Cobra hardware \$417
 Bracket for after market headrest \$180

Quantum® 6000 Series Group 3

Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

TILT POWER POSITIONING ELECTRONICS

Select Electronics to be used

PG Remote Plus

Q-Logic (TB Tilt only)

Electronics for Tilt Only

Through Toggle Box²⁰ NC

Left Right (mounts to side of armrest)

Through Single Switch Tilt Module^{21,22} E2310 \$500

Tilt or Recline Only Power Positioning

Single Motor Controller:

Through Remote Plus Controller E2310 \$1175

Through Q-Logic Controller E2310 \$1175

Through 2-Directional Control Box E2310
w/ dual Switch Jacks^{21,22} \$1175

Multiple Motor Controller:²³

Through Remote Plus Controller E2311 \$2375

Through Toggle¹⁴ (Select one) E2311 \$2375

4-way Toggle

4-way Button Toggle

Quad Push-Button Control

Toggle box mounts to side of armrest.

Left Right

TILT/RECLINE POWER POSITIONING

All options require motor control module.

PG Electronics- Multiple Motor controller

Single Function:

Operation through PG E2310 \$1175

Operation through toggle NC
(Select one)

4-way Toggle

4-way Button Toggle

Quad Push-Button Control

Toggle box mounts to side of armrest.

Left Right

PG Electronics- Multiple Motor controller

Multiple Function:

Operation through PG E2311 \$2375

Operation through toggle E2311 \$2375
(Select one)

4-way Toggle

4-way Button Toggle

Quad Push-Button Control

Toggle box mounts to side of armrest.

Left Right

ACCESSORIES

Rear Basket^{6,7,12} \$103

Cane/Crutch Holder^{6,7,12} \$98

Oxygen Holder^{6,7,12,13} E2208 \$103

Oxygen Holder^{6,7,10,13} E2208 \$150

(TRU-Balance Tilt & Synergy Seat only)

Stroller handle \$85

(Synergy Manual Tilt only)

Accessory Mounting Bracket \$50

(1" adaptor tube for use w/ TB Power Recline)

Walker Holder^{2,6,7,12} \$92

Flag Holder (SOF on rehab seating) \$17

Weather Cover, Large, Black \$202

Saddle bag(contoured seating only) \$28

Cup holder \$17

Left Right

N/A on desk length arm with S/A
mount or toggle box.

Lap Belt Std.

Pediatric 38" N/A w/ limited recline seating

50"⁸ 60" 70" 80" 90"

Vent battery storage kit¹¹ \$250

(Includes extended battery cover, strap &
bracket, for use with third battery in base)

Can accommodate up to NF-22 battery-not included

N/A w/ footplatform, high mount footplatform
or super low system

Footnotes for pages 5 & 6:

- STF will increase.
- S/A joystick mount N/A.
- Measured from seat pan to top of back pan w/ 3-4" gap on tilt & 7-8" gap on recline to accommodate the cushion.
- N/A w/ flip-back, cane mounted arms.
- Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
- Includes mounting bracket and hardware.
- N/A on reclining seats.
- N/A on manual recline high-back.
- S/A laterals interfere w/ flip-back arms.
- Manual Tilt w/ Shape Back and oxygen holder requires the use of a pushbar.
- Vent Tray is SOF due to vent and seating type.
- N/A w/ tilt systems.
- N/A w/ flip-back arms.
- N/A w/ Attendant manual recline.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- Comes with swing-away / removable hardware; HCPCS code is E1028.
- These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
- N/A w/ waterfall arm pads.
- Consult specialty control order form for a wide variety of switch types and sizes-- 2 required.
- Does not include switch.
- Multiple motor controller options are only available w/ a PES.
- N/A w/ Shape Back.
- Includes price for removable, adjustable, quick-height arms & 2-post arms.
- Only available on TB Tilt Systems.
- Solid plate back is N/A w/manual recline.

Lowest seat-to-floor height

Specialty Seat:	16.875"
Synergy Seat:	16.625"
TB 2 Tilt:	18.125"
(w/ Ang-Adj towers)	
TB 2 Tilt:	17.375"
(w/ fixed towers towers)	
TB 2 Recline:	16.625"
TB 2 Tilt/Recline:	18.125"
(w/ Ang-Adj towers)	
TB 2 Tilt/Recline:	17.375"
(w/ fixed towers towers)	
High-back Seat:	21.5"

Quantum Custom Delivery Service

Includes: liftgate service and delivery of your product(s) inside your facility
Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.

All measurements +/- 1/8" with 18" x 18" seat, fixed towers (please see Tilt information above), 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Increased widths may increase seat-to-floor height. Contact customer service for details. Listed measurements are in the lowest position. It does not include elevating systems.