

# QUANTUM® 4FRONT<sup>2</sup>



## 4FRONT<sup>2</sup>

shown in Military Green with  
TRU Balance<sup>®</sup> 4 power positioning,  
Q-Logic 3 environmental  
controls and Stealth<sup>®</sup>  
TRU Comfort 2 cushion



FDA Class II  
Medical Device

The 4Front<sup>®</sup> 2 is the next generation of the 4Front<sup>®</sup> Power Chair! New Smart Traction Control (STC), along with Smooth Ride Suspension (SRS), provide stability for smooth obstacle transitions and unmatched driving performance. Additional new features include 4-pole motors and redesigned caster arms for improved durability and aesthetics and reduced weight. TRU Balance<sup>®</sup> 4 Power Positioning with anterior tilt and memory seating.

## SMART TRACTION CONTROL



| [quantumrehab.com](http://quantumrehab.com)

## Specifications

<b>DRIVE WHEELS</b>	Standard: 14" flat free Optional: 14" knobby 14" pneumatic
<b>CASTER WHEELS</b>	Rear: 9" flat free
<b>ANTI-TIP WHEELS</b>	Front: 4" flat free
<b>MAXIMUM SPEED <sup>1</sup></b>	Up to 6.25 mph (3.5 mph when elevated)
<b>GROUND CLEARANCE</b>	3.14" (battery tray)
<b>TURNING RADIUS <sup>2</sup></b>	25.9"
<b>OVERALL LENGTH <sup>2</sup></b>	39.9"
<b>BASE WIDTH</b>	24"
<b>SEATING SIZES</b>	
<b>TRU Balance® 4 Anterior Tilt</b>	Width: 16-22" Depth: 14-22"
<b>TRU Balance® 3 Static Seat</b>	Width: 12-22" Depth: 12-22"
<b>TRU Balance® 3 Power Tilt TRU Balance® 4 Power Tilt</b>	Width: 12-22" Depth: 12-22"
<b>TRU Balance® 3 Power Lift &amp; Tilt</b>	Width: 12-22" Depth: 12-22"
<b>TRU Balance® Flex Seating</b>	Width: 15-22" Depth: 15-20"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>TRU Balance® 4</b>	17.75-18.75"
<b>TRU Balance® 3 Static Seat <sup>5</sup></b>	15-17" (not available with Synergy®)
<b>TRU Balance® 3 Power Tilt <sup>5</sup> TRU Balance® 4 Power Tilt <sup>5</sup></b>	17.75-19.75"
<b>TRU Balance® 3 Power Lift &amp; Tilt <sup>5</sup></b>	17.75-19.75"
<b>TRU Balance® Flex Seating</b>	16.5-18.5"
<b>BATTERY SIZE <sup>3</sup></b>	
<b>BATTERY WEIGHT <sup>7</sup></b>	Group 24: 51.5 lbs. each Group 34: 44.8 lbs. each
<b>AVAILABLE ELECTRONICS <sup>8</sup></b>	120A Q-Logic 3 75A Q-Logic NES
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	4-pole
<b>WEIGHT CAPACITY <sup>4</sup></b>	300 lbs.
<b>BASE WEIGHT</b>	130 lbs.
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard (K0848/K0849) Single Power: (K0856/K0857) Multiple Power: (K0861)

- Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- Without front riggings
- Battery capacity may change when adding power positioning systems.
- Weight capacity is dictated by the selection of seating.
- All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- Battery weight may vary +/- 3 lbs. based on manufacturer.
- All Quantum Rehab® units can be ordered with specialty controls.
- Call for availability
- Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front 2's suspension system, seat-to-floor heights are approximately 0.25" lower when the unit is loaded.

Quantum Rehab® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

## Standard Features

- 4-pole motors
- Smart Traction Control (STC)
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions
- Redesigned caster arms for improved durability and aesthetics and reduced weight
- Q-Logic 3e controller
- Increased speeds up to 6.25 mph
- LED front and rear fender lights
- Quiet ride technology

## Options

- 10" of power adjustable safe seat elevation at speeds of up to 3.5 mph
- Depth-adjustable vent tray<sup>8</sup>
- Compact vent tray<sup>8</sup>
- Quantum® backup camera

## Colors



## Hydrographics

