

QUANTUM® 4FRONT²



4FRONT²
shown in Military Green with
TRU Balance[®] 4 power positioning,
Q-Logic 3 environmental
controls and Stealth[®]
TRU Comfort 2
cushion.

AVAILABLE WITH 3.5 MPH
10"
SAFE SEAT ELEVATION



The 4Front[®] 2 is the next generation of the 4Front[®] Power Chair! New Smart Traction Control (STC), along with Smooth Ride Suspension (SRS), provide stability for smooth obstacle transitions and unmatched driving performance. Additional new features include 4-pole motors and redesigned caster arms for improved durability and aesthetics and reduced weight. TRU Balance[®] 4 Power Positioning with anterior tilt and memory seating.

SMART TRACTION CONTROL



QUANTUM[®]
#1 FOR REHAB POWER



quantumrehab.com

Specifications

DRIVE WHEELS	Standard: 14" flat free Optional: 14" knobby 14" pneumatic
CASTER WHEELS	Rear: 9" flat free
ANTI-TIP WHEELS	Front: 4" flat free
MAXIMUM SPEED ¹	Up to 6.25 mph (3.5 mph when elevated)
GROUND CLEARANCE	3.14" (battery tray)
TURNING RADIUS ²	25.9"
OVERALL LENGTH ²	39.9"
BASE WIDTH	24"
SEATING SIZES	
TRU Balance® 4 Anterior Tilt	Width: 16-22" Depth: 14-22"
TRU Balance® 3 Static Seat	Width: 12-22" Depth: 12-22"
TRU Balance® 3 Power Tilt TRU Balance® 4 Power Tilt	Width: 12-22" Depth: 12-22"
TRU Balance® 3 Power Lift & Tilt	Width: 12-22" Depth: 12-22"
TRU Balance® Flex Seating	Width: 15-22" Depth: 15-20"
SEAT-TO-FLOOR HEIGHTS	
TRU Balance® 4	17.75-18.75"
TRU Balance® 3 Static Seat ⁵	15-17" (not available with Synergy®)
TRU Balance® 3 Power Tilt ⁵ TRU Balance® 4 Power Tilt ⁵	17.75-19.75"
TRU Balance® 3 Power Lift & Tilt ⁵	17.75-19.75"
TRU Balance® Flex Seating	16.5-18.5"
BATTERY SIZE ³	
BATTERY WEIGHT ⁷	Group 24 (2) or Group 34 (2)
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 75A Q-Logic NES
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-pole
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	130 lbs.
MEDICARE CODE COVERAGE	Group 3: Standard (K0848/K0849) Single Power: (K0856/K0857) Multiple Power: (K0861)

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 7) All Quantum Rehab® units can be ordered with specialty controls.
- 8) Call for availability
- 9) Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front 2's suspension system, seat-to-floor heights are approximately 0.25" lower when the unit is loaded.

Standard Features

- 4-pole motors
- Smart Traction Control (STC)
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions
- Redesigned caster arms for improved durability and aesthetics and reduced weight
- Q-Logic 3e controller
- Increased speeds up to 6.25 mph
- LED front and rear fender lights
- Quiet ride technology


Options

- 10" of power adjustable safe seat elevation at speeds of up to 3.5 mph
- Depth-adjustable vent tray⁸
- Compact vent tray⁸
- Quantum® backup camera

Colors



Hydrographics



QUANTUM®

#1 FOR REHAB POWER

(US) 866-800-2002

| (Canada) 888-570-1113

| quantumrehab.com

