



12" x 12"
seat sizes now
available with
iLevel®
-SOF ONLY-



Q6 edge[®] 2.0

With features like standard 6 mph motors and advanced responsiveness, the Q6 Edge[®] 2.0 is engineered to meet the performance needs of the most active user. The Q6 Edge 2.0 accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel[®] technology.

quantumrehab.com

QA-Q6EDGE2-022316



Beamsville, Ontario

888-570-1113

The Rehab Power Chair
Solutions Company

Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ⁴	Up to 6 mph
GROUND CLEARANCE	2.7" (battery tray)
TURNING RADIUS ⁵	20.5"
OVERALL LENGTH	35.5" without foot riggings
BASE WIDTH	24.4"
SEATING SIZES	
Synergy®/TRU-Balance® 3 Static Seat	Width: 10-22" Depth: 10-22"
TRU-Balance® 3 Power Tilt	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Lift & Tilt	Width: 16-22" Depth: 14-22"
SEAT-TO-FLOOR HEIGHTS	
Synergy®/TRU-Balance® 3 Static Seat ⁸	16.5-18.5"
TRU-Balance® 3 Power Tilt ⁸	17.5-19.5"
TRU-Balance® 3 Lift & Tilt ⁸	17.5-19.5"
RECLINE ¹⁰	Yes
BATTERY SIZE ⁶	NF-22
BATTERY WEIGHT ¹⁴	38 lbs.
AVAILABLE ELECTRONICS ¹²	75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole Mid-Wheel 6 [®]
WEIGHT CAPACITY ⁷	300 lbs.
BASE WEIGHT	138.6 lbs.

Colors

SEAT BACKS for units without TRU-Balance® 3 are only available in *Back in Black*.




Features

- Standard 6 mph motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics
- iLevel® option available with up to 10" of seat elevation (in 16 seconds) at walking speed (up to 3.5 mph)
- iLevel is retrofittable on the Q6 Edge® 2.0
- Compact mid-wheel drive base for enhanced maneuverability and access
- ATX Suspension for optimum ride and enhanced obstacle climbing
- Wide range of standard and rehab seating systems to meet any user's needs
- Low impact OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups

Options

- TRU-Balance® 3 10" Power Adjustable Seat with iLevel® option allows 10" of seat elevation in 16 seconds while at walking speed (3.5 mph)
- TRU-Balance® 3 Seating System
- TRU-Comfort 2 Seating
- Q-Logic 2 EX Drive Controls
- Trilogy vent tray⁹
- Compact vent tray⁹
- LED full lighting package
- 4" wide flat-free knobby extended traction tires⁹
- New LED fender lights
- 12" x 12" seat sizes with iLevel (SOF only)
- Wide range of specialty seating options (See Stealth Products order forms)
- Wide range of specialty drive controls (See the Q-Logic controls and Stealth Products order forms)
- Wide range of optional accessories (See order form)



NEW LED FENDER LIGHTS

- Standard
- Optional
- Requires HD Suspension
- Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- Without front riggings
- Battery capacity may change when adding power positioning systems.
- Weight capacity is dictated by the selection of seating.
- All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- Call for availability.
- Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- Requires Group 34 batteries
- All Quantum Rehab® units can be ordered with specialty controls
- Required with high-speed motors
- Battery weight may vary +/- 3 lbs. based on manufacturer.



Beamsville, Ontario • 888-570-1113 • www.quantumrehab.com • www.ilevel.rehab

