

# EDGE 3

EDGE 3 shown in Blue Jean  
with iLevel® power seat elevation  
and Q-Logic 3 with Bluetooth®

AVAILABLE WITH  
12" iLEVEL® POWER  
SEAT ELEVATION AT  
**4.5 MPH**



**PROUDLY MADE IN**  
**AMERICA**   
OF FOREIGN AND DOMESTIC PARTS



[quantumrehab.com](http://quantumrehab.com)

  
**QUANTUM®**



## Specifications

<b>DRIVE WHEELS</b>	Standard: 14" Solid Optional: 14" Pneumatic
<b>CASTER WHEELS</b>	Front: 6" Solid Rear: 6" Solid
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED <sup>1</sup></b>	Up to 6 mph Up to 4.5 mph w/iLevel activated
<b>GROUND CLEARANCE</b>	2.7" (battery tray)
<b>TURNING RADIUS <sup>2</sup></b>	20.5"
<b>OVERALL LENGTH <sup>2</sup></b>	32.23" without front riggings
<b>BASE WIDTH</b>	24.4"
<b>SEATING SIZES</b>	
<b>Synergy®/TRU-Balance® 3 Static Seat</b>	Width: 10-22" Depth: 10-22"
<b>TRU-Balance® 3 Power Tilt</b>	Width: 12-22" Depth: 12-22"
<b>TRU-Balance® 3 Tilt and iLevel®</b>	Width: 12-22" Depth: 12-22"
<b>TRU-Balance® Flex Seating</b>	Width: 15-22" Depth: 12-20"
<b>SEAT-TO-FLOOR HEIGHTS <sup>10</sup></b>	
<b>Synergy®/TRU-Balance® 3 Static Seat <sup>5</sup></b>	16.5-18.5"
<b>TRU-Balance® 3 Power Tilt <sup>5</sup></b>	17.5-19.5"
<b>TRU-Balance® 3 Tilt &amp; iLevel® <sup>5</sup></b>	17.5-19.5"
<b>TRU-Balance® Flex Seating</b>	16.5" - 19.5"
<b>MANUAL RECLINE <sup>6</sup></b>	Yes
<b>BATTERY SIZE <sup>3</sup></b>	NF-22 (2)
<b>BATTERY WEIGHT <sup>7</sup></b>	38 lbs. each
<b>AVAILABLE ELECTRONICS <sup>8</sup></b>	120A Q-Logic 3 EX 75A Q-Logic 2 NE 75A Q-Logic 2 NE+
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-Pole, Mid-Wheel 6®
<b>WEIGHT CAPACITY <sup>4</sup></b>	300 lbs.
<b>BASE WEIGHT</b>	138.6 lbs.

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.

2) Without front riggings

3) Battery capacity may change when adding power positioning systems.

4) Weight capacity is dictated by the selection of seating.

5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.

6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.

7) Battery weight may vary +/- 3 lbs. based on manufacturer.


8) All Quantum Rehab® units can be ordered with specialty controls.

9) Call for availability

10) 16.5" seat-to-floor height is not available with power tilt or iLevel®. 19.5" seat-to-floor height is not available with static seat or power recline.

## Features


- Upgraded SRS (Smooth Ride Suspension)
- 4.5 mph with iLevel®
- Fender lights and USB charger included
- Bold new colors
- Made in the USA of foreign and domestic parts
- XLR USB Charger



## 7 New Colors\*



BLUE  
JEAN



CHERRY  
BOMB



LEMON  
CRUSH



GREEN  
MACHINE



RASPBERRY  
BERET



EVOLUTION  
ORANGE



SUGAR  
PLUM

## HYDROGRAPHICS



CARBON  
FIBER



DIGITAL  
CAMO

\* Visit [quantumrehab.com](http://quantumrehab.com) to see all 19 available colors.