



Q610

The compact Q610 is the ultimate choice for clients who demand the superb stability of Mid-Wheel 6® Design and tight-quarter maneuverability. Highly adaptable, the Q610 accepts a wide range of seating and electronics options while its powerful drivetrain delivers high-performance torque.

quantumrehab.com

CA-Q610-01116

Beamsville, Ontario

Q
QUANTUM®

888-570-1113

The Rehab Power Chair
Solutions Company

Specifications

DRIVE WHEELS	10" solid (standard) or 10" pneumatic (optional)
OMNI-CASTERS	Front: 5" solid Rear: 6" solid
SUSPENSION	ATX Suspension (Active-Trac® with extra stability)
MAXIMUM SPEED ¹	Up to 4 mph
BRAKES	Intelligent Braking (electronic, regenerative disc brakes)
GROUND CLEARANCE	2.5" (frame)
TURNING RADIUS ²	22"
OVERALL SIZE ²	Length: 34.74" (without foot riggings) Width: 23.19"
SEATING	Synergy® Seat: Adjustable Width: 10"-20" Adjustable Depth: 10"-20" TRU-Balance® 2 Power Positioning Systems: Adjustable Width: 14"-20" Adjustable Depth: 14"-20" TRU-Balance® 2 Lift & Tilt: Adjustable Width: 16"-20" Adjustable Depth: 16"-20" High-back Seat (20" x 20" max. dim.) Medium-back Seat (20" x 20" max. dim.)
DRIVETRAIN	2-motor, Mid-Wheel 6®
BATTERIES	Two 12 volt, deep cycle U-1
RANGE PER CHARGE ¹	Up to 20 miles
BATTERY CHARGER	5A, off-board
AVAILABLE ELECTRONICS	70A PG VR2 70A Q-Logic 2 EX
WEIGHT CAPACITY	300 lbs. maximum (user weight) ³
BASE WEIGHT	100.2 lbs.
BATTERY WEIGHT	23.4 lbs. (each)
WARRANTY	Lifetime limited warranty on frame 2-year limited warranty on electronics and drive motor

Features

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® 2 Power Positioning
- Side-mounted, easily accessible freewheel levers
- ATX Suspension (Active-Trac® with extra stability) for enhanced performance over more varied terrain
- Easy front access to batteries
- Powerful GLIDE motors for added power and torque

Options

- TRU-Balance® 2 Lift & Tilt
- TRU-Comfort Seating System
- Q-Logic 2 EX Drive Controls
- LED full lighting package
- Power elevating seat
- Compact vent tray

Colors



FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
- 3) 250 lbs. with TRU-Balance® Power Positioning Systems



Beamsville, Ontario • 888-570-1113 • www.quantumrehab.com • www.ilevel.rehab

