



# RIVAL®

The Rival® rear-wheel drive power chair features innovative design that provides exceptional outdoor performance and handling while delivering outstanding tight space maneuverability.

The Rival accepts a complete range of rehab seating and electronics options.

[quantumrehab.com](http://quantumrehab.com)

CA-RIVAL-020516

**Q**  
QUANTUM®

Beamsville, Ontario

888-570-1113

The Rehab Power Chair  
Solutions Company

# Specifications

<b>DRIVE WHEELS</b>	<b>Standard:</b> 14" solid <b>Optional:</b> pneumatic
<b>CASTER WHEELS</b>	<b>Front:</b> 9" pneumatic
<b>SUSPENSION</b>	Active-Trac® ATX with Long Trail Arm
<b>MAXIMUM SPEED</b>	Up to 12 km/h (7.5 mph) <sup>6</sup>
<b>BRAKES<sup>1</sup></b>	Intelligent braking (electronic, regenerative disc brakes)
<b>GROUND CLEARANCE</b>	3.25" (frame base)
<b>TURNING RADIUS<sup>2</sup></b>	35.75"
<b>OVERALL SIZE<sup>2</sup></b>	Length: 36" without foot riggings Width: 25"
<b>SEATING</b>	<b>Synergy® / TRU-Balance® 3 Static Seat:</b> Adjustable Width: 12"-22" Adjustable Depth: 12"-22" <b>TRU-Balance® 3 Power Positioning Systems:</b> Adjustable Width: 12"-22" Adjustable Depth: 12"-22" <b>High-back Seat</b> (22" x 22" max. dim.) <b>Medium-back Seat</b> (22" x 22" max. dim.)
<b>DRIVETRAIN</b>	2-motor, rear-wheel
<b>BATTERIES</b>	Two 12 volt, deep cycle, Group 24
<b>RANGE PER CHARGE<sup>1,3</sup></b>	Up to 19.5 miles
<b>BATTERY CHARGER</b>	8A, off-board
<b>AVAILABLE ELECTRONICS</b>	75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX
<b>WEIGHT CAPACITY</b>	300 lbs.
<b>BASE WEIGHT</b>	129 lbs.
<b>BATTERY WEIGHT<sup>4</sup></b>	52 lbs.
<b>WARRANTY</b>	Lifetime limited warranty on frame; 2-year warranty on electronics; 18-month warranty on drive motor

### FOOTNOTES:

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
- 3) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4, standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, -5%.
- 4) Battery weight may vary based on manufacturer.
- 5) Optional.
- 6) With optional high-speed motor package

# Features

- Compact size for outstanding maneuverability in tight spaces
- Maximum speed up to 7.5 mph<sup>6</sup>
- High-performance, high-efficiency motor package
- Accepts a wide range of seating and positioning systems
- Sliding battery tray for easy access
- Two suspension packages

# Options

- TRU-Balance® 3 Power Tilt
- TRU-Balance® 3 Power Recline
- TRU-Balance® 3 (10") Lift
- Manual and Power Elevating Leg Rests
- TRU-Comfort 2 Seating
- Q-Logic 2 Controller
- 9" solid casters
- 10" pneumatic casters
- Full Lighting Package
- Centermount Adjustable Foot Platform
- Accu-Trac Motors
- Power Articulating Foot Platform
- Fender light package (coming soon)
- 8" solid casters (coming soon)

# Colors

SEAT BACKS for units without TRU-Balance® 3 are only available in *Back in Black*.



Beamsville, Ontario • 888-570-1113 • [www.quantumrehab.com](http://www.quantumrehab.com) • [www.ilevel.rehab](http://www.ilevel.rehab)

