

# EDGE® 3

EDGE® 3 shown in Blue Jean with captain's seating

**NOW AVAILABLE**

**MAXIMUM SPEED UP TO 7 MPH**  
with new optional 4-pole motors



FDA Class II  
Medical Device\*

Upgraded SRS (Smooth Ride Suspension) provides a smooth, comfortable ride and advanced stability. Optional iLevel® is available at 4.5 mph with up to 12" of power adjustable seat height. Standard features include front and rear LED fender lights, a USB charger, and drive wheel color accents.



**BLUE JEAN**



**QUANTUM®**  
#1 FOR REHAB POWER



[quantumrehab.com](http://quantumrehab.com)

## Specifications

<b>DRIVE WHEELS</b>	<b>Standard:</b> 14" Solid <b>Optional:</b> 14" Pneumatic; 14" Knobby
<b>CASTER WHEELS</b>	<b>Front:</b> 6" Solid <b>Rear:</b> 6" Solid
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 7 mph <sup>6</sup> (4.5 mph with iLevel®)
<b>GROUND CLEARANCE</b>	2.7" (battery tray)
<b>TURNING RADIUS</b> <sup>2</sup>	20.5"
<b>OVERALL LENGTH</b> <sup>2</sup>	35.5"
<b>BASE WIDTH</b>	24.4"
<b>SEATING SIZES</b>	
<b>CAPTAIN'S SEATING</b>	<b>Widths:</b> 16", 18", 20" <b>Depths:</b> 16", 18", 20", 22"
<b>SEAT-TO-FLOOR HEIGHTS</b> <sup>10</sup>	
<b>CAPTAIN'S SEATING</b>	Lowest, middle, highest
<b>BATTERY SIZE</b> <sup>3</sup>	NF-22 (2)
<b>BATTERY WEIGHT</b> <sup>7</sup>	40 lbs. each
<b>AVAILABLE ELECTRONICS</b> <sup>8</sup>	120A Q-Logic 3 75A Q-Logic NE 75A Q-Logic NE+
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	Standard: 2-Pole Optional: 4-Pole
<b>WEIGHT CAPACITY</b> <sup>4</sup>	300 lbs.
<b>BASE WEIGHT</b>	138.6 lbs.
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard (K0848/K0849)

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.

2) Without front riggings

3) Battery capacity may change when adding power positioning systems.

4) Weight capacity is dictated by the selection of seating.

5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.  
Listed measurements are in the lowest position; does not include elevating system unless specified.

6) 7 mph speed requires optional 4-pole motors and Q-Logic 3 EX or Q-Logic 3e controller.

7) Battery weight may vary +/- 3 lbs. based on manufacturer.

8) All Quantum Rehab® units can be ordered with specialty controls.

**PROUDLY MADE IN**  
**AMERICA**  
**OF FOREIGN AND DOMESTIC PARTS**

\*The Quantum® Edge 3 is an FDA Class II Medical Device designed to aid individuals with mobility impairments.

**QUANTUM®**

#1 FOR REHAB POWER

(US) 866-800-2002

(CAN) 888-570-1113

quantumrehab.com




## Features

- Upgraded SRS (Smooth Ride Suspension)
- 7 mph speed with new 4-pole motor option
- Fender lights and USB charger included
- Wide array of color choices
- Made in America of foreign and domestic parts
- Back depth adjustable
- Limited recline hinges from 85-115 degrees
- Can be folded flat for transport
- Armrests offer 3" of height adjustment and 10 degrees of angle adjustment

## Options

- iLevel® option allows 12" of seat elevation while at walking speed (3.5 mph)
- Q-Logic 3 EX Drive Controls
- Cup holder
- Oxygen holder
- Cane/crutch holder
- Rear basket
- WC-19 Transit
- Accu-Trac Advanced Tracking Technology

## Colors



## Hydrographics

