

EDGE® 3

EDGE® 3 shown in *Blue Jean* with captain's seating

NOW AVAILABLE

MAXIMUM SPEED UP TO 7 MPH
with new optional 4-pole motors



Upgraded SRS (Smooth Ride Suspension) provides a smooth, comfortable ride and advanced stability. Optional iLevel® is available at 4.5 mph with up to 12" of power adjustable seat height. Standard features include front and rear LED fender lights, a USB charger, and drive wheel color accents. The Edge® 3 is available in a total of 20 color options, including seven bold colors.



BLUE JEAN

Q
QUANTUM®
#1 FOR REHAB POWER



quantumrehab.com

Specifications

| | |
|--|---|
| DRIVE WHEELS | Standard: 14" Solid Optional: 14" Pneumatic; 14" Knobby |
| CASTER WHEELS | Front: 6" Solid Rear: 6" Solid |
| ANTI-TIP WHEELS | N/A |
| MAXIMUM SPEED ¹ | Up to 7 mph ¹¹ Up to 4.5 mph w/iLevel activated |
| GROUND CLEARANCE | 2.7" (battery tray) |
| TURNING RADIUS ² | 20.5" |
| OVERALL LENGTH ² | 32.23" without front riggings |
| BASE WIDTH | 24.4" |
| SEATING SIZES | |
| CAPTAIN'S SEATING | Width: 16", 18", 20" Depth: 16-18", 18-20" |
| SEAT-TO-FLOOR HEIGHTS ¹⁰ | |
| CAPTAIN'S SEATING | Lowest, middle, highest |
| MANUAL RECLINE ⁶ | No |
| BATTERY SIZE ³ | NF-22 (2) |
| BATTERY WEIGHT ⁷ | 38 lbs. each |
| AVAILABLE ELECTRONICS ⁸ | 120A Q-Logic 3 EX 120A Q-Logic 3e 75A Q-Logic NE 75A Q-Logic NE+ |
| BATTERY CHARGER | 8A Off-board |
| MOTOR PACKAGES | Standard: 2-Pole, Mid-Wheel 6® Optional: 4-Pole, Mid-Wheel 6® |
| WEIGHT CAPACITY ⁴ | 300 lbs. |
| BASE WEIGHT | 138.6 lbs. |
| MEDICARE CODE COVERAGE | Group 3: Standard (K0848/K0849) |

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability.
- 10) 16.5" seat-to-floor height is not available with power tilt or iLevel. 19.5" seat-to-floor height is not available with static seat or power recline.
- 11) 7 mph speed requires optional 4-pole motors and Q-Logic 3 EX or Q-Logic 3e controller.

Features

- Upgraded SRS (Smooth Ride Suspension)
- 7 mph speed with new 4-pole motor option
- Fender lights and USB charger included
- 20 bold, contemporary color choices
- Made in America of foreign and domestic parts
- Back depth adjustable
- Limited recline hinges from 85-115 degrees
- Can be folded flat for transport
- Armrests offer 3" of height adjustment and 10 degrees of angle adjustment

Options

- iLevel® option allows 12" of seat elevation while at walking speed (3.5 mph)
- Q-Logic 3 EX Drive Controls
- Cup holder
- Oxygen holder
- Cane/crutch holder
- Rear basket
- WC-19 Transit
- Accu-Trac Advanced Tracking Technology

Colors

Visit quantumrehab.com to see all 20 colors.



BLUE JEAN CHERRY BOMB LEMON CRUSH GREEN MACHINE



RASPBERRY BERET EVOLUTION ORANGE SUGAR PLUM

HYDROGRAPHICS (Not available on side shrouds)



CARBON FIBER DIGITAL CAMO FADED GLORY