

EDGE® 3

EDGE® 3 shown in *Blue Jean* with TRU-Balance® 3 seating system

NOW AVAILABLE

MAXIMUM SPEED UP TO 7 MPH
with new optional 4-pole motors



Upgraded SRS (Smooth Ride Suspension) provides a smooth, comfortable ride and advanced stability. Optional iLevel® is available at 4.5 mph with up to 12" of power adjustable seat height. Standard features include front and rear LED fender lights, a USB charger, and drive wheel color accents. The Edge 3 is available in a total of 20 color options, including seven bold colors.

QUANTUM



BLUE JEAN

Q
QUANTUM®
#1 FOR REHAB POWER



quantumrehab.com

Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic; 14" Knobby
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Up to 7 mph ¹¹ Up to 4.5 mph w/iLevel activated
GROUND CLEARANCE	2.7" (battery tray)
TURNING RADIUS ²	20.5"
OVERALL LENGTH ²	32.23" without front riggings
BASE WIDTH	24.4"
SEATING SIZES	
Synergy®/TRU-Balance® 3 Static Seat	Width: 10-22" Depth: 10-22"
TRU-Balance® 3 Power Tilt	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Tilt and iLevel®	Width: 12-22" Depth: 12-22"
TRU-Balance® Flex Seating	Width: 15-22" Depth: 15-20"
SEAT-TO-FLOOR HEIGHTS ¹⁰	
Synergy®/TRU-Balance® 3 Static Seat ⁵	Static: 16.5 - 18.5" Synergy: 15.875 - 20.875"
TRU-Balance® 3 Power Tilt ⁵	17.5-19.5"
TRU-Balance® 3 Tilt & iLevel® ⁵	17.5-19.5"
TRU-Balance® Flex Seating	16.5 - 18.25"
MANUAL RECLINE ⁶	No
BATTERY SIZE ³	NF-22 (2)
BATTERY WEIGHT ⁷	38 lbs. each
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 EX 120A Q-Logic 3e 75A Q-Logic NE 75A Q-Logic NE+
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole, Mid-Wheel 6®
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	138.6 lbs.
MEDICARE CODE COVERAGE	Group 3: Standard (K0848/K0849), Single Power (K0856/K0857), Multiple Power (K0861)



1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
 2) Without front riggings
 3) Battery capacity may change when adding power positioning systems.
 4) Weight capacity is dictated by the selection of seating.
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
 8) All Quantum Rehab® units can be ordered with specialty controls.
 9) Call for availability
 10) 16.5" seat-to-floor height is not available with power tilt or iLevel®. 19.5" seat-to-floor height is not available with static seat or power recline.
 11) 7 mph speed requires optional 4-pole motors and Q-Logic 3 EX or Q-Logic 3e controller.

Features

- Upgraded SRS (Smooth Ride Suspension)
- 7 mph speed with new 4-pole motor option
- 4.5 mph with iLevel®
- Fender lights and USB charger included
- Bold new colors
- Made in America of foreign and domestic parts
- Available with optional WC19 crash-test approved occupied transit kit
- Available with optional Accu-Trac Advanced Tracking Technology

Colors

Visit quantumrehab.com to see all 20 colors.



HYDROGRAPHICS (Not available on side shrouds)



PROUDLY MADE IN AMERICA
OF FOREIGN AND DOMESTIC PARTS