

EDGE® 3

EDGE® 3 shown in *Blue Jean* with
TRU-Balance® 3 seating system

NOW AVAILABLE

MAXIMUM SPEED UP TO 7 MPH
with new optional 4-pole motors



FDA Class II
Medical Device*

Upgraded SRS (Smooth Ride Suspension) provides a smooth, comfortable ride and advanced stability. Optional iLevel® is available at 4.5 mph with up to 12" of power adjustable seat height. Standard features include front and rear LED fender lights, a USB charger, and drive wheel color accents. TRU-Balance® 4 Power Positioning with anterior tilt and memory seating is now available!

QUANTUM



BLUE JEAN

Q
QUANTUM®
#1 FOR REHAB POWER



| quantumrehab.com

Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic; 14" Knobby
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Standard: Up to 6.25 mph ⁶ Optional: Up to 7 mph ⁹ (4.5 mph with iLevel [®])
GROUND CLEARANCE	2.7" (battery tray)
TURNING RADIUS ²	20.5"
OVERALL LENGTH ²	35.5"
BASE WIDTH	24.4"
SEATING SIZES	
TRU-Balance[®] 4	Width: 16-22" Depth: 14-22"
Synergy[®] Seat	Width: 10-22" Depth: 10-22"
TRU-Balance[®] 3 Static Seat	Width: 12-22" Depth: 12-22"
TRU-Balance[®] 3 Power Tilt	Width: 12-22" Depth: 12-22"
TRU-Balance[®] 3 Tilt and iLevel[®]	Width: 12-22" Depth: 12-22"
TRU-Balance[®] Flex Seating	Width: 15-22" Depth: 15-20"
SEAT-TO-FLOOR HEIGHTS ¹⁰	
TRU-Balance[®] 4	17.5-19.25"
Synergy[®] Seat	15.875-20.875"
TRU-Balance[®] 3 Static Seat ⁵	16.5-18.5"
TRU-Balance[®] 3 Power Tilt ⁵	17.5-19.5"
TRU-Balance[®] 3 Tilt & iLevel[®] ⁵	17.5-19.5"
TRU-Balance[®] Flex Seating	16.5 - 18.25"
BATTERY SIZE ³	NF-22 (2)
BATTERY WEIGHT ⁷	40 lbs. each
AVAILABLE ELECTRONICS ⁸	120A Q-Logic 3 75A Q-Logic NE 75A Q-Logic NE+
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-pole and 4-pole
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	138.6 lbs.
MEDICARE CODE COVERAGE	Group 3: Standard (K0848/K0849), Single Power (K0856/K0857), Multiple Power (K0861)


- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
2) Without front riggings
3) Battery capacity may change when adding power positioning systems.
4) Weight capacity is dictated by the selection of seating.
5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.
Listed measurements are in the lowest position; does not include elevating system unless specified.
6) 6.25 mph speed available only with Q-Logic 3.
7) Battery weight may vary +/- 3 lbs. based on manufacturer.
8) All Quantum Rehab[®] units can be ordered with specialty controls.
9) Up to 7 mph speed available only with 4-pole motor option and Q-Logic 3 or Q-Logic 3e.
10) Up to 7 mph speed available only with 4-pole motor option and Q-Logic 3 or Q-Logic 3e.

*The Quantum[®] Edge 3 is an FDA Class II Medical Device designed to aid individuals with mobility impairments.


Features

- Upgraded SRS (Smooth Ride Suspension)
- 7 mph speed with new 4-pole motor option
- Up to 4.5 mph with iLevel[®]
- Fender lights and USB charger included
- Made in America of foreign and domestic parts
- Available with optional WC19 crash-test approved occupied transit kit
- Available with optional Accu-Trac Advanced Tracking Technology

Colors



Hydrographics



PROUDLY MADE IN
AMERICA
OF FOREIGN AND DOMESTIC PARTS

QUANTUM[®]

#1 FOR REHAB POWER

(US) 866-800-2002

(CAN) 888-570-1113

quantumrehab.com

