

# EDGE 3 *Stretto*®



EDGE 3 STRETTO®  
shown in *Sugar Plum* with  
captain's seating

The Edge 3 Stretto® with Captain's seating features a narrow width for outstanding maneuverability and is equipped with independent SRS (Smooth Ride Suspension), making it ideal for children, teenagers and small adults. The Stretto is available with 12.5" or 14" drive wheels and optional 12" of iLevel® power adjustable seat height.

QUANTUM



Narrowest  
Base in  
America

**QUANTUM**®  
#1 FOR REHAB POWER



[quantumrehab.com](http://quantumrehab.com)

# Specifications

<b>DRIVE WHEELS</b>	14" or 12.5" pneumatic or flat free tires
<b>CASTER WHEELS</b>	<b>Front:</b> 5" flat free <b>Rear:</b> 6" flat free
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 6.25 mph Up to 3.5 mph w/iLevel activated
<b>GROUND CLEARANCE</b>	2.9"
<b>TURNING RADIUS</b> <sup>2</sup>	21"
<b>OVERALL LENGTH</b> <sup>2</sup>	32.5" without front riggings
<b>BASE WIDTH</b>	21.75" with 14" drive wheels 20.47" with 12.5" drive wheels
<b>SEATING SIZES</b>	
<b>Captain's seating</b>	<b>Width:</b> 16", 18", 20" <b>Depth:</b> 16", 18", 20"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>No iLevel® with standard seat cushion</b> <sup>5</sup>	18.5" - 20.5"
<b>No iLevel® without standard seat cushion</b> <sup>5</sup>	15.5" - 17.5"
<b>iLevel® with standard seat cushion</b> <sup>5</sup>	20.5" - 22.5"
<b>iLevel® without standard seat cushion</b> <sup>5</sup>	17.5" - 19.5"
<b>MANUAL RECLINE</b> <sup>6</sup>	No
<b>BATTERY SIZE</b> <sup>3</sup>	40 AH (2) or 55 AH (2)
<b>BATTERY WEIGHT</b> <sup>7</sup>	40 AH: 31 lbs. each 55 AH: 40 lbs. each
<b>AVAILABLE ELECTRONICS</b> <sup>8</sup>	120A Q-Logic 3 EX 120A Q-Logic 3e 75A Q-Logic NE 75A Q-Logic NE+
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-Pole, Mid-Wheel 6®
<b>WEIGHT CAPACITY</b> <sup>4</sup>	300 lbs.
<b>BASE WEIGHT</b>	126.4 lbs.
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard (K0848/K0849)

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.  
 2) Without front riggings  
 3) Battery capacity may change when adding power positioning systems.  
 4) Weight capacity is dictated by the selection of seating.  
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.  
 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.  
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.  
 8) All Quantum Rehab® units can be ordered with specialty controls.  
 9) Call for availability.

# Features

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- Fender lights and USB charger included
- 19 bold, contemporary color choices
- Made in America of foreign and domestic parts
- Back depth adjustable
- Limited recline hinges from 85-115 degrees
- Can be folded flat for transport
- Armrests offer 3" of height adjustment and 10 degrees of angle adjustment

# Options

- iLevel® option allows 12" of seat elevation while at walking speed (3.5 mph)
- Q-Logic 3 EX Drive Controls
- Cup holder
- Oxygen holder
- Cane/crutch holder
- Rear basket
- WC-19 transit kit
- Accu-Trac Advanced Tracking Technology

# Colors



# Hydrographics

