

EDGE 3 *Stretto*®



EDGE 3 STRETTO®
shown in Green Machine with
captain's seating

FDA Class II
Medical Device *

The Edge 3 Stretto® with Captain's seating features a narrow width for outstanding maneuverability and is equipped with independent SRS (Smooth Ride Suspension), making it ideal for children, teenagers and small adults. The Stretto is available with 12.5" or 14" drive wheels and optional 12" of iLevel® power adjustable seat height. TRU-Balance® 4 Power Positioning with anterior tilt and memory seating is now available!

QUANTUM



Narrowest
Base in
America


QUANTUM®



| quantumrehab.com

Specifications

DRIVE WHEELS	14" or 12.5" pneumatic or flat free tires Optional: 14" knobby
CASTER WHEELS	Front: 5" flat free Rear: 6" flat free
MAXIMUM SPEED ¹	Up to 6.25 mph (3.5 mph with iLevel®)
GROUND CLEARANCE	2.9"
TURNING RADIUS ²	21"
OVERALL LENGTH ²	32.5"
BASE WIDTH	21.5" with 14" drive wheels 20.47" with 12.5" drive wheels
SEATING SIZES	
Captain's seating	Width: 16", 18", 20" Depth: 16", 18", 20"
SEAT-TO-FLOOR HEIGHTS	
Lowest, Middle, Highest	
BATTERY SIZE ³	40 AH (2) or NF-22 (2)
BATTERY WEIGHT ⁷	40 AH: 31 lbs. each NF-22: 40 lbs. each
AVAILABLE ELECTRONICS ⁶	120A Q-Logic 3 75A Q-Logic NE 75A Q-Logic NE+
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	126.4 lbs.
MEDICARE CODE COVERAGE	Group 3: Standard (K0848/K0849)

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
 2) Without front riggings
 3) Battery capacity may change when adding power positioning systems.
 4) Weight capacity is dictated by the selection of seating.
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.
 Listed measurements are in the lowest position; does not include elevating system unless specified.
 6) All Quantum Rehab® units can be ordered with specialty controls.
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.

Features

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- Fender lights and USB charger included
- Made in America of foreign and domestic parts
- Back depth adjustable
- Limited recline hinges from 85-115 degrees
- Can be folded flat for transport
- Armrests offer 3" of height adjustment and 10 degrees of angle adjustment

Options

- iLevel® option allows 12" of seat elevation while at walking speed (3.5 mph)
- Q-Logic 3 Drive Controls
- Cup holder
- Oxygen holder
- Cane/crutch holder
- Rear basket
- WC-19 transit kit
- Accu-Trac Advanced Tracking Technology

Colors



Hydrographics

