

QUANTUM® **EDGE 4**

THE LEADING EDGE
IN GROUP 3 JUST TOOK
A **QUANTUM** LEAP.



FDA Class II
Medical Device*

FEATURING 8 NEW
COLOR OPTIONS!



“ THE EDGE 4 ISN'T JUST MOBILITY,
IT'S MOMENTUM. IT CARRIES YOU
FORWARD WITH POWER, GRACE, AND THE
KIND OF CONFIDENCE THAT MAKES
THE WORLD FEEL WIDE OPEN. ”

STEWART LUNDY



QUANTUM®



quantumrehab.com

Specifications

DRIVE WHEELS	14" pneumatic or flat free tires Optional: 14" knobby Pneumatic knobby tires not available
CASTER WHEELS	7" flat free
MAXIMUM SPEED ¹	Up to 7 mph (4.5 mph with iLevel®)
RANGE PER CHARGE ^{1,9}	Up to 18.1 miles
GROUND CLEARANCE	2.87"
TURNING RADIUS	22.8"
LENGTH	38.27"
BASE WIDTH	24.02"
SEATING SIZES	
TRU Balance® 4 Anterior Tilt ⁵	Width: 16-22" Depth: 14-22"
TRU Balance® 3 Static Seat	Width: 12-22" Depth: 12-20"
TRU Balance® 3 Power Tilt TRU Balance® 4 Power Tilt	Width: 12-22" Depth: 12-22"
TRU Balance® 3 Tilt and iLevel®	Width: 12-22" Depth: 12-22"
TRU Balance® Flex Seating	Width: 10-22" Depth: 10-20"
SEAT-TO-FLOOR HEIGHTS	
TRU Balance® 4 Anterior Tilt ⁵	18.25-21.25"
TRU Balance® 3 Power Tilt ⁵ TRU Balance® 4 Power Tilt ⁵	17.75-19.75"
TRU Balance® 3 Tilt & iLevel® ⁵	17.5-19.5"
TRU Balance® Flex Seating	16.25-18.25"
BATTERY SIZE ³	
BATTERY WEIGHT ⁷	C5 = 65 Ah C20 = 75 Ah
AVAILABLE ELECTRONICS ⁶	120A Q-Logic 3 90A Q-Logic NES
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-Pole
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	154 lbs. without batteries
MEDICARE CODE COVERAGE	Group 3: Standard (K0848 ⁸ /K0849 ⁸) Single Power (K0857 ⁸ /K0856) Multiple Power (K0861)

- 1) Varies with user weight, terrain type, battery charge, battery condition and tire pressure.
 2) Without front riggings
 3) Battery capacity may change when adding power positioning systems.
 4) Weight capacity is dictated by the selection of seating.
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.
 Listed measurements are in the lowest position; does not include elevating system unless specified.
 6) All Quantum Rehab® units can be ordered with specialty controls.
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
 8) Pending PDAC approval
 9) Range shown is measured using Group 24 batteries

* FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

Features

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 4.5 mph)
- Standard accent lights and USB charger included
- Wide array of color choices
- Made in America of foreign and domestic parts



Options

- iLevel® option allows 12" of seat elevation while at walking speed (4.5 mph)
- Q-Logic 3 Drive Controls
- TRU Balance® 3 & 4 Seating
- Depth-adjustable vent tray
- 12x12" seat sizes with iLevel
- WC19 crash-test approved occupied transit kit (pending test results)
- Accu-Trac Advanced Tracking Technology (coming soon)

Colors



QUANTUM®