

# QUANTUM® **EDGE 4**

THE LEADING EDGE  
IN GROUP 3 JUST TOOK  
A **QUANTUM** LEAP.



FDA Class II  
Medical Device\*

FEATURING 8 NEW  
COLOR OPTIONS!



“ THE EDGE 4 ISN'T JUST MOBILITY,  
IT'S MOMENTUM. IT CARRIES YOU  
FORWARD WITH POWER, GRACE, AND THE  
KIND OF CONFIDENCE THAT MAKES  
THE WORLD FEEL WIDE OPEN. ”

STEWART LUNDY



QUANTUM®



quantumrehab.com

# Specifications

<b>DRIVE WHEELS</b>	14" pneumatic or flat free tires <b>Optional:</b> 14" knobby Pneumatic knobby tires not available
<b>CASTER WHEELS</b>	7" flat free
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 7 mph (4.5 mph with iLevel®)
<b>RANGE PER CHARGE</b> <sup>1,8</sup>	Up to 18.1 miles
<b>GROUND CLEARANCE</b>	2.87"
<b>TURNING RADIUS</b>	22.8"
<b>LENGTH</b>	38.27"
<b>BASE WIDTH</b>	24.02"
<b>SEATING SIZES</b>	
<b>TRU Balance® 4 Anterior Tilt</b> <sup>5</sup>	<b>Width:</b> 16-22" <b>Depth:</b> 14-22"
<b>TRU Balance® 3 Static Seat</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-20"
<b>TRU Balance® 3 Power Tilt</b> <b>TRU Balance® 4 Power Tilt</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU Balance® 3 Tilt and iLevel®</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU Balance® Flex Seating</b>	<b>Width:</b> 10-22" <b>Depth:</b> 10-20"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>TRU Balance® 4 Anterior Tilt</b> <sup>5</sup>	18.25-21.25"
<b>TRU Balance® 3 Power Tilt</b> <sup>5</sup> <b>TRU Balance® 4 Power Tilt</b> <sup>5</sup>	17.75-19.75"
<b>TRU Balance® 3 Tilt &amp; iLevel®</b> <sup>5</sup>	17.5-19.5"
<b>TRU Balance® Flex Seating</b>	16.25-18.25"
<b>BATTERY SIZE</b> <sup>3</sup>	
	C5 = 65 Ah C20 = 75 Ah
<b>BATTERY WEIGHT</b> <sup>7</sup>	47.8 lbs. each
<b>AVAILABLE ELECTRONICS</b> <sup>6</sup>	120A Q-Logic 3 90A Q-Logic NES
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	4-Pole
<b>WEIGHT CAPACITY</b> <sup>4</sup>	300 lbs.
<b>BASE WEIGHT</b>	154 lbs. without batteries
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard (K0848/K0849) Single Power (K0857 / K0856) Multiple Power (K0861)

1) Varies with user weight, terrain type, battery charge, battery condition and tire pressure.  
 2) Without front riggings  
 3) Battery capacity may change when adding power positioning systems.  
 4) Weight capacity is dictated by the selection of seating.  
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.  
 Listed measurements are in the lowest position; does not include elevating system unless specified.  
 6) All Quantum Rehab® units can be ordered with specialty controls.  
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.  
 8) Range shown is measured using Group 24 batteries

\* FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

# Features

- Independent SRS (Smooth Ride Suspension)
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 4.5 mph)
- Standard accent lights and USB charger included
- Wide array of color choices
- Made in America of foreign and domestic parts



# Options

- iLevel® option allows 12" of seat elevation while at walking speed (4.5 mph)
- Q-Logic 3 Drive Controls
- TRU Balance® 3 & 4 Seating
- Depth-adjustable vent tray
- 12x12" seat sizes with iLevel
- WC19 crash-test approved occupied transit kit (pending test results)
- Accu-Trac Advanced Tracking Technology (coming soon)

# Colors



**QUANTUM®**

(US) 866-800-2002 • quantumrehab.com • Follow us on social media for our latest innovations:

