



J4 HD shown in *Electric Blue*  
with Q4 HD Power Tilt Seating  
and NE+ controller

# J4<sup>HD</sup>

With features like Q4 seating, Mid-Wheel 6<sup>®</sup> Drive Design and ATX Suspension, the J4 HD is engineered for tight-quarter maneuverability. The J4 HD delivers rehab capability and performance at an excellent value.

QUANTUM<sup>®</sup>

## Specifications

<b>DRIVE WHEELS</b>	Standard: 14" Flat Free
<b>CASTER WHEELS</b>	Front: 6" Solid Rear: 6" Solid
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED <sup>1</sup></b>	Up to 5 mph
<b>GROUND CLEARANCE</b>	3"
<b>TURNING RADIUS <sup>2</sup></b>	21"
<b>OVERALL LENGTH</b>	32.3"
<b>BASE WIDTH</b>	24.25"
<b>SEATING SIZES</b>	
<b>Q4 Seating</b> (in 2" increments)	<b>Width:</b> 20-24" <b>Depth:</b> 20-24"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>Q4 Seating</b>	17.25-19.25"
<b>BATTERY SIZE <sup>3</sup></b>	
<b>NF-22</b>	
<b>BATTERY WEIGHT <sup>7</sup></b>	
<b>38 lbs.</b>	
<b>AVAILABLE ELECTRONICS <sup>6</sup></b>	
<b>75A Q-Logic NE+</b>	
<b>BATTERY CHARGER</b>	
<b>8A Off-board</b>	
<b>MOTOR PACKAGES</b>	
<b>Two Motor, Mid-Wheel 6<sup>®</sup></b>	
<b>WEIGHT CAPACITY <sup>4</sup></b>	
<b>400 lbs.</b>	
<b>BASE WEIGHT</b>	
<b>131.6 lbs.</b>	
<b>MEDICARE CODE COVERAGE</b>	
<b>Group 2:</b> <b>Single Power HD (K0837)</b>	

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) All Quantum Rehab<sup>®</sup> units can be ordered with specialty controls
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.

## Colors



Red Corvette



Electric Blue



Back in Black

## Q4 HD Power Tilt Features



**QUANTUM<sup>®</sup>**

#1 FOR REHAB POWER

(US) 866-800-2002

(Canada) 888-570-1113

quantumrehab.com

