

# 1 F O R R E H A B P O W E R



FDA Class II  
Medical Device

**J4 HD** shown in *Electric Blue*  
with Q4 Power Tilt Seating  
and NE+ controller

# J4<sup>HD</sup>

With features like Q4 seating, Mid-Wheel 6® Drive Design and  
ATX Suspension, the J4 HD is engineered for tight-quarter maneuverability.  
The J4 HD delivers rehab capability and performance at an excellent value.

**QUANTUM®**

## Specifications

<b>DRIVE WHEELS</b>	14" Flat Free
<b>CASTER WHEELS</b>	<b>Front:</b> 6" Solid <b>Rear:</b> 6" Solid
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED <sup>1</sup></b>	Up to 5 mph
<b>GROUND CLEARANCE</b>	2.9" (battery tray)
<b>TURNING RADIUS <sup>2</sup></b>	21"
<b>OVERALL LENGTH</b>	35.5"
<b>BASE WIDTH</b>	24.5"
<b>SEATING SIZES</b>	
<b>Q4 Seating</b> (in 2" increments)	<b>Width:</b> 20-24" <b>Depth:</b> 20-24"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>Q4 Seating</b>	17.25-19.25"
<b>BATTERY SIZE <sup>3</sup></b>	
<b>BATTERY WEIGHT <sup>7</sup></b>	40 lbs. each.
<b>AVAILABLE ELECTRONICS <sup>6</sup></b>	75A Q-Logic NE+
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	Two Motor
<b>WEIGHT CAPACITY <sup>4</sup></b>	400 lbs.
<b>BASE WEIGHT</b>	131.6 lbs.
<b>MEDICARE CODE COVERAGE</b>	<b>Group 2:</b> Single Power HD (K0837)

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) All Quantum Rehab® units can be ordered with specialty controls
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.

Quantum Rehab® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

## Colors



Red  
Corvette



Electric  
Blue



Back in  
Black

## Q4 Power Tilt Features



**QUANTUM®**

#1 FOR REHAB POWER

(US) 866-800-2002

(Canada) 888-570-1113

quantumrehab.com

