



J4 shown in Red Corvette
with Q4 Power Tilt Seating

FDA Class II
Medical Device*

J4

With features like Q4 seating, Mid-Wheel 6® Drive Design and ATX Suspension, the J4 is engineered for tight-quarter maneuverability. The J4 delivers rehab capability and performance at an excellent value.

QUANTUM®

Specifications

DRIVE WHEELS	14" Flat Free
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Up to 5 mph
GROUND CLEARANCE	2.9" at battery tray
TURNING RADIUS ²	21"
OVERALL LENGTH	35.5"
BASE WIDTH	24.5"
SEATING SIZES	
Q4 Seating	Width: 16-22" (in 2" increments) Depth: 16-22" (in 2" increments)
TRU-Balance® Flex Seating	Width: 15-22" Depth: 15-20"
SEAT-TO-FLOOR HEIGHTS	
Q4 Seating	17.25-19.25"
TRU-Balance® Flex Seating	16.5-18.25"
BATTERY SIZE ³	
	NF-22
BATTERY WEIGHT ⁷	
	40 lbs.
AVAILABLE ELECTRONICS ⁶	
	75A Q-Logic NE 75A Q-Logic NE+ 75A Q-Logic NES
BATTERY CHARGER	
	8A Off-board
MOTOR PACKAGES	
	Two Motor
WEIGHT CAPACITY ⁴	
	300 lbs.
BASE WEIGHT	
	131.6 lbs.
MEDICARE CODE COVERAGE	
	Group 2: Single Power (K0835)

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) All Quantum Rehab® units can be ordered with specialty controls
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.

Colors



Red Corvette



Electric Blue



Back in Black

*The Quantum® J4 is an FDA Class II Medical Device designed to aid individuals with mobility impairments.

QUANTUM®

#1 FOR REHAB POWER

(US) 866-800-2002

(CAN) 888-570-1113

quantumrehab.com



Features

- Standard 5 mph motors offer advanced performance, precise control and responsiveness
- Bold aesthetics
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters for enhanced performance over more varied terrain

Options

- Q4 Power Tilt
- Power Articulating Foot Platform



J4 shown in *Red Corvette* with Q4 Power Tilt seating

FDA Class II Medical Device*