


# J4E HD

Group 2: Single Power (K0837)



FDA Class II Medical Device<sup>6</sup>

## Features

- 400 lb. weight capacity
- Standard 5 mph motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics
- Available with up to 10" of seat elevation
- Safe Seat Elevation speed of up to 1.25 mph

## Seating Options

- Q4 50° Tilt
- Power articulating foot platform

## Specifications

<b>DRIVE WHEELS</b>	14" Solid
<b>CASTER WHEELS</b>	<b>Front:</b> 6" Solid <b>Rear:</b> 6" Solid
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 5 mph
<b>RANGE</b>	Up to 16.4 miles
<b>GROUND CLEARANCE</b>	2.7" (battery tray)
<b>TURNING RADIUS</b> <sup>2</sup>	23"
<b>OVERALL LENGTH</b> <sup>2</sup>	35.5"
<b>BASE WIDTH</b>	24.4"
<b>SEATING SIZES</b> (Q4 LIFT TILT)	<b>Width:</b> 20-24" <b>Depth:</b> 20-24"
<b>SEAT-TO-FLOOR HEIGHTS</b> <sup>4</sup> (Q4 LIFT TILT) <sup>4</sup>	17.5-19.5"
<b>BATTERY SIZE</b> <sup>3</sup>	NF-22
<b>BATTERY WEIGHT</b> <sup>5</sup>	38 lbs.
<b>AVAILABLE ELECTRONICS</b>	120A Q-Logic 3e 75A Q-Logic NES
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	4-Pole motors
<b>WEIGHT CAPACITY</b>	400 lbs.
<b>BASE WEIGHT</b>	138.6 lbs.
<b>MEDICARE CODE COVERAGE</b>	<b>Group 2:</b> Single Power (K0837)

## Available Colors



Red Corvette



Electric Blue



Back in Black

1) Speed varies with user weight, terrain type, battery charge/condition, and tire pressure.

2) Without front riggings

3) Battery capacity may change when adding power positioning systems.

4) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 400 lbs. weight capacity and solid tires. Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.

5) Battery weight may vary +/- 3 lbs. based on manufacturer.

6) Quantum Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.



quantumrehab.com

866.800.2002

QUANTUM®