

J4^E

Group 2: Single Power (K0835)



Features

- Standard 5 mph motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics
- Available with up to 10 inches of seat elevation

Seating Options

- Q4 50° Tilt
- Power articulating foot platform
- 10 inch seat elevation

Specifications

DRIVE WHEELS	14" Solid
CASTER WHEELS	Front: 6" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ¹	Up to 5 mph
RANGE	Up to 16.1 miles
GROUND CLEARANCE	3" (battery tray)
TURNING RADIUS ²	20.5"
OVERALL LENGTH ²	35.5"
BASE WIDTH	24.4"
SEATING SIZES (Q4 LIFT TILT)	Width: 16", 18", 20" or 22" Depth: 16-22"
SEAT-TO-FLOOR HEIGHTS ⁵ (Q4 LIFT TILT) ⁵	17.5-19.5"
BATTERY SIZE ³	NF-22
BATTERY WEIGHT ⁷	38 lbs.
AVAILABLE ELECTRONICS	120A Q-Logic 3e 75A Q-Logic NE+
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	2-Pole Mid-Wheel 6 [®]
WEIGHT CAPACITY ⁴	300 lbs.
BASE WEIGHT	138.6 lbs.
MEDICARE CODE COVERAGE	Group 2: Single Power (K0835)

Available Colors



Red Corvette



Electric Blue



Back in Black

- 1) Speed varies with user weight, terrain type, battery charge/condition, and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires. Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.



quantumrehab.com

866.800.2002

QUANTUM[®]
#1 FOR REHAB POWER