



J6
shown in Red Corvette with
TRU-Balance® 3 Power Tilt

FDA Class II
Medical Device*

J6

The compact J6 is the ultimate choice for clients who demand the superb stability of Mid-Wheel 6® Design and tight-quarter maneuverability. Highly adaptable, the J6 accepts a wide range of seating and electronics options while its powerful drivetrain delivers high-performance torque.

QUANTUM®

Specifications

DRIVE WHEELS	Standard: 10" Solid Optional: 10" Pneumatic
CASTER WHEELS	Front: 5" Solid Rear: 6" Solid
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ⁴	Up to 4 mph
GROUND CLEARANCE	2.5" (frame)
TURNING RADIUS ⁵	22"
OVERALL LENGTH ⁵	34.74"
BASE WIDTH	23.19"
SEATING SIZES	
Synergy[®]/Static Seat	Width: 10-22" Depth: 10-22"
TRU-Balance[®] 3 Power Tilt	Width: 14-20" Depth: 14-20"
SEAT-TO-FLOOR HEIGHTS	
Synergy[®]/Static Seat ⁸	Synergy [®] : 15.875-17.875" Static: 16.75"
TRU-Balance[®] 3 Power Tilt ⁸	16.75"
BATTERY SIZE ⁶	U-1
BATTERY WEIGHT ¹⁰	23.4 lbs.
AVAILABLE ELECTRONICS ³	70A Q-Logic NE+ 70A Q-Logic NE
BATTERY CHARGER	5A Off-board
MOTOR PACKAGES	2-Pole
WEIGHT CAPACITY ⁷	300 lbs. ⁹
BASE WEIGHT	100.2 lbs.
MEDICARE CODE COVERAGE	Group 2: Standard (K0822/K0823), Single Power (K0835)

- 1) Standard
- 2) Optional
- 3) All Quantum Rehab[®] units can be ordered with specialty controls
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems.
- 7) Weight capacity is dictated by the selection of seating.
- 8) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 9) 250 lbs. weight capacity with TRU-Balance[®] Systems
- 10) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 11) Requires Group 34 batteries
- 12) All Quantum Rehab[®] units can be ordered with specialty controls

Features

- Mid-Wheel 6[®] allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance[®] 3 Power Positioning
- Side-mounted, easily accessible freewheel levers
- ATX Suspension (Active-Trac[®] with extra stability) incorporates enhanced performance over more varied terrain
- Easy front access to batteries

Options

- Swing-away joystick mount
- Multi-axis joystick mount
- Elevating leg rests
- Oxygen tank holder
- Neoprene leg rest sleeves
- Gel adductor buttons
- Power elevating seat
- Angle-adjustable footplates
- 70° swing-away leg rests

Colors



Red
Corvette



Electric
Blue



Back in
Black



Evolution
Orange



Pretty in
Pink

*The Quantum[®] J6 is an FDA Class II Medical Device designed to aid individuals with mobility impairments.

QUANTUM[®]

#1 FOR REHAB POWER

(US) 866-800-2002

(CAN) 888-570-1113

quantumrehab.com

