

P E D I A T R I C P O W E R C H A I R S O L U T I O N S

KOZMO

shown in *Rose Quartz*
with Pediatric Seating



KOZMO™

The all new Kozmo™ is engineered to meet the needs of even the most active child. The Kozmo's highly compact size provides easy access to kid-sized spaces, while kid-friendly color choices add a fun, personal touch. Whether on the playground or in the home, your child will have a blast with the Quantum® Kozmo.

QUANTUM®

Specifications

| | |
|---|---|
| DRIVE WHEELS | 9" Solid |
| CASTER WHEELS | Front: 6" Solid |
| ANTI-TIP WHEELS | 3" Solid |
| MAXIMUM SPEED ¹ | Up to 3.6 mph |
| GROUND CLEARANCE ² | 1.6" at anti-tip bracket 1.9" at motor |
| TURNING RADIUS ² | 25.45" with foot platform |
| OVERALL LENGTH | 33.09" with front rigging |
| BASE WIDTH ⁵ | 22.6" |
| SEAT TYPE | Kozmo™ pediatric, foldable, static seat |
| SEAT-TO-FLOOR HEIGHT | 15.75"-16.75" |
| BATTERY SIZE ^{3,4} | (2) 18A, 12 V, deep cycle |
| ELECTRONICS | 40A, PG GC3 (non-programmable) |
| BATTERY CHARGER | Off-board, 3.5A |
| WEIGHT CAPACITY | 200 lbs. |
| BASE WEIGHT | 67 lbs. |
| STANDARD SEAT WEIGHT | 30.6 lbs. |
| BATTERY PACK WEIGHT ⁴ | 31 lbs. |
| MEDICARE CODE COVERAGE | Group 2: Standard (K0822) |

1) Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition tire type and condition. These specifications can be subject to a variance of +10%, -5%.

2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

3) AGM or Gel-Cell type required. *See "Batteries and Charging."

4) Battery weight may vary based on manufacturer.

5) Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA and ISO 7176 series standards.



Features

- Choice of 6 colors
- 200 lbs. weight capacity
- Maximum speed up to 3.6 mph
- 18 amp batteries
- Compact turning radius of 25.45"
- Pediatric, height-adjustable footplate
- 38" lap belt


Options

- Cup Holder
- Weather Cover
- Additional Battery Pack
- Rear Mount Bracket Cane/Crutch Holder
- Glove Box
- Quantum Backpack
- Swingaway Pediatric Legrests
- Oxygen Holder
- Flag
- Back Pack
- Tablet Holder
- Cell/Media Holder

Colors




Ruby Red



Sapphire Blue



Rose Quartz



Citrine Yellow



Amber Orange



Pearl White