

PEDIATRIC POWER CHAIR SOLUTIONS



The new Kozmo™ Portable features one-hand, feather-touch disassembly, for quick and convenient transport or storage.



KOZMO™ **PORTABLE**

The Kozmo™ Portable is engineered to meet the needs of even the most active child. The Kozmo Portable's highly compact size provides easy access to kid-sized spaces, while kid-friendly color choices add a fun, personal touch. One-hand, feather-touch disassembly makes transport fast and convenient for parents on the go.

QUANTUM®

Specifications

DRIVE WHEELS	9" Solid
CASTER WHEELS	Front: 6" Solid
ANTI-TIP WHEELS	3" Solid
MAXIMUM SPEED ¹	Up to 3.6 mph
GROUND CLEARANCE ²	1.6" at anti-tip bracket 1.9" at motor
TURNING RADIUS ²	25.45" with foot platform
OVERALL LENGTH	33.09" with front rigging
BASE WIDTH ⁵	22.6"
SEAT TYPE	Kozmo™ pediatric, foldable, static seat
SEAT-TO-FLOOR HEIGHT	15.75"-16.75"
BATTERY SIZE ^{3,4}	(2) 18A, 12 V, deep cycle
ELECTRONICS	40A, PG GC3 (non-programmable)
BATTERY CHARGER	Off-board, 3.5A
MOTOR PACKAGES	Two-motor, rear-wheel drive
WEIGHT CAPACITY	200 lbs.
BASE WEIGHT	67 lbs.
STANDARD SEAT WEIGHT	48.5 lbs.
BATTERY PACK WEIGHT ⁴	31 lbs.
MEDICARE CODE COVERAGE	Group 2: Standard (K0822)

1) Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition tire type and condition. These specifications can be subject to a variance of +10%, -5%.

2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

3) AGM or Gel-Cell type required. *See "Batteries and Charging."

4) Battery weight may vary based on manufacturer.

5) Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA and ISO 7176 series standards.



Features

- Choice of 6 colors
- 200 lbs. weight capacity
- Maximum speed up to 3.6 mph
- 18 amp batteries
- Compact turning radius of 25.45"
- Pediatric, height-adjustable footplate
- One-hand feather-touch disassembly
- 38" lap belt

Options

- Cup Holder
- Weather Cover
- Additional Battery Pack
- Rear Mount Bracket Cane/Crutch Holder
- Glove Box
- Quantum Backpack
- Swingaway Pediatric Legrests
- Oxygen Holder
- Flag
- Tablet Holder
- Cell/Media Holder

Colors




Ruby Red




Sapphire Blue



Rose Quartz



Citrine Yellow



Amber Orange



Pearl White

QUANTUM®

#1 FOR REHAB POWER

(US) 866-800-2002

(CAN) 888-570-1113

quantumrehab.com

