

**EDGE 3 STRETTO®
CONFIGURATION WEIGHTS:
FULLY LOADED**

CONFIGURATION	WEIGHT
iLevel, tilt, recline, AFP, 22NF batteries, 12.5" flat free tires, Redesigned back	408 lbs
iLevel, tilt, recline, AFP, 22NF batteries, 12.5" flat free tires, Redesigned back, Q-Back pan	406 lbs
iLevel, tilt, recline, AFP, 40AH batteries, 12.5" flat free tires, redesigned back	392 lbs
iLevel, tilt, recline, AFP, 40AH batteries, 12.5" flat free tires, redesigned back, Q-back pan	390 lbs
iLevel, tilt, recline, AFP, 22NF, 14" flat free tires, Redesigned back	413 lbs
iLevel, tilt, recline, AFP, 22NF batteries, 14" flat free tires, redesigned back, Q-Back pan	411 lbs
iLevel, tilt, recline, AFP, 40AH batteries, 14" flat free tires, redesigned back	397 lbs
iLevel, tilt, recline, AFP 40AH batteries, 14" flat free tires, redesigned back, Q-Back pan	395 lbs

**EDGE 3
3 STRETTO®
CONFIGURATION
WEIGHTS: TILT,
RECLINE, AFP**

CONFIGURATION	WEIGHT
Tilt, recline, AFP, 22NF batteries, 14" flat free tires, redesigned back	363 lbs
Tilt, recline, AFP, 22NF batteries, 14" flat free tires, Q-Back	361 lbs
Tilt, recline, AFP, 40AH batteries, 14" flat free tires, redesigned back	347 lbs
Tilt, recline, AFP, 40AH batteries, 14" flat free tires, Q-Back	345 lbs


Note: All configuration weights are taken with the headrest, seat cushion, and back cushion removed.

EDGE 3
Stretto®

THE ULTIMATE PEDIATRIC POWER BASE



COLORS*



BLUE JEAN

GREEN MACHINE

LEMON CRUSH

RASPBERRY BERET

EVOLUTION ORANGE

* Visit quantumrehab.com to see all colors.



QUANTUM®
#1 FOR REHAB POWER



866-800-2002 • quantumrehab.com • lifeatlevel.com



Q
QUANTUM®
#1 FOR REHAB POWER

EDGE 3 STRETTO® FAQs

WHY WAS THE EDGE 3 STRETTO® CREATED?

Clinicians and providers have been asking Quantum for a true pediatric power base. While existing Quantum bases were able to accommodate pediatric seating, the Stretto has a narrower width. The Stretto has an overall width of 20.47" with 12.5" drive wheels, making it the narrowest power base in America. It is a milestone for pediatric power!

WHAT IS SRS?

SRS (*Smooth Ride Suspension*) is a superior power chair suspension developed after two years of research and consumer feedback.

Fun Fact...

Stretto is the Italian word for *narrow*.

Edge 3 Stretto® featuring independent SRS technology

WHAT IS SRS? On the Edge 3 Stretto®, independent SRS has an increased shock absorption due to shock angle. The suspension creates a 30 percent to 45 percent improvement in ride as found through vibration testing. The spring rates have been optimized for maximum obstacle climbing and overall comfort. The Stretto is the first mid-wheel drive Quantum base engineered from the ground up with true independent suspension without using an articulating beam.

IS THE EDGE 3 STRETTO® CODED?

Yes. It is coded under Group 3: Standard (K0848/K0849), Single Power (K0856) and Multiple Power (K0861).

What is the maximum speed? Up to 6.25 mph

How fast can iLevel drive with it? Up to 3.5 mph

How many inches of power seat elevation does iLevel have on it? Up to 12 inches

Compatible with EZ Lock? Yes

WC 19 crash tested? Available on all configurations

Available electronics? Q-Logic 3 EX, Q-Logic 3e, Q-Logic NE and Q-Logic NE+

Maximum range? Up to 15.9 miles with 40AH batteries; up to 20.3 miles with 55 AH batteries

Base weight? 126.4 lbs.

Base width? 21.75" with 14" drive wheels; 20.47" with 12.5" drive wheels

Weight capacity? 300 lbs.

