

**Q6 EDGE 3 STRETTO™  
CONFIGURATION WEIGHTS:  
FULLY LOADED**

CONFIGURATION	WEIGHT
iLevel, tilt, recline, AFP, 22NF batteries, 12.5" flat free tires, Redesigned back	<b>408 lbs</b>
iLevel, tilt, recline, AFP, 22NF batteries, 12.5" flat free tires, Redesigned back, Q-Back pan	<b>406 lbs</b>
iLevel, tilt, recline, AFP, 40AH batteries, 12.5" flat free tires, redesigned back	<b>392 lbs</b>
iLevel, tilt, recline, AFP, 40AH batteries, 12.5" flat free tires, redesigned back, Q-back pan	<b>390 lbs</b>
iLevel, tilt, recline, AFP, 22NF, 14" flat free tires, Redesigned back	<b>413 lbs</b>
iLevel, tilt, recline, AFP, 22NF batteries, 14" flat free tires, redesigned back, Q-Back pan	<b>411 lbs</b>
iLevel, tilt, recline, AFP, 40AH batteries, 14" flat free tires, redesigned back	<b>397 lbs</b>
iLevel, tilt, recline, AFP 40AH batteries, 14" flat free tires, redesigned back, Q-Back pan	<b>395 lbs</b>

**Q6 EDGE  
3 STRETTO™  
CONFIGURATION  
WEIGHTS: TILT,  
RECLINE, AFP**

CONFIGURATION	WEIGHT
Tilt, recline, AFP, 22NF batteries, 14" flat free tires, redesigned back	<b>363 lbs</b>
Tilt, recline, AFP, 22NF batteries, 14" flat free tires, Q-Back	<b>361 lbs</b>
Tilt, recline, AFP, 40AH batteries, 14" flat free tires, redesigned back	<b>347 lbs</b>
Tilt, recline, AFP, 40AH batteries, 14" flat free tires, Q-Back	<b>345 lbs</b>

*Note:* All configuration weights are taken with the headrest, seat cushion, and back cushion removed.

**Q6 EDGE 3  
Stretto™**

THE ULTIMATE PEDIATRIC POWER BASE



**COLORS\***



- BLUE JEAN
- CHERRY BOMB
- GREEN MACHINE
- LEMON CRUSH
- RASPBERRY BERET
- EVOLUTION ORANGE
- SUGAR PLUM



866-800-2002 • [quantumrehab.com](http://quantumrehab.com) • [lifeatlevel.com](http://lifeatlevel.com)



*Every other chair  
is at least  
2.5" wider!*

\* Visit [quantumrehab.com](http://quantumrehab.com) to see all 17 colors.

## EDGE 3 STRETTO™ FAQs

### WHY WAS THE EDGE 3 STRETTO™ CREATED?

Clinicians and providers have been asking Quantum for a true pediatric power base. While existing Quantum bases were able to accommodate pediatric seating, the Stretto has a narrower width. The Stretto has an overall width of 20.75" with 12.5" drive wheels, making it the narrowest power base in America. It is a milestone for pediatric power!

### WHAT IS SRS?

SRS (*Smooth Ride Suspension*) is a superior power chair suspension developed after two years of research and consumer feedback.

### Fun Fact...

Stretto is the Italian word for *narrow*.

Edge 3 Stretto™ featuring independent SRS technology

**WHAT IS SRS?** On the Edge 3 Stretto™, independent SRS has an increased shock absorption due to shock angle. The suspension creates a 30 percent to 45 percent improvement in ride as found through vibration testing. The spring rates have been optimized for maximum obstacle climbing and overall comfort. The Edge 3 Stretto™ is the first mid-wheel drive Quantum base engineered from the ground up with true independent suspension without using an articulating beam.

### IS THE EDGE 3 STRETTO™ CODED?

Yes. It is coded under Group 3: Standard (K0848), Single Power (K0856), and Multiple Power (K0861). Captain's seat (K0849) is pending PDAC approval.

**What is the maximum speed?** Up to 6.25 mph

**How fast can iLevel drive with it?** Up to 3.5 mph

**How many inches of power seat elevation does iLevel have on it?** Up to 12 inches

**Compatible with EZ Lock?** Shortly after launch

**WC 19 crash tested?** Available on all configurations

**Available electronics?** Q-Logic 3, Q-Logic NE, and Q-Logic NE+

**Maximum range?** Up to 15.9 miles with 40AH batteries; up to 20.3 miles with 55 AH batteries

**Base weight?** 126.4 lbs.

**Base width?** 21.75" with 14" drive wheels; 20.75" with 12.5" drive wheels

**Weight capacity?** 300 lbs.

