



Q6 EDGE® Z  
shown in *Back in Black*  
with iLevel® technology

# Q6 edge® Z

The Q6 Edge® Z utilizes powerful 4-pole motors for maximum torque and Mid-Wheel 6® design for aggressive outdoor performance and precision indoor maneuverability. The Q6 Edge® Z accepts a broad range of exclusive seating and electronics options for unmatched adaptability and rehab capability, and is available with our iLevel® technology.

**QUANTUM®**

## Specifications

<b>DRIVE WHEELS</b>	Standard: 14" Run flat Optional: 14" Pneumatic 14" Knobby
<b>CASTER WHEELS</b>	<b>Front:</b> 6" Solid <b>Rear:</b> 6" Solid
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 6 mph (3.5 mph with iLevel®)
<b>GROUND CLEARANCE</b>	3" (base frame)
<b>TURNING RADIUS</b> <sup>2</sup>	19.75"
<b>OVERALL LENGTH</b> <sup>2</sup>	35.6"
<b>BASE WIDTH</b>	26.25"
<b>SEATING SIZES</b>	
<b>SYNERGY®</b>	<b>Width:</b> 10-24" <b>Depth:</b> 10-24"
<b>TRU-BALANCE® 3 STATIC</b>	<b>Width:</b> 12-24" <b>Depth:</b> 12-24"
<b>TRU-BALANCE® 3 POWER TILT</b>	<b>Width:</b> 12-24" <b>Depth:</b> 12-24"
<b>TRU-BALANCE® 3 LIFT</b>	<b>Width:</b> 12-24" <b>Depth:</b> 12-24"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>SYNERGY®/TRU-BALANCE® 3 STATIC</b> <sup>5</sup>	16.5-18.5"
<b>TRU-BALANCE® 3 POWER TILT</b> <sup>5</sup>	18.5-20.5"
<b>TRU-BALANCE® 3 iLEVEL® &amp; TILT</b> <sup>5</sup>	18.5-20.5"
<b>BATTERY SIZE</b> <sup>3</sup>	
	Group 24 (2)
<b>BATTERY WEIGHT</b> <sup>7</sup>	
	52 lbs. each
<b>AVAILABLE ELECTRONICS</b> <sup>8</sup>	
	120A Q-Logic 3 75A Q-Logic NE+ 75A Q-Logic NE
<b>BATTERY CHARGER</b>	
	8A Off-board
<b>MOTOR PACKAGES</b>	
	4-Pole
<b>WEIGHT CAPACITY</b>	
	300 lbs.
<b>BASE WEIGHT</b>	
	143.5 lbs.

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Part of Community Use Package with 6 mph motors. Increases width by 1.5".
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Depths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Seat backs for units without TRU-Balance® 3 seating are only available in Back in Black.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability

## Features

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® 3 Power Positioning Systems
- Low impact OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters and semi-independent rear caster beam for enhanced performance over more varied terrain
- iLevel® option, featuring patented Extra Stability Technology®, enables up to 12" of seat elevation while at walking speed (up to 3.5 mph)
- Standard switch-operated LED fender lights

## Options

- iLevel® power adjustable seat height
- TRU-Balance® 3 Seating
- Q-Logic 3 Drive Controls
- Accu-Trac Advanced Tracking Technology<sup>9</sup>
- Trilogy vent tray<sup>9</sup>
- Compact vent tray<sup>9</sup>
- 4" wide flat-free or pneumatic knobby extended traction tires<sup>9,4</sup>

## Colors

Available on Q-Logic 3 controllers and units with TRU-Balance® 3 (base and back)<sup>6</sup>



Pretty in Pink



Green Machine



Lemon Crush



Deep Purple



Evolution Orange



Electric Blue



Red Corvette



Midnight Blue



Grey Street



Back in Black



Night Train

**QUANTUM®**

#1 FOR REHAB POWER

(US) 866-800-2002

| (Canada) 888-570-1113

| quantumrehab.com

