

NEW!

QUANTUM®

Remote Stop Switch

Designed to be used as a training tool, the Quantum® Remote Stop Switch helps new power chair users familiarize themselves with their seating system and to navigate confidently.

The Quantum® Remote Stop Switch consists of a receiver that mounts to the power chair and a remote fob that, when activated by a caregiver, stops the chair.

Available on the following units...

- STRETTO
- EDGE® 2
- EDGE® 3
- 4FRONT® 2
- 4FRONT® 2 HD
- R-TRAK
- EDGE® HD
- Q1450
- J6

FDA Class II
Medical Device



RECEIVER

REMOTE FOB



POINTS OF INTEREST

- The remote fob operates up to 20 feet behind the power chair
- The fob and receiver are paired prior to shipping; additional fobs can be paired to allow multiple fobs to be paired to one receiver
- When activated, the power chair stops and is put into a *drive inhibit* until the fob is pressed again, resetting the system into a normal drive status
- Available on most Quantum® Power Chairs
- Retrofittable in the field with Q-Logic electronics
- Basic operating instructions can be quickly accessed by scanning the QR code on the back of the fob

QUANTUM®

quantumrehab.com • 866.800.2002



Quantum® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

032025