



RIVAL®  
shown in *Evolution Orange*  
with TRU-Balance® 3  
Power Positioning

# RIVAL®

The Rival® rear-wheel drive power chair features innovative design that provides exceptional outdoor performance and handling while delivering outstanding tight space maneuverability.

The Rival accepts a complete range of rehab seating and electronics options.

QUANTUM®

# Specifications

<b>DRIVE WHEELS</b>	<b>Standard:</b> 14" Solid <b>Optional:</b> 14" Pneumatic; 14" Knobby
<b>CASTER WHEELS</b>	<b>Front:</b> 9" Solid <b>Optional:</b> 9" Pneumatic
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED</b> <sup>4</sup>	Up to 6 mph (Optional: 8 mph)
<b>GROUND CLEARANCE</b>	3.25" (from motor)
<b>TURNING RADIUS</b> <sup>5</sup>	25.5"
<b>OVERALL LENGTH</b> <sup>5</sup>	36"
<b>BASE WIDTH</b>	25"
<b>SEATING SIZES</b>	
<b>Synergy® / Static Seat</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 12-22" <sup>3</sup> <b>Depth:</b> 12-22" <sup>3</sup>
<b>TRU-Balance® 3 Lift &amp; Tilt</b>	<b>Width:</b> 12-22" <sup>3</sup> <b>Depth:</b> 12-22" <sup>3</sup>
<b>TRU-Balance® Flex Seat</b>	<b>Width:</b> 15-22" <b>Depth:</b> 15-20"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>Synergy®/Static Seat</b> <sup>8</sup>	Static: 16.5 - 18.1" Synergy: 15.875 - 17.875"
<b>TRU-Balance® 3 Power Tilt</b> <sup>8</sup>	18-19.6"
<b>TRU-Balance® 3 Lift &amp; Tilt</b> <sup>8</sup>	18.5-20.1"
<b>TRU-Balance® Flex Seat</b>	16.5"
<b>BATTERY SIZE</b> <sup>6</sup>	Group 24
<b>BATTERY WEIGHT</b> <sup>9</sup>	51 lbs.
<b>AVAILABLE ELECTRONICS</b> <sup>3</sup>	120A Q-Logic 3 100A Q-Logic 3 75A Q-Logic NE+ 75A Q-Logic NE
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	4-Pole
<b>WEIGHT CAPACITY</b> <sup>7</sup>	300 lbs.
<b>BASE WEIGHT</b>	129 lbs.
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)

- 1) Standard
- 2) Optional
- 3) All Quantum Rehab® units can be ordered with specialty controls
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems.
- 7) Weight capacity is dictated by the selection of seating.
- 8) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 9) Battery weight may vary +/- 3 lbs. based on manufacturer.

# Features

- Compact size for outstanding maneuverability in tight spaces
- Maximum speed up to 8 mph<sup>2</sup>
- High-performance, high-efficiency motor package
- Accepts a wide range of seating and positioning systems
- Sliding battery tray for easy access
- Two suspension packages

# Options

- TRU-Balance® 3 Power Tilt
- TRU-Balance® 3 Power Recline
- TRU-Balance® 3 (10") Lift
- Manual and Power Elevating Leg Rests
- TRU-Comfort 2 Seating
- Q-Logic 2 Controller
- Full Lighting Package
- Centermount Adjustable Foot Platform
- Accu-Trac Advanced Tracking Technology
- Power Articulating Foot Platform

# Colors



Pretty in Pink



Green Machine



Deep Purple



Evolution Orange



Electric Blue



Red Corvette



Grey Street



Midnight Blue



Night Train



Back in Black