

REHAB POWER CHAIR SOLUTIONS



RIVAL®
shown in *Orange Crush*
with TRU-Balance® 3
Power Positioning

RIVAL®

The Rival® rear-wheel drive power chair features innovative design that provides exceptional outdoor performance and handling while delivering outstanding tight space maneuverability.

The Rival accepts a complete range of rehab seating and electronics options.

QUANTUM®

Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic 14" Knobby
CASTER WHEELS	Front: 9" Pneumatic
ANTI-TIP WHEELS	N/A
MAXIMUM SPEED ⁴	Up to 6 mph (Optional: 8 mph)
GROUND CLEARANCE	3.25" (frame base)
TURNING RADIUS ⁵	35.75"
OVERALL LENGTH	36" without foot riggings
BASE WIDTH	25"

SEATING SIZES	
Synergy® / Static Seat	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Power Tilt	Width: 12-22" ³ Depth: 12-22" ³
TRU-Balance® 3 Lift & Tilt	Width: 12-22" ³ Depth: 12-22" ³

SEAT-TO-FLOOR HEIGHTS	
Synergy®/Static Seat ⁸	16.5-18.1"
TRU-Balance® Flex Seat	16.5-18.1"
TRU-Balance® 3 Power Tilt ⁸	18-19.6"
TRU-Balance® 3 Lift & Tilt ⁸	18.5-20.1"

MANUAL TILT/RECLINE ¹⁰	Yes
BATTERY SIZE ⁶	Group 24
BATTERY WEIGHT ¹⁴	52 lbs.
AVAILABLE ELECTRONICS ¹²	75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-Pole 6 mph ¹ or 8 mph ²
WEIGHT CAPACITY ⁷	300 lbs.
BASE WEIGHT	129 lbs.
MEDICARE CODE COVERAGE	Group 3: Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)

- 1) Standard
- 2) Optional
- 3) Requires HD Suspension
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems.
- 7) Weight capacity is dictated by the selection of seating.
- 8) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 9) 250 lbs. weight capacity with TRU-Balance® Systems
- 10) The Synergy® Manual Tilt weight capacity is 75 lbs. and maximum dimensions are 14" W x 14" D. Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 11) Requires Group 34 batteries
- 12) All Quantum Rehab® units can be ordered with specialty controls
- 13) Required with high-speed motors
- 14) Battery weight may vary +/- 3 lbs. based on manufacturer.

Features

- Compact size for outstanding maneuverability in tight spaces
- Maximum speed up to 8 mph²
- High-performance, high-efficiency motor package
- Accepts a wide range of seating and positioning systems
- Sliding battery tray for easy access
- Two suspension packages

Options

- TRU-Balance® 3 Power Tilt
- TRU-Balance® 3 Power Recline
- TRU-Balance® 3 (10") Lift
- Manual and Power Elevating Leg Rests
- TRU-Comfort 2 Seating
- Q-Logic 2 Controller
- Full Lighting Package
- Centermount Adjustable Foot Platform
- Accu-Trac Advanced Tracking Technology
- Power Articulating Foot Platform

Colors

