









	TRU-Comfort 2 SPP	Essence SPP	Solution® SPP	Spectrum® Foam/Gel SPP	Premiere P	Glacial SP	Zen SP	Simplicity G
DESCRIPTION¹	 Multi-layered foam Deep contour, CoolCore® cushion cover has "pockets" on sides	 Low profile Tru-Comfort 2, lateral and pre-ischial shelf Contoured, HD molded foam base, coccyx and sacral relief	 Moderate contour base foam, viscoelastic foam topper with twin cell aqueous viscous gel insert in ischial well	 Low contour base foam with sacral relief, blue XXSoft FRG foam topper (sizes 22 wide and wider, green MedSoft FRG foam topper) with twin cell ultra-lightweight viscous gel insert in ischial well.	 Moderate contour base foam with high density positioning wedges built in for offloading, gray XSoft FRG foam topper.	 New SIM G base, fluid resistant overlay, top and sides TPU heat-sealed, low profile, lightweight cooling gel topper	 Low profile, lightweight base foam, blue XXSoft FRG foam topper.	 High-density foam Stealth-Tek cover standard More contour to sacral/coccyx relief well
TOP FOAM	ViscoElastic	ViscoElastic	ViscoElastic	FRG top foam of Visco	FRG top foam or Visco	N/A	FRG top foam or Visco	N/A
WEIGHT CAPACITY	14-15" W: 175 lbs. 16-21" W: 350 lbs. 22+" W: 450 lbs.	10-15" W: 175 lbs. 16-21" W: 350 lbs. 22" W: 450 lbs.	10-15" W: 175 lbs. 16-20" W: 350 lbs. 22+" W: 500 lbs.	14" W: 175 lbs. Foam & Gel 16-20" W: 350 lbs. Foam & Gel 22+" W: 650 lbs. Foam; 500 lbs. Gel	14-15" W: 175 lbs. 16-21" W: 350 lbs. 22" W: 500 lbs.	12-15" W: 175 lbs. 16-21" W: 350 lbs. 22+" W: 400 lbs.	12-15" W: 175 lbs. 16-21" W: 350 lbs. 22+": 400 lbs.	12-15" W: 175 lbs. 16-21" W: 350 lbs. 22+" W: 400 lbs.
WEIGHT⁵	3.9 lbs.	5.5 lbs.	4.9 lbs.	Foam: 3.75 lbs. Gel: 4.3 lbs.	3.55 lbs.	4.3 lbs.	2.5 lbs.	2.2 lbs.
HEIGHT³,⁴ (AT HIGHEST POINT OF CONTOUR)	4.9"	3.7"	3"	3"	3.7"	2.9"	2.9"	2.9"
DIMENSIONS²,³,⁴	14x14" up to 24x24" Odd sizes in 15-21" W SOF available	10x10" up to 22x24" Odd sizes in 15-21" W SOF available	10x10" up to 24x24" Odd sizes in 15-19" W SOF available	Foam: 14x18" up to 32x32" Gel: 14x18" up to 24x26" Odd sizes in 17" and 19" W SOF available	14x13" up to 22x22" Odd sizes up to 21" W SOF available	12x12" up to 24x24" Odd sizes in 15-21" W No SOF	12x12" up to 24x24" Odd sizes in 15-21" W SOF available	14x16" up to 24x24" Odd sizes in 15-21" W No SOF
MSRP	\$441 - \$551	\$441 - \$656	\$414 - \$475	\$507 - \$585 (\$651 for Spectrum Foam SPP larger than 24")	\$358 - \$524	\$364 - \$485	\$348 - \$468	\$122 - \$182
HCPCS	E2607/E2608	E2607/E2608	E2607/E2608	E2607/E2608	E2605/E2606	E2603/E2604	E2603/E2604	E2601/E2602
ALLOWANCES⁶	\$340.87 / \$409.35	\$340.87 / \$409.35	\$340.87 / \$409.35	\$340.87 / \$409.35	\$316.55/ \$493.84	\$178.26 / \$221.59	\$178.26 / \$221.59	\$71.92 / \$140.41

- 1) CoolCore® is standard on all cushions except the Simplicity G.
2) For additional dimensions or special order sizes, contact sales for assistance. Special order sizes available for all cushions except the Glacial SP and the Simplicity G.
3) Due to manufacturing tolerances and continued product improvement the Length and Depth dimensions are subject to a variance of +or- 5%
4) Due to manufacturing tolerances and continued product improvement any Height dimension is subject to a variance of +or- 15%
5) Each weight corresponds to an 18"x18" cushion.
6) Allowances listed are estimations for cushions billed in conjunction with Group 3 power chair bases using the KU modifier and may vary by state.

QUANTUM®
#1 FOR REHAB POWER

Rev B/November 2024

