



	Simplicity G	Glacial SP	Zen SP	Premiere P	Spectrum® Foam/Gel SPP	Solution® SPP	TRU-Comfort 2 SPP	Essence SPP
DESCRIPTION¹	 High-density foam Stealth-Tek cover standard More contour to sacral/coccyx relief well	 New SIM G base, fluid resistant overlay, top and sides TPU heat-sealed, low profile, lightweight cooling gel topper	 New SIM G base, SunMate FRG XXSoft (blue topper), low profile, lightweight	 Gray X-Soft SunMate FRG foam topper, high density positioning wedges built in for offloading.	 SunMate FRG foam top (SunMate FRG XXSoft (blue); Medsoft (green)) with multi-layered foam insert or twin cell ultra-lightweight viscous gel insert in ischial well, sacral relief	 Standard viscoelastic foam top with twin cell aqueous viscous gel (heavier gel) insert in ischial well, mid-level contour	 Multi-layered foam Deep contour, CoolCore® cushion cover has "pockets" on sides	 Low profile Tru-Comfort 2, lateral and pre-ischial shelf Contoured, HD molded foam base, coccyx and sacral relief
WEIGHT CAPACITY	400 lbs.	400 lbs.	400 lbs.	500 lbs.	Foam: 650 lbs. Gel: 550 lbs.	500 lbs.	450 lbs.	450 lbs.
OVERALL WEIGHT	2 lbs.	3.5 lbs.	2.5 lbs.	3.5 lbs.	Foam: 3 lbs. Gel: 3.5 lbs.	5.4 lbs.	4 lbs.	3.5 - 5 lbs.
HEIGHT	2"	2.5"	3"	3.5"	3"	3"	4.5"	3.5"
DIMENSIONS²	10x10" up to 24x24" Odd sizes up to 21" SOF available	10x10" up to 24x24" Odd sizes from 15" up to 21" NO SOF	10x10" up to 24x24" Odd sizes available SOF available	14x13" up to 22x22" Odd sizes available SOF available	Foam: 14x18" up to 32x28" Gel: 14x18" up to 24x26" Odd sizes in 17" and 19" SOF available	10x10" up to 24x24" Odd sizes in 17" and 19" SOF available	14x14" up to 24x24" Odd sizes in 17" and 19" SOF available	15x15" up to 22x24" Odd sizes in 15" SOF available
MSRP	\$110 - \$165	\$330 - \$440	\$315 - \$425	\$325 - \$475	\$460 - \$530 (\$590 for Spectrum Foam SPP larger than 24")	\$375 - \$430	\$400 - \$500	\$440 - \$595
HCPCS	E2601/E2602	E2603/E2604	E2603/E2604	E2605/E2606	E2607/E2608	E2607/E2608	E2607/E2608	E2607/E2608
ALLOWANCES³	\$61.24 / \$119.55	\$151.77 / \$188.66	\$151.77 / \$188.66	\$269.53 / \$420.48	\$290.23 / \$348.54	\$290.23 / \$348.54	\$290.23 / \$348.54	\$290.23 / \$348.54

1) CoolCore® is standard on all cushions except the Simplicity G.
 2) Special order sizes available for all cushions except the Glacial SP. Contact sales for additional assistance.
 3) Allowances listed are estimations for cushions billed in conjunction with Group 3 power chair bases using the KU modifier and may vary by state.