

TRU BALANCE® 4

POWER POSITIONING SYSTEMS
WITH ANTERIOR TILT


TRU-Balance® 4
with Anterior Tilt
shown on Stretto



CHEST BAR



KNEE BLOCKS



DUAL-ACTUATOR AFP

FDA Class II
Medical Device*



TRU-Balance® 4
with Anterior Tilt
shown on Stretto

FDA Class II
Medical Device*

Anterior tilt options:

10° | 20° | 30°

Up to 8 customizable memory
seating positions

Articulating Foot Platform lowers
to assist with transfers

Available on the Stretto,
Edge® 3 and 4Front® 2.

Specifications:

Seating sizes: 16"W x 14"D - 22"W x 22"D

Anterior tilt range: 0-30° with ranges of 0-10°, 0-20° and 0-30°

Posterior tilt range: 0-50°

Recline range: 85-165°

AFP articulation range: 0-7.25"

AFP elevation range: 85-165°

Drive speeds in TRU-Balance 4 anterior tilt:

Up to 10° of anterior tilt, up to full lift: 2 mph

10-20° of anterior tilt, up to 6" of lift: 0.5 mph

10-20° of anterior tilt, over 6" of lift: drive lockout

Over 20° of anterior tilt: drive lockout

*Quantum Rehab® power chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

QUANTUM®
#1 FOR REHAB POWER

866-800-2002 • quantumrehab.com • lifeatilevel.com •

