

# TRU BALANCE<sup>®</sup>

## FLEX



FDA Class II  
Medical Device

TRU-Balance<sup>®</sup> Flex Seating is the no-charge static seat option available on the Stretto, Edge<sup>®</sup> 3, Q6 Edge<sup>®</sup> 2.0, J4 and R-TRAK power chairs.



FDA Class II  
Medical Device

  
**QUANTUM<sup>®</sup>**  
#1 FOR REHAB POWER



## Specifications

<b>WEIGHT CAPACITY</b>	300 lbs.
<b>SEAT DIMENSIONS</b>	Width: 15-22" Depth: 15-20"
<b>FLIP-FORWARD BACK CANE</b>	Angle can be set from 85° to 115° (in 5° increments)
<b>BACK OPTIONS</b>	TRU-Comfort 2, sling back, Synergy® back, and Sport back
<b>SEAT-TO-FLOOR HEIGHT</b>	16.5" - 18.25"
<b>ARMS</b>	Single-post armrests
<b>ARM HEIGHT RANGES</b>	9"-14" (adjustable in 1" increments)



## Features

- Easily adjustable, this versatile design can grow with your client
- Same ease of width and depth adjustment as TRU-Balance® 3
- No-charge, TRU-Balance® Flex, height-adjustable, single-post armrests with clothing guard standard
- Accepts all rail-mounted TRU-Balance 3 accessories
- Angle-adjustable fold-down back canes
- Single seat pan covers all sizes
- Height and angle-adjustable seat towers
- 4-post system on Edge 3 Stretto®, Edge® 3 and Q6 Edge® 2.0
- WC19 crash-test approved occupied transit kit option on TRU-Balance® Flex seating with Stretto, Edge® 3, Q6 Edge® 2.0 and J4



## Optional Components

- TRU-Balance® 3 style back canes with the following options:
  - *Standard angle adjustable hatchet assembly*
  - *Quick fold-down cane assembly*
  - *Endo and Meso offset brackets*
  - *Optional height from 14"-25"*
- Pediatric Flex Tilt available in 10-14" W x 10-14" D (Stretto only)
- Accepts aftermarket back rests
- Flex drop-in "T" armrests (desk or full length), flip-back cane mounted armrests, and dual-post armrests
- Cane mounted flip-back arms – available on TRU-Balance canes
- TRU-Balance #8 swing-away legrests
- Elevating legrests
- Center Mount Foot Platform

**QUANTUM®**

#1 FOR REHAB POWER

Follow us



• 866.800.2002 • [quantumrehab.com](http://quantumrehab.com)