

TRU BALANCE[®]

FLEX



TRU-Balance[®] Flex Seating is a low-cost static seat option available on the Edge 3, Q6 Edge[®] 2.0 and Rival[®].

TRU-Balance[®] Flex Seating shown on Edge 3 power base in Grey Street with NE+ controller





Specifications

WEIGHT CAPACITY	300 lbs.
SEAT DIMENSIONS	Width: 15-22" Depth: 15-20"
FLIP-FORWARD BACK CANE	Angle can be set from 85° to 115° (in 5° increments)
BACK OPTIONS	TRU-Comfort 2, sling back, Synergy® back, and Sport back
SEAT-TO-FLOOR HEIGHT	16.5" on Edge 3, Q6 Edge® 2.0, and Rival®
ARMS	Single-post armrests
ARM HEIGHT RANGES	9"-14" (adjustable in 1" increments)



Optional Components

- TRU-Balance® 3 style back canes with the following options:
 - *Standard angle adjustable hatchet assembly*
 - *Quick fold-down cane assembly*
 - *Endo and Meso offset brackets*
 - *Optional height from 14"-25"*
- Accepts aftermarket back rests
- Flex drop-in "T" armrests (desk or full length), flip-back cane mounted armrests, and dual-post armrests
- Cane mounted flip-back arms – available on TRU-Balance canes
- TRU-Balance #8 swing-away legrests
- Elevating legrests
- Center Mount Foot Platform



Features

- Easily adjustable, this versatile design can grow with your client
- Same ease of width and depth adjustment as TRU-Balance® 3
- No-charge, TRU-Balance® Flex, height-adjustable, single-post armrests with clothing guard standard
- Accepts all rail-mounted TRU-Balance 3 accessories
- Angle-adjustable fold-down back canes
- Single seat pan covers all sizes
- Height and angle-adjustable seat towers