



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 2.0 Order Form with Synergy Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

☐ Quantum Q6 Edge 2.0 3S-SS\$8,704
Part: Q6 EDGE 2.0 3S-SS. HCPCS: K0848
Solid seat pan standard with no power option.

☐ Quantum Q6 Edge 2.0 3SP-SS\$8,399
Part: Q6 EDGE 2.0 3SP-SS. HCPCS: K0856
Single power option for Power Elevating Seat.

2. DRIVE MOTORS

REQUIRED.

☐ 6 mph Drive Motors Standard
Part: MOT147433

☐ Accu-Trac Motor Technology \$1,387
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

Account Number: _____ Date: _____ Quote Name: _____

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

- | | | | |
|---|--|--|--|
| ☐  Red Corvette
Part: DGN153201-PRD1 | ☐  Night Train
Part: DGN153201-PBK3 | ☐  Deep Purple
Part: DGN153201-PPR5 | ☐  Lemon Crush
Part: DGN153201-PYL4 |
| ☐  Electric Blue
Part: DGN153201-PBU9 | ☐  Pretty in Pink
Part: DGN153201-PPK3 | ☐  Grey Street
Part: DGN153201-PSL3 | ☐  Green Machine
Part: DGN153201-PGN11 |
| ☐  Back in Black
Part: DGN153201-PBK1 | ☐  Evolution Orange
Part: DGN153201-POR7 | ☐  Midnight Blue
Part: DGN153201-PBU10 | ☐  Military Green
Part: DGN153201-PGN8 |

4. POWER ADJUSTABLE SEAT HEIGHT

Optional. The seat to floor height is 15.875" with 6" stroke and 16.875" with 7" stroke. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit. *Please verify the appropriate billing code to be used for this item with all third-party payors.

- | | |
|---|---|
| ☐ Power Elevating Seat with 7" Stroke \$1,990
Part: XRF000131. HCPCS: E2298/E2300^ | ☐ Power Elevating Seat with 6" Stroke \$1,990
Part: XRF000132. HCPCS: E2298/E2300^ |
|---|---|

5. BASE MOUNTED LIGHTS

Optional. Select only one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

5.1a. Accent Lights

Accent lights operate by pressing included single switch button.


- ☐ **Switch Accent Lights (without iLevel 12" Lift System)**..... \$268
Part: XRF000177

5.1b. Base Mounted Light

- ☐ **Base Mounted, LED, Full Light Package** \$1,066
Part: XRF000129
Requires QLNES 10 Key, Q-Logic 3 Joystick for Full Light Package, or Q-Logic 3e Joystick.



Base Mounted Lights



Accent Lights

6. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

- ☐ **Quantum Backup Camera** \$375
Part: MEC1910365

1. Optional Camera Mounting Bracket

- | | |
|---|--|
| ☐ Backup Camera Headrest Bracket \$21
Part: MEC1910373
Works with .5" x .5" vertical supports only. | ☐ Backup Camera Track Mount Bracket \$28
Part: MEC2104227 |
|---|--|

2. Optional Display Mounting Bracket

- ☐ **Backup Camera Q-Logic Joystick Mount Bracket** \$107
Part: MEC2106123



Backup Camera

7. WHEEL OPTIONS

REQUIRED. Select a type of drive wheels and accent color if applicable. Knobby tires increase overall base width by 1.5".

7.1. Drive Wheels

- | | |
|---|---|
| ☐ 3"x14" Flat Free Black Tires (No Accent Color) Standard
Part: PTOINDV3525 (pair) | ☐ 3"x14" Pneumatic Black Tires (No Accent Color) No Charge
Part: WHL150034 x2 |
| ☐ 3"x14" Flat Free Black Tires for Wheel Accents No Charge
Part: WHL1611297 x2 | ☐ 3"x14" Pneumatic Black Tires for Wheel Accents No Charge
Part: WHL1611296 x2 |
| ☐ 4"x14" Knobby, Flat Free Black Tires \$592
Part: XRF000368 | ☐ 4"x14" Knobby, Pneumatic Black Tires \$592
Part: XRF000367 |

7.2. Drive Wheel Accent Color

Optional. Not available with knobby tires.

- | | | |
|--|---|---|
| ☐ Back in Black Accent Color \$64
Part: XRF000179 | ☐ Deep Purple Accent Color \$64
Part: XRF000184 | ☐ Electric Blue Accent Color \$64
Part: XRF000188 |
| ☐ Red Corvette Accent Color \$64
Part: XRF000180 | ☐ Pretty in Pink Accent Color \$64
Part: XRF000185 | ☐ Lemon Crush Accent Color \$64
Part: XRF000189 |
| ☐ Night Train Accent Color \$64
Part: XRF000181 | ☐ Evolution Orange Accent Color \$64
Part: XRF000186 | ☐ Green Machine Accent Color \$64
Part: XRF000190 |
| ☐ Grey Street Accent Color \$64
Part: XRF000182 | ☐ Midnight Blue Accent Color \$64
Part: XRF000187 | ☐ Military Green Accent Color \$64
Part: XRF000191 |

8. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

- | | |
|--|---|
| ☐ NF-22, Gel Cell, Interceptor 2 Battery \$362 ea.
Part: BAT2004376. HCPCS: E2361 | ☐ NF-22, Gel Cell, MK Battery \$441 ea.
Part: BATGEL1008. HCPCS: E2361 |
| ☐ NF-22, AGM, MK Battery \$441 ea.
Part: BAT121346. HCPCS: E2361 | |

9. SEAT TOWERS**Optional.** Not available with Power Adjustable Seat Height.

- ☐ **Extended Seat Towers** No Charge
 Part: MEC120301
 Increases seat to floor height.

10. ELECTRONICS**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**10.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

- ☐ **QLNES 6 Key Joystick** Standard ☐ **QLNES 10 Key Joystick** \$429
 Part: CTL2205379 Part: CTL2205380
 Controls up to 2 actuators. Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

**10.1b. Expandable Controls**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

- ☐ **Q-Logic 3 EX Controller Module** \$745 ☐ **Q-Logic 3 Expandable Controller Harness** \$467
 Part: XRF000002. HCPCS: E2377 Part: XRF000264. HCPCS: E2313

2. Q-Logic 3 Joystick

- ☐ **Q-Logic 3 EX Joystick** \$857 ☐ **Q-Logic 3 EX Joystick for Full Lights** \$857
 Part: CTL167724 Part: CTL167725
☐ **Q-Logic 3e Joystick** \$735 ☐ **No Joystick, Setup for Q-Logic**
 Part: CTL1910600 Part: XRF000010
 Can operate full light package. Not available with alternative or specialty controls. Requires selection of Enhanced Display

**10.2. Q-Logic Magnetic Charger**

Available with QLNES. Included with Q-Logic 3e Joystick.

- ☐ **Magnetic to XLR Charger Adapter with Wall Holder** \$61
 Part: ELE2103635

10.3. Q-Logic 3 Joystick Armor**REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

- | | | |
|--|---|---|
| ☐ Black, Soft Touch Standard
Part: DGN124307
Also matches Night Train shrouds. | ☐ Deep Purple \$64
Part: DGN124336-PPR5 | ☐ Electric Blue \$64
Part: DGN124336-PBU9 |
| ☐ Red Corvette \$64
Part: DGN124336-PRD1 | ☐ Pretty in Pink \$64
Part: DGN124336-PPK3 | ☐ Lemon Crush \$64
Part: DGN124336-PYL4 |
| ☐ Back in Black \$64
Part: DGN124336-PBK1 | ☐ Evolution Orange \$64
Part: DGN124336-POR7 | ☐ Green Machine \$64
Part: DGN124336-PGN11 |
| ☐ Grey Street \$64
Part: DGN124336-PSL3 | ☐ Midnight Blue \$64
Part: DGN124336-PBU10 | ☐ Military Green \$64
Part: DGN124336-PGN8 |


11. JOYSTICK HANDLES**Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

- | | |
|---|--|
| ☐ Mushroom Handle (B) \$155
Part: ACCASMB1860. HCPCS: E2323 | ☐ Extended I-Handle (G) \$174
Part: ACCASMB1857. HCPCS: E2323 |
| ☐ Chin Cup (C) \$155
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control) | ☐ Bodypoint 3" Goal Post Handle \$179
Part: PLSKNCP1057. HCPCS: E2323 |
| ☐ T-Bar (D) \$195
Part: ACCASMB1858. HCPCS: E2323 | ☐ Bodypoint 4" Goal Post Handle \$179
Part: PLSKNCP1058. HCPCS: E2323 |
| ☐ Medium Ball Handle (E) \$155
Part: ACCASMB1859. HCPCS: E2323 | ☐ Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308 |
| ☐ Soft Ball Handle (F) \$174
Part: ACCASMB1856. HCPCS: E2323 | |




12. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**12.1a. Fixed, In-line Joystick Bracket**☐ Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617**12.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

☐ Swing-away, Q-Logic 3 \$325
Part: ACC1908434. HCPCS: E1032**12.1c. Basic Swing Away Joystick Mount**☐ Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1032☐ Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1032**12.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

☐ Retractable Mount \$375
Part: XRF000003. HCPCS: E1032
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.
 ☐ Multi-Axis Mount \$129
Part: XRF000004
Allows height adjustment of joystick.
**12.2. Joystick Mount Location**
☐ Left Side Mount
Part: PTOINDV1225
 ☐ Right Side Mount
Part: PTOINDV1226
13. SYNERGY SEATING SYSTEM**Optional.** Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 22"W x 22"D.**13.1. Synergy Frame Selection****REQUIRED.**☐ Synergy Seating System \$465
Part: PTOINDV2995**13.2. Synergy Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. Omitting a backrest selection requires a provider to facilitate an aftermarket backrest that best fits the consumer's needs. The third-party backrest if applicable must be occupied transit approved per ISO/RESNA standards. When modifying a power chair, including the addition of third-party componentry, to accommodate the medical needs of the power chair occupant, a risk assessment in conformance with ISO 14971, as outlined in ISO 7176, should be performed. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>☐ WC-19 Occupied Transit Option \$295
Part: KITASMB1830**13.3. Seat Width****REQUIRED.** Maximum seat width is 22". Part number is SETWDTH plus number following # below.
☐ 10" W (#1000) ☐ 11" W (#1001) ☐ 12" W (#1002) ☐ 13" W (#1003) ☐ 14" W (#1004) ☐ 15" W (#1005) ☐ 16" W (#1006) ☐ 17" W (#1007) ☐ 18" W (#1008) ☐ 19" W (#1009) ☐ 20" W (#1010) ☐ 21" W (#1011) ☐ 22" W (#1012)
13.4. Seat Width Package**REQUIRED.** Maximum seat width is 22".
☐ 10"-14" Width Package
Part: FRMASMB18011
 ☐ 14"-18" Width Package
Part: FRMASMB18012
 ☐ 16"-20" Width Package
Part: FRMASMB18013
 ☐ 20"-24" Width Package
Part: FRMASMB18014
13.5. Seat Depth**REQUIRED.** Maximum seat depth is 22". Part number is SETDPTH plus number following # below.
☐ 10" D (#1000) ☐ 11" D (#1001) ☐ 12" D (#1002) ☐ 13" D (#1003) ☐ 14" D (#1004) ☐ 15" D (#1005) ☐ 16" D (#1006) ☐ 17" D (#1007) ☐ 18" D (#1008) ☐ 19" D (#1009) ☐ 20" D (#1010) ☐ 21" D (#1011) ☐ 22" D (#1012)
13.6. Seat Depth Package**REQUIRED.** Maximum seat depth is 22".
☐ 10"-14" Depth Package
Part: FRMASMB6981
 ☐ 14"-18" Depth Package
Part: FRMASMB17823
 ☐ 15"-20" Depth Package
Part: FRMASMB7732
 ☐ 20"-24" Depth Package
Part: FRMASMB7733

Low	Med	High
Synergy Seat		
15.875"	16.875"	17.875"
Extended Towers		
Add 3 inches to above		
Seat to Floor Height		

13.7. Synergy Seat to Floor**REQUIRED.**
☐ Lowest Seat to Floor
Part: PTOINDV3558
 ☐ Middle Seat to Floor
Part: PTOINDV3559
 ☐ Highest Seat to Floor
Part: PTOINDV3560

Synergy seating section continued on following pages.

13. SYNERGY SEATING SYSTEM CONTINUED**13.8. Stealth Seat Cushion**

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

Cushion Model ####-###XXYY	Width	Depth
Cushion Top Foam V=Visco Top Foam: 1" thick, XXSoft viscoelastic foam topper. F=FRG Top Foam: Under 22"W, blue XXSoft FRG foam topper. 22"W and over, green MedSoft FRG foam topper. N=No Top Foam	Cushion Liner U=Urethane Liner: 2Way Stretch. P=Polyurethane Liner: 4Way Stretch. N=No Liner	Cushion Cover C=CoolCore: Chemical free, 4Way stretch material that offers moisture wicking, moisture transportation, and moisture evaporation properties. M=Mesh: Superior air circulation, non-slip. T=Stealth Tek: Waterproof with stretch capability.

1. Cushion Type

- | | |
|---|---|
| ☐ TRU-Comfort 2 SPP
Reference part: XRF000116. See cushion order form.
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608 | ☐ Premiere P
Reference part: XRF000122. See cushion order form.
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606 |
| ☐ Essence SPP
Reference part: XRF000118. See cushion order form.
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608 | ☐ Glacial SP
Reference part: XRF000114. See cushion order form.
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604 |
| ☐ Solution SPP
Reference part: XRF000113. See cushion order form.
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608 | ☐ Zen SP
Reference part: XRF000115. See cushion order form.
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604 |
| ☐ Spectrum Foam SPP
Reference part: XRF000112. See cushion order form.
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608 | ☐ Simplicity G
Reference part: XRF000111. See cushion order form.
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602 |
| ☐ Spectrum Gel SPP
Reference part: XRF000117. See cushion order form.
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608 | |

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

- | | | |
|---|--|---|
| ☐ Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller | ☐ Medium Wedge Positioning Kit\$72
Part: POP1909006-M
14"W and smaller | ☐ Large Wedge Positioning Kit \$105
Part: POP1909006-L
14"W and smaller |
| ☐ Small Wedge Positioning Kit\$39
Part: POP1712340-S
15"W and larger | ☐ Medium Wedge Positioning Kit\$72
Part: POP1712340-M
15"W and larger | ☐ Large Wedge Positioning Kit \$105
Part: POP1712340-L
15"W and larger |

3. Cushion Top Foam

- | | |
|---|--|
| ☐ Visco Top Foam
Reference part: XRF000380 | ☐ No Top Foam
Reference part: XRF000382 |
| ☐ FRG Top Foam
Reference part: XRF000381 | |

4. Cushion Liner

- | | |
|---|---|
| ☐ Stealth Urethane Liner
Reference part: XRF000383 | ☐ No Liner
Reference part: XRF000385 |
| ☐ Stealth Polyurethane Liner
Reference part: XRF000384 | |

5. Cushion Cover

- | | |
|---|--|
| ☐ CoolCore Cover
Reference part: XRF000386 | ☐ Stealth Tek Cover
Reference part: XRF000388 |
| ☐ Mesh Cover
Reference part: XRF000387 | |

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

- ☐
- ABS Flat Rigidizer**
-
- Reference part: XRF000389
-
- MSRP: \$48-\$69

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner					
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

Synergy seating section continued on following pages.

13. SYNERGY SEATING SYSTEM CONTINUED**7. Additional Cushion Liner**

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

☐ Additional Stealth Urethane Liner
Reference Part: XRF000390☐ Additional Stealth Polyurethane Liner
Reference Part: XRF000391**8. Additional Cushion Cover**

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

☐ Additional CoolCore Cover
Reference Part: XRF000392☐ Additional Stealth Tek Cover
Reference Part: XRF000394☐ Additional Mesh Cover
Reference Part: XRF000393**13.9. Synergy Seating Back Cane Selection****REQUIRED.**☐ Synergy Back Cane No Charge
Part: FRMASMB15845☐ 2" Endomorph (1" per side) Back Cane \$209
Part: FRMASMB12115☐ Flip Forward Back Cane \$188
Part: FRMASMB15846
Not available with a vent tray, oxygen holder, or cane mounted and cantilever armrests.☐ 4" Endomorph (2" per side) Back Cane \$209
Part: FRMASMB12139☐ 2" Mesomorph (1" per side) Back Cane \$205
Part: FRMASMB12067**13.10. Synergy Seating Back Cane Height****REQUIRED with Synergy seating EXCEPT if a TRU-Comfort back is desired.** 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles.☐ 14" High Canes
Part: FRMASMB15832☐ 17" High Canes
Part: FRMASMB15835☐ 20" High Canes
Part: FRMASMB15838☐ 23" High Canes
Part: FRMASMB15841☐ 15" High Canes
Part: FRMASMB15833☐ 18" High Canes
Part: FRMASMB15836☐ 21" High Canes
Part: FRMASMB15839☐ 24" High Canes
Part: FRMASMB15842☐ 16" High Canes
Part: FRMASMB15834☐ 19" High Canes
Part: FRMASMB15837☐ 22" High Canes
Part: FRMASMB15840☐ 25" High Canes
Part: FRMASMB15843**13.11. Synergy Seating Back Cushion****REQUIRED.** Select one type of back.**1a. Sport Back Selection**☐ Sport Back Cushion, 21" Wide or Less \$402
Part: XRF000361☐ Sport Back Cushion, 22" Wide or More \$616
Part: XRF000362**1b. TRU-Comfort 4 Back Selection**

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

☐ TRU-Comfort Back, Back Cane Option, Less than 22" Wide \$759
Part: TC4V10W21WBC. HCPCS: E2620☐ TRU-Comfort Back, Back Cane Option, 22" Wide or Greater ... \$872
Part: TC4V22W32WBC. HCPCS: E2621**1c. Q-Back 2.0, AL Series, Active Contour**

Select a back type and the specific overall back height in subsequent section.

☐ Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL14AS. HCPCS: E2613☐ Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL18AT. HCPCS: E2613☐ Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL14AT. HCPCS: E2613☐ Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL18AX. HCPCS: E2613☐ Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL16AS. HCPCS: E2613☐ Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18AXX. HCPCS: E2613☐ Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL16AT. HCPCS: E2613☐ Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL20AS. HCPCS: E2613☐ Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL16AX. HCPCS: E2613☐ Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL20AT. HCPCS: E2613☐ Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL16AXX. HCPCS: E2613☐ Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL20AX. HCPCS: E2613☐ Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL18AS. HCPCS: E2613☐ Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20AXX. HCPCS: E2613

Synergy seating section continued on following pages.

13. SYNERGY SEATING SYSTEM CONTINUED**1d. Q-Back 2.0, AL Series, Deep Contour**

Select a back type and the specific overall back height in subsequent section.

- | | |
|--|--|
| ☐ Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL14DS. HCPCS: E2620 | ☐ Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL18DT. HCPCS: E2620 |
| ☐ Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL14DT. HCPCS: E2620 | ☐ Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL18DX. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL16DS. HCPCS: E2620 | ☐ Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18DXX. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL16DT. HCPCS: E2620 | ☐ Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL20DS. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL16DX. HCPCS: E2620 | ☐ Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL20DT. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL16DXX. HCPCS: E2620 | ☐ Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL20DX. HCPCS: E2620 |
| ☐ Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL18DS. HCPCS: E2620 | ☐ Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20DXX. HCPCS: E2620 |

1e. Q-Back 2.0, AL Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|--|
| ☐ Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL14ES. HCPCS: E2620 | ☐ Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL18EX. HCPCS: E2620 |
| ☐ Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$703
Part: QB2AL16ES. HCPCS: E2620 | ☐ Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18EXX. HCPCS: E2620 |
| ☐ Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$703
Part: QB2AL16ET. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL20ES. HCPCS: E2620 |
| ☐ Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL16EX. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL20ET. HCPCS: E2620 |
| ☐ Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL18ES. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL20EX. HCPCS: E2620 |
| ☐ Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL18ET. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20EXX. HCPCS: E2620 |

1f. Q-Back 2.0, AFT Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|---|
| ☐ AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT16AS. HCPCS: E2613 | ☐ AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT18AS. HCPCS: E2613 |
|---|---|

1g. Q-Back 2.0, AFT Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|---|
| ☐ AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$643
Part: QB2AFT16DS. HCPCS: E2620 | ☐ AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$643
Part: QB2AFT18DT. HCPCS: E2620 |
| ☐ AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$643
Part: QB2AFT16DT. HCPCS: E2620 | ☐ AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AFT20DT. HCPCS: E2620 |
| ☐ AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AFT18DS. HCPCS: E2620 | |

1h. Q-Back 2.0, AFT Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|--|
| ☐ AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline..... \$703
Part: QB2AFT16EXX. HCPCS: E2620 | |
|---|--|

1i. Omit Back Selection

- | | |
|--|--|
| ☐ Omit Back..... No Charge
Part: PTOINDV1287 | |
|--|--|

Synergy seating section continued on following pages.

13. SYNERGY SEATING SYSTEM CONTINUED**13.12. Synergy Seating Overall Back Height****REQUIRED.**

- | | |
|---|---|
| ☐ 14" Overall Back Height
Part: XRF000136 | ☐ 22" Overall Back Height
Part: XRF000144 |
| ☐ 15" Overall Back Height
Part: XRF000137 | ☐ 23" Overall Back Height
Part: XRF000145 |
| ☐ 16" Overall Back Height
Part: XRF000138 | ☐ 24" Overall Back Height
Part: XRF000146 |
| ☐ 17" Overall Back Height
Part: XRF000139 | ☐ 25" Overall Back Height
Part: XRF000147 |
| ☐ 18" Overall Back Height
Part: XRF000140 | ☐ 26" Overall Back Height
Part: SETHGHT1039 |
| ☐ 19" Overall Back Height
Part: XRF000141 | ☐ 27" Overall Back Height
Part: SETHGHT1040 |
| ☐ 20" Overall Back Height
Part: XRF000142 | ☐ 28" Overall Back Height
Part: SETHGHT1041 |
| ☐ 21" Overall Back Height
Part: XRF000143 | ☐ 29" Overall Back Height
Part: SETHGHT1042 |

		Back Width											
		10"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
Overall Height	14"	P	T	P	T	P	T	-	-	-	-	-	-
	15"	P	T	P	T	P	T	P	T	P	T	P	T
	16"	P	T	P	T	P	T	P	T	P	T	P	T
	17"	P	T	P	T	P	T	P	T	P	T	P	T
	18"	P	T	P	T	P	T	P	T	P	T	P	T
	19"	P	T	P	T	P	T	P	T	P	T	P	T
	20"	-	-	P	T	P	T	P	T	P	T	P	T
	21"	-	-	P	T	P	T	P	T	P	T	P	T
	22"	-	-	P	T	P	T	P	T	P	T	P	T
	23"	-	-	P	T	P	T	P	T	P	T	P	T
	24"	-	-	P	T	P	T	P	T	P	T	P	T
	25"	-	-	P	T	P	T	P	T	P	T	P	T
	26"	-	-	P	T	P	T	P	T	P	T	P	T
	27"	-	-	P	T	P	T	P	T	P	T	P	T
	28"	-	-	P	T	P	T	P	T	P	T	P	T
	29"	-	-	P	T	P	T	P	T	P	T	P	T

☐ Sport Back ☐ TRU-Comfort Back

Seat Back Availability Chart

13.13. Synergy Armrests**REQUIRED with Synergy seating.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. Flip Back, Cane Mounted Armrests**

May interfere with swing-away laterals.

- | | |
|--|---|
| ☐ ☐ Full Length, Flip Back Armrests..... \$63
L R Left Part: FRMASMB7038. Right part: FRMASMB7039 | ☐ ☐ Desk Length, Flip Back Armrests \$63
L R Left part: FRMASMB7040. Right part: FRMASMB7041 |
|--|---|

1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

- | | |
|--|---|
| ☐ ☐ Full Length, Two-Post Armrests..... \$161 ea.
L R Left part: XRF000051. Right part: XRF000052. HCPCS: E0973 | ☐ ☐ Desk Length, Two-Post Armrests \$161 ea.
L R Left part: XRF000053. Right part: XRF000054. HCPCS: E0973 |
|--|---|

1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

- | | |
|--|--|
| ☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H \$92 ea.
L R Left Part: XRF000245. Right part: XRF000246. HCPCS: E0973 | ☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H..... \$92 ea.
L R Left part: XRF000247. Right part: XRF000248. HCPCS: E0973 |
|--|--|

1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

- | | |
|---|--|
| ☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H..... \$92 ea.
L R Left Part: XRF000249. Right part: XRF000250. HCPCS: E0973 | ☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14"H \$92 ea.
L R Left part: XRF000251. Right part: XRF000252. HCPCS: E0973 |
|---|--|

1e. Pediatric, Removable, Height Adjustable Armrests

Only available with desk length armpads.

- | |
|---|
| ☐ ☐ Pediatric, Removable Armrests \$81 ea.
L R Left Part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973 |
|---|

1f. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

- | | |
|---|--|
| ☐ ☐ Full Length, Cantilever Armrests \$105 ea.
L R Left Part: FRMASMB16188. Right part: FRMASMB16189 | ☐ ☐ Desk Length, Cantilever Armrests..... \$105 ea.
L R Left part: FRMASMB16190. Right part: FRMASMB16191 |
|---|--|

1g. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

- | | |
|---|---|
| ☐ ☐ Full Length, Single Post Armrests \$79 ea.
L R Left Part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973 | ☐ ☐ Desk Length, Single Post Armrests \$79 ea.
L R Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973 |
|---|---|

1h. Drop-In, Height Adjustable Armrests (T-Arm)

Includes flat armpads.

- | | |
|--|--|
| ☐ ☐ Full Length, Drop-In Armrests..... No Charge
L R Left part: ARM1803253. Right part: ARM156742 | ☐ ☐ Desk Length, Drop-In Armrests No Charge
L R Left part: ARM1803254. Right part: ARM156743. |
|--|--|

13.14. Synergy Armpads**REQUIRED with armrest selection above.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

- | | |
|---|--|
| ☐ Straight Armpads..... Standard
Part: XRF000148 | ☐ Waterfall Armpads..... No Charge
Part: XRF000149 |
| ☐ No Armpad Selection No Charge
Part: PTOINDV1373 | |

13.15. Joystick Bracket Receiver**REQUIRED.**

- | | |
|---|---|
| ☐ Standard Fixed Receiver Standard
Part: FRMASMB16062 | ☐ Stealth Height Adjustable Receiver..... \$236
Part: XRF000064. HCPCS: K0108 |
|---|---|

Synergy seating section continued on following pages.

14. HEADRESTS**Optional.** Select a headrest pad and bracket.**14.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads, Single Point Adjustable**

- ☐ Stealth 8" Comfort Plus, Cool Core (CP08CCSPM)\$179
Part: CP08CCSPM. HCPCS: E0955
- ☐ Stealth 14" Comfort Plus, Cool Core (CP14CCSPM)\$179
Part: CP14CCSPM. HCPCS: E0955
- ☐ Stealth 10" Comfort Plus, Cool Core (CP10CCSPM)\$179
Part: CP10CCSPM. HCPCS: E0955
- ☐ Stealth 6", Comfort Plus, Cool Core (CP06CCSPM)\$179
Part: CP06CCSPM. HCPCS: E0955

2. Headrest Brackets

- ☐ Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033
- ☐ Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033

14.2. Flip-down Headrest Bracket

- ☐ Stealth Flip-down Headrest Bracket (FDM350)..... \$170
Part: INDPART2219

The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

14.3. Q Cam Lock**Optional.**

- ☐ Q Cam Lock..... \$51
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

15. LAP BELTS**REQUIRED.** Select one lap belt below.**15.1a. Standard Lap Belts**

- ☐ 38" Lap Belt
Part: ACCBELT1006
- ☐ 50" Lap Belt
Part: ACCBELT1000
- ☐ 60" Lap Belt
Part: ACCBELT1002
- ☐ 70" Lap Belt
Part: ACCBELT1003
- ☐ 80" Lap Belt
Part: ACCBELT1004
- ☐ 90" Lap Belt
Part: ACCBELT1001

15.1b. Optional Lap Belts

- ☐ Stealth Positioning 60" Single Pull, Push Button, Medium..... \$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q
- ☐ Stealth Positioning 60" Single Pull, Push Button, Extra Large.... \$50
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q
- ☐ Stealth Positioning 60" Single Pull, Push Button, Large \$50
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q
- ☐ 60" Sure-Lok Lap Belt (Non-retractable)..... \$120
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016

16. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

16.1a. Base Mounted Foot Platform

May interfere with seat mounted legrests. Not available with center mount foot platforms.

- ☐ Base Mounted Foot Platform Standard
Part: FRMASMB16077
- ☐ Omit base mounted foot platform No Charge
Part: PTOINDV3674
- ☐ Bell Shaped, Base Mounted Foot Platform \$252
Part: RIGASMB7120041

16.1b. Basic Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

1. Foot Platform Assembly

- ☐ Basic Center Mount..... No Charge
Part: XRF000323

2. Basic Center Mount Extension and Angle Adjustable Footplates

- ☐ Short Extension, 5"x9" Footplates \$184
Part: XRF000312. HCPCS: K0040
- ☐ Long Extension, 5"x9" Footplates \$184
Part: XRF000314. HCPCS: K0040
- ☐ Medium Extension, 5"x9" Footplates \$184
Part: XRF000313. HCPCS: K0040

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

- ☐ Small Calf Pads, Basic Center Mount.....\$161
Part: RIG2104369
- ☐ Large Calf Pads, Basic Center Mount\$161
Part: RIG2202997
Only available with Long Length Extension.

Legrests section continued on following page.

16. LEGRESTS CONTINUED**4. Basic Center Mount Additional Accessories**

- ☐ Basic Center Mount Footplate Connector \$16
Part: RIG2104240

16.1c. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform Assembly

- ☐ Power Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

2. AFP Lower Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility. Do not select with Basic Center Mount.

- ☐ Medium Extension, 5"x8" Footplates
Part: RIG2404016
- ☐ Long Extension, 5"x8" Footplates
Part: RIG2404017
- ☐ Medium Extension, 5"x11" Footplates
Part: RIG2404018
- ☐ Long Extension, 5"x11" Footplates
Part: RIG2404019

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

- ☐ Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018
- ☐ Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019

3b. Center Mount Calf Panel

- ☐ Small Calf Panel, 18"Wx8"H \$188
Part: FRMASMB7110012
- ☐ Large Calf Panel, 24"Wx8"H \$242
Part: RIG123825
Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

- ☐ Width Extension for 8" \$70
Part: RIG2404919
- ☐ Width Extension for 11" \$70
Part: RIG2404920
- ☐ Outside Edge \$70
Part: RIG2404921
Not available with width extensions.
- ☐ Outside Edge for Width Extension \$70
Part: RIG2404922
Available with footplates and width extension.
- ☐ Front and Rear Cap \$16
Part: RIG2404923

5. Additional Accessories

- ☐ Small Foam Calf Pad Cover (pair) \$67
Part: ST-CSCS0775-Q
- ☐ Large Foam Calf Pad Cover (pair) \$83
Part: ST-CSCS0985-Q
- ☐ Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-Q
- ☐ Large Foam Panel Cover (pair) \$282
Part: ST-CSCP0927-Q
- ☐ Footplate Connector \$16
Part: RIGASMB7120023

6. Stealth Gel Footplate Covers

- ☐ 5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K0108
- ☐ 5"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K0108
- ☐ 5"x11" Tapered Gel Footplate Covers \$270
Part: ST-AFPG5011T-Q. HCPCS: K0108
- ☐ 5"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-AFPG5011E-Q. HCPCS: K0108

16.1d. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1a. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

- ☐ ☐ ☐ 60°-90° Swing-Away Legrest, Static Seat \$88 ea.
Pair part: XRF000082. Left: XRF000083. Right: XRF000084
Not available with 11"-13" extensions.
- ☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short, Static Seat \$88 ea.
Pair part: XRF000085. Left: XRF000086. Right: XRF000087
For use only with 11"-13" lower extensions.

Swing-Away Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Extruded, 13" - 16.5"	!	✓	✓
Extruded, 15.5" - 19.5"	!	!	!
Angle Adjustable, 13" - 16.5"	!	✓	✓
Angle Adjustable, 15.5" - 19.5"	!	!	!
Multi-Axis, 15.5" - 19.5"	!	!	!

☒ This rigging configuration causes no interference issues at this seat-to-floor height.

☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

☒ This rigging configuration does not work with this seat-to-floor height.

Legrests section continued on following pages.

16. LEGRESTS CONTINUED**2. Lower Extension**

☐ ☐ ☐ **11"-13" Length**
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.

☐ ☐ ☐ **13"-16.5" Length**
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.

☐ ☐ ☐ **15.5"-19.5" Length**
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ **5"Wx8"D Footplates with Heel Loops** \$92 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

☐ ☐ ☐ **6"Wx11"D¹ Footplates with Heel Loops** \$92 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

☐ ☐ ☐ **6"Wx8"D¹ Footplates with Heel Loops** \$92 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

☐ ☐ ☐ **8"Wx11"D² Footplates with Heel Loops** \$92 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

☐ **Therafin 31515 Multi-Axis Amputee Support** \$259 ea.
 QTY Part: INDPART2822. HCPCS: E1020

☐ **Therafin 32387 Amputee Support Gel Cover** \$183 ea.
 QTY Part: INDPART2823. HCPCS: K0108

16.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ **Elevating Legrest Assembly** \$171 ea.
 Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

☐ ☐ ☐ **11"-14" Length**
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

☐ ☐ ☐ **14"-17" Length**
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

☐ ☐ ☐ **8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length**
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

☐ ☐ ☐ **13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length**
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Standard, 11"-14"	✓	✓	✓
Standard, 14"-17"	!	!	!
Angle Adjustable, 12.5"-16.5"	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!

☒ This rigging configuration causes no interference issues at this seat-to-floor height.
☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2b. Residual Limb Support

Do not select a footplate or accessory.

☐ **Left Residual Limb Support with Calf Pad** \$342
 Part: FRMASMB5301. HCPCS: E1020

☐ **Right Residual Limb Support with Calf Pad** \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ **7.25"Wx6"D Reinforced Footplate** Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ **8.5"Wx6"D Reinforced Footplate** Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18".

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ **5"Wx8"D Footplates with Heel Loops** \$92 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

☐ ☐ ☐ **6"Wx11"D¹ Footplates with Heel Loops** \$92 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

☐ ☐ ☐ **6"Wx8"D¹ Footplates with Heel Loops** \$92 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

☐ ☐ ☐ **8"Wx11"D² Footplates with Heel Loops** \$92 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

☐ **Falcon Gel Padded Knee Adductor Buttons (KBE) (pair)** \$134
 Part: INDPART2242. HCPCS: E0953

16.2. Legrest Accessories

☐ **Heel Loops (pair)** \$54
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

☐ **Wheel Bumper** \$38
 Part: XRF000100

Legrests section continued on following pages.

16. LEGRESTS CONTINUED**16.3. Set Legrest Extension Length****REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.****1a. Centermount Extension Length**

- | | | | | | | | | | |
|--|--|--|---|---|---|---|---|---|---|
| ☐ 2" L
(#3989) | ☐ 5" L
(#3992) | ☐ 8" L
(#3680) | ☐ 9.5" L
(#3683) | ☐ 11" L
(#3686) | ☐ 12.5" L
(#3689) | ☐ 14" L
(#3692) | ☐ 15.5" L
(#3695) | ☐ 17" L
(#3698) | ☐ 18.5" L
(#3701) |
| ☐ 3" L
(#3990) | ☐ 6" L
(#3993) | ☐ 8.5" L
(#3681) | ☐ 10" L
(#3684) | ☐ 11.5" L
(#3687) | ☐ 13" L
(#3690) | ☐ 14.5" L
(#3693) | ☐ 16" L
(#3696) | ☐ 17.5" L
(#3699) | ☐ 19" L
(#3702) |
| ☐ 4" L
(#3991) | ☐ 7" L
(#3994) | ☐ 9" L
(#3682) | ☐ 10.5" L
(#3685) | ☐ 12" L
(#3688) | ☐ 13.5" L
(#3691) | ☐ 15" L
(#3694) | ☐ 16.5" L
(#3697) | ☐ 18" L
(#3700) | |

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

- | | | | | | | | | | | |
|--|--|---|---|---|---|---|---|---|---|---|
| ☐ 7" L
(#3470) | ☐ 8.5" L
(#3473) | ☐ 10" L
(#3476) | ☐ 11.5" L
(#3479) | ☐ 13" L
(#3482) | ☐ 14.5" L
(#3485) | ☐ 16" L
(#3488) | ☐ 17.5" L
(#3491) | ☐ 19" L
(#3494) | ☐ 20.5" L
(#3497) | ☐ 22" L
(#3669) |
| ☐ 7.5" L
(#3471) | ☐ 9" L
(#3474) | ☐ 10.5" L
(#3477) | ☐ 12" L
(#3480) | ☐ 13.5" L
(#3483) | ☐ 15" L
(#3486) | ☐ 16.5" L
(#3489) | ☐ 18" L
(#3492) | ☐ 19.5" L
(#3495) | ☐ 21" L
(#3498) | ☐ 22.5" L
(#3670) |
| ☐ 8" L
(#3472) | ☐ 9.5" L
(#3475) | ☐ 11" L
(#3478) | ☐ 12.5" L
(#3481) | ☐ 14" L
(#3484) | ☐ 15.5" L
(#3487) | ☐ 17" L
(#3490) | ☐ 18.5" L
(#3493) | ☐ 20" L
(#3496) | ☐ 21.5" L
(#3600) | |

17. SYNERGY POWER OPTION ELECTRONICS**REQUIRED for Synergy Seat with Power Legrests or Elevating Seat. Select the desired control type for the seating system.**

- | | |
|---|---|
| ☐ Single Actuator Through Controller Interface\$2,206
Part: XRF000335. HCPCS: E2310 | ☐ Multiple Actuator Through Controller and Interface..... \$3,101
Part: XRF000336. HCPCS: E2311 |
| ☐ Single Actuator Through Controller and Interface (AAM)..... \$3,101
Part: XRF000068. HCPCS: E2310 | |

18. SEAT MOUNTED LED LIGHTS**Optional.** Requires QLNE5 10 Key, Q-Logic 3 Joystick for Full Light Package, or Q-Logic 3e Joystick.

- | |
|--|
| ☐ Seat Mounted, LED, Full Light Package \$1,045
Part: XRF000130
Not available with Base Mounted, LED, Full Light Package. |
|--|

19. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

- | | |
|--|--|
| ☐ Synergy Seat Cup Holder No Charge
Part: ACCASMB2634 | ☐ Grip Drink Holder (20085) \$15.50
Part: ACC152993 |
| ☐ Oxygen Holder for Synergy Seat.....\$161
Not available with Synergy Manual Recline.
Part: FRMASMB7030. HCPCS: E2208 | ☐ Accessory Bar\$132
Part: XRF000128 |
| ☐ Walker Holder \$111
Synergy part: FRMASMB7031 | ☐ Large, Black, Weather Cover \$202
Part: ACCCOVER1021 |
| ☐ Cane/Crutch Holder\$118
Synergy part: FRMASMB4822. HCPCS: E2207 | ☐ Privacy Bag.....\$43
Part: ACC1907570 |
| ☐ Rear Basket\$126
Synergy part: FRMASMB4821 | ☐ Personal Item Pouch\$86
Part: ACC2101017 |
| ☐ Flag Holder\$55
Synergy/Flex part: ACCASMB2415. | ☐ Transfer Handles (pair) \$214
Part: XRF000396. HCPCS: E0705
Available with Flip-Up armrests only. |
| ☐ Saddle Bag.....\$40
Part: DWR1020D010 | ☐ Remote Stop Switch \$1,321
Part: XRF000379 |
| ☐ Stealth Cell Phone/Media Device Holder \$255
Part: ACC7120006 | |
| ☐ Stealth Large Tablet Holder \$625.28
Part: ST-URPKTAB0000-Q | |
| ☐ Quantum Backpack.....\$171
Part: ACC1905853 | |

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

- | | |
|---|--|
| ☐ Aftermarket Accessory Mounting Bracket, Slotted \$54 ea.
QTY Part: ACC2407046 | ☐ Aftermarket Accessory Mounting Bracket, Non-slotted .. \$54 ea.
QTY Part: ACC2405496 |
|---|--|

20. BRANCH PATCHES

Optional.

- | | | |
|---|---|--|
| ☐ Air Force Patch\$21
Part: ACC140240 | ☐ Coast Guard Patch\$21
Part: ACC140239 | ☐ Navy Patch\$21
Part: ACC140244 |
| ☐ Army Patch\$21
Part: ACC140241 | ☐ Marines Patch\$21
Part: ACC140243 | ☐ Purple Heart Patch\$21
Part: ACC140238 |



21. POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

21.1. Stealth Thoracic Laterals

Swing-away laterals may interfere with 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

- ☐ Swing-away, UniLink, Back Cane, pair (ULPK404104B00)\$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034

2. Pads

- | | |
|--|--|
| ☐ Curved Pads, 3"x5", pair (TWBL-3X5)..... \$250
Part: INDPART2663. HCPCS: E0956 | ☐ Flat Pads, 3"x5", pair (ULP846050320)..... \$250
Part: ST-ULP846050320-Q. HCPCS: E0956 |
| ☐ Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250
Part: INDPART2664. HCPCS: E0956 | ☐ Flat Pads, 4"x6", pair (ULP835050320)..... \$250
Part: ST-ULP835050320-Q. HCPCS: E0956 |

3. Installation

- ☐ Install Laterals
Part: PTOINDV3499
- ☐ Do Not Install
Part: PTOINDV3500

21.2. Stealth Unilink Pelvic/Thigh Guide

1. UniLink Pelvic/Thigh Hardware

^ Please verify the appropriate billing code depending upon the placement of the support pad.

- | | |
|---|---|
| ☐ UniLink Flip Down Hip/Thigh Support, Left..... \$332 ea.
Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^ | ☐ UniLink Flip Down Hip/Thigh Support, Right \$332 ea.
Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^ |
| ☐ Flip Down with LEFT Full Surface Contact \$281
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^ | ☐ Flip Down with RIGHT Full Surface Contact..... \$281
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^ |
| ☐ UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^ | ☐ UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^ |

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used for this item with all third-party payors.

- | | |
|---|--|
| ☐ UniLink, 4" x 4" Pad \$106 ea.
Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^ | ☐ UniLink, 4" x 12" Pad \$106 ea.
Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^ |
| ☐ UniLink, 4" x 8" Pad..... \$106 ea.
Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^ | |

3. UniLink Pelvic/Thigh Additional Links

- | | |
|--|---|
| ☐ UniLink, 1.5" Link \$81 ea.
Part: ST-UIL001150A000-Q | ☐ UniLink, 2" Link..... \$81 ea.
Part: ST-UIL001200A000-Q |
|--|---|

4. UniLink Pelvic/Thigh Guide Installation

- ☐ Install on Left
Part: PTOINDV3585
- ☐ Install on Right
Part: PTOINDV3586
- ☐ Do Not Install
Part: PTOINDV3719

21.4. Stealth Positioning Chest Harness

- | | |
|--|---|
| ☐ Upper Classic, Medium (IFPS1103) \$154 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960 | ☐ Upper Contour, Medium (IFPS2103)..... \$154 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960 |
| ☐ Upper Classic, Large (IFPS1104) \$154 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960 | ☐ Upper Contour, Large (IFPS2104) \$154 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960 |



Positioning components section continued on following page.

21. POSITIONING COMPONENTS CONTINUED**21.5. Stealth Positioning Chest Straps**

- ☐ Chest Strap, Medium w/ Belt Anchor (IFCS4113).. \$97 ea.
Part: ST-IFCS4113-Q. HCPCS: E0960
- ☐ Chest Strap, Large w/ Belt Anchor (IFCS4114)..... \$97 ea.
Part: ST-IFCS4114-Q. HCPCS: E0960

**21.6. Stealth Swing-Away Headrest Mounts**


These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

- ☐ UniLink Facial Lateral, Swing Away, Right..... \$318 ea.
Part: ST-SAR2A-R-Q. HCPCS: E1028
- ☐ UniLink Facial Lateral, Swing Away, Left \$318 ea.
Part: ST-SAR2A-L-Q. HCPCS: E1028
- ☐ UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028
- ☐ UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028
- ☐ UniLink Facial Lateral, Dual Swing Away \$434 ea.
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028

**2. Facial Lateral Pads and Accessories**

- ☐ Egg Switch Mount (ESM) (A)\$57 ea.
Part: INDPART2815
- ☐ Facial Pad, 4" 5x2", Cool Core (900M-Q)\$61 ea.
Part: ST-900M-Q. HCPCS: K0108
- ☐ Facial Pad, 6"x2", Cool Core (920M-Q).....\$61 ea.
Part: ST-920M-Q. HCPCS: K0108
- ☐ Swing-Away Joystick Mount,..... \$276 ea.
18" Rod (C)
Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest

**21.7. Stealth Elbow Support****1. Elbow Support Bracket**

- ☐ Left, Swing-Away, Elbow Support Bracket\$199 ea.
Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028
- ☐ Right, Swing-Away, Elbow Support Bracket.....\$199 ea.
Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

- ☐ 3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K0108
- ☐ 4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.
Part: POS124067. HCPCS: K0108

21.8. Stealth Foot Boxes

- ☐ Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954
- ☐ Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954
- ☐ Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

**21.9. Stealth Positioning Shoe Holder**

- ☐ Small Shoe Holder, \$87 ea.
3/4" Webbing (IFSH4211)
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952
- ☐ Large Shoe Holder,\$87 ea.
3/4" Webbing (IFSH4213)
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952
- ☐ Medium Shoe Holder, \$87 ea.
3/4" Webbing (IFSH4212)
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952



22. LAP TRAY

Optional.

22.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

- ☐ Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q, HCPCS: E0950
- ☐ Simple Solutions Large Tray with Cut Out..... \$309
Part: ST-SSLXL010001171-Q, HCPCS: E0950
- ☐ Simple Solutions Medium Tray with Cut Out..... \$309
Part: ST-SSLXM010001171-Q, HCPCS: E0950

2. Mounting Hardware

- ☐ Swing Away Mounting Hardware.....\$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	S	M	M	M	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	S	M	M	M	L	L	
Drop-In, T-Arm	S	S	S	S	S	S	M	M	M	L	L	
Two-Post, Flip-Back	S	S	S	S	S	S	M	M	M	L	L	
Cantilever	S	S	S	S	S	M	M	M	M	L	L	
Q-Captain, Flip-Up						S	M	M		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	M	M	M	M	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

:Compatible

:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

23. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Select a type of specialty controls below or to refer to the specialty controls order form.

23.1. Enhanced Display

- ☐ Q-Logic 3 Enhanced Display Kit..... \$2,020
Part: CTL167726
Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

- ☐ Black, Soft Touch Standard ☐ Night Train..... \$64 ☐ Pretty in Pink \$64 ☐ Lemon Crush \$64
Part: DGN137728 Part: DGN137729-PBK3 Part: DGN137729-PPK3 Part: DGN137729-PYL4
- ☐ Red Corvette..... \$64 ☐ Grey Street \$64 ☐ Midnight Blue \$64 ☐ Green Machine \$64
Part: DGN137729-PRD1 Part: DGN137729-PSL3 Part: DGN137729-PBU10 Part: DGN137729-PGN11
- ☐ Back in Black..... \$64 ☐ Deep Purple \$64 ☐ Electric Blue \$64 ☐ Military Green \$64
Part: DGN137729-PBK1 Part: DGN137729-PPR5 Part: DGN137729-PBU9 Part: DGN137729-PGN8

**23.2. Q-Logic 3 SAJ**

- ☐ Q-Logic 3 EX SAJ \$857
Part: CTL143980
Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.

**23.3. Head Array**

- ☐ i-Drive 4.0 Bluetooth Tri-array Series,..... \$4,774
3 Switch with Egg Switch
Part: ST-IDHBT200-U-Q, HCPCS: E2330
- ☐ i-Drive 4.0 Bluetooth Pro Series,..... \$4,825
3 Switch with Egg Switch
Part: ST-IDHBT300-U-Q, HCPCS: E2330
- ☐ Adult, Curved, i-Drive Tri-Array with \$4,474
Mounting Hardware
Part: ST-IDHBT400-U-Q, HCPCS: E2330
- ☐ Junior, Curved, i-Drive Tri-Array with \$4,474
Mounting Hardware
Part: ST-IDHBT450-U-Q, HCPCS: E2330



Alternative controls and electronics section continued on following page.

23. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

23.4. Egg Switches

- | | |
|---|--|
| ☐ <u>QTY</u> Red Egg Switch (ESRED).....\$95 ea.
Part: SWTMCRO1077 | ☐ <u>QTY</u> Green Egg Switch (ESGRN)\$95 ea.
Part: SWTMCRO1078 |
| ☐ <u>QTY</u> Blue Egg Switch (ESBLU)\$95 ea.
Part: SWTMCRO1079 | ☐ <u>QTY</u> Yellow Egg Switch (ESYEL)\$95 ea.
Part: SWTMCRO1080 |
| ☐ <u>QTY</u> Black Egg Switch (ESBLK)\$95 ea.
Part: SWTMCRO1081 | |



23.5. Switch Adapter

- ☐ QTY 2 to 1 Switch Adapter..... \$20
Part: ELE2501711
Connects 2 single switches to a stereo jack to act as a dual switch.

23.6. Environmental Controls

- ☐ QTY Q-Logic 3 ECU Module Kit..... \$2,570
Part: CTL168057

