



Q6 Edge® 2.0 X Order Form with Synergy Seating

Account Number: _____ Date: _____
 Provider Name: _____
 Contact: _____
 Phone: _____ Fax: _____
 Email: _____
 PO Number: _____
 Quote Name/Marked for: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____

Q6 edge[®]
2.0x



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
 B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
 C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a power base. This power base has a 300 lb. weight capacity with standard 6 mph programmable motors, 8 amp off board charger, fender lights and USB mobile device charger. Dimensions of base without legrests: 24"W x 35.5"L

Quantum Q6 Edge 2.0 X 3S-SS \$10,195
 Part: Q6EDGE 2.0X 3S-SS. HCPCS: K0848
 Solid seat pan standard with no power option.

Quantum Q6 Edge 2.0 X 3SP-SS \$10,195
 Part: Q6EDGE 2.0X 3SP-SS. HCPCS: K0856
 Single power option for Power Elevating Seat for Synergy Seat.

2. DRIVE MOTORS













REQUIRED.

6 mph Drive Motors Standard
 Part: MOT147433

Accu-Trac Motor Technology \$1,360
 Part: MOT147881. HCPCS: K0108
 Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

- | | | | |
|--|---|---|---|
|  Red Corvette
Part: DGN153201-PRD1 |  Night Train
Part: DGN153201-PBK3 |  Deep Purple
Part: DGN153201-PPR5 |  Lemon Crush
Part: DGN153201-PYL4 |
|  Electric Blue
Part: DGN153201-PBU9 |  Pretty in Pink
Part: DGN153201-PPK3 |  Grey Street
Part: DGN153201-PSL3 |  Green Machine
Part: DGN153201-PGN11 |
|  Back in Black
Part: DGN153201-PBK1 |  Evolution Orange
Part: DGN153201-POR7 |  Midnight Blue
Part: DGN153201-PBU10 |  Military Green
Part: DGN153201-PGN8 |

4. POWER ADJUSTABLE SEAT HEIGHT

Optional. The seat to floor height is 15.875" with 6" stroke and 16.875" with 7" stroke. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit.

- | | |
|--|--|
| Power Elevating Seat with 7" Stroke \$1,990
Part: XRF000131. HCPCS: E2300 | Power Elevating Seat with 6" Stroke \$1,990
Part: XRF000132. HCPCS: E2300 |
|--|--|

5. QUANTUM BACKUP CAMERA

Optional. Requires Q-Logic 3 electronics. The backup camera may interfere with armrest mounted accessories such as media device holder.

- Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

- | | |
|---|--|
| Backup Camera Headrest Bracket \$21
Part: MEC1910373
Works with .5" x .5" vertical supports only. | Backup Camera Track Mount Bracket \$27
Part: MEC2104227 |
|---|--|

2. Optional Display Mounting Bracket

- | | |
|---|--|
| Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123 | Backup Camera NE/NE+ Joystick Mount Bracket \$105
Part: MEC2106132 |
|---|--|



6. WHEEL OPTIONS

REQUIRED. Select a type of drive wheels and accent color if applicable. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

- | | |
|---|---|
| 3"x14" Flat Free Black Tires (No Accent Color) Standard
Part: PTOINDV3525 (pair) | 3"x14" Pneumatic Black Tires (No Accent Color) No Charge
Part: WHL150034 x2 |
| 3"x14" Flat Free Black Tires for Wheel Accents No Charge
Part: WHL1611297 x2 | 3"x14" Pneumatic Black Tires for Wheel Accents No Charge
Part: WHL1611296 x2 |
| 4"x14" Knobby, Flat Free Black Tires (No Accent Color) \$290 ea.
Part: WHL2006647 x2 | 4"x14" Knobby, Pneumatic Black Tires (No Accent Color) \$290 ea.
Part: WHL2006822 x2 |

6.2. Drive Wheel Accent Color

Optional. Not available with knobby tires.

- | | | |
|--|---|---|
| Back in Black Accent Color \$63
Part: XRF000179 | Deep Purple Accent Color \$63
Part: XRF000184 | Electric Blue Accent Color \$63
Part: XRF000188 |
| Red Corvette Accent Color \$63
Part: XRF000180 | Pretty in Pink Accent Color \$63
Part: XRF000185 | Lemon Crush Accent Color \$63
Part: XRF000189 |
| Night Train Accent Color \$63
Part: XRF000181 | Evolution Orange Accent Color \$63
Part: XRF000186 | Green Machine Accent Color \$63
Part: XRF000190 |
| Grey Street Accent Color \$63
Part: XRF000182 | Midnight Blue Accent Color \$63
Part: XRF000187 | Military Green Accent Color \$63
Part: XRF000191 |

7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

- | | |
|--|---|
| NF-22, Gel Cell, Interceptor Battery \$341 ea.
Part: BAT2004376. HCPCS: E2361 | NF-22, Gel Cell, MK Battery \$415 ea.
Part: BATGEL1008. HCPCS: E2361 |
| NF-22, AGM, MK Battery \$415 ea.
Part: BAT121346. HCPCS: E2361 | |

8. SEAT TOWERS

Optional. Not available with Power Adjustable Seat Height.

- Extended Seat Towers No Charge
Part: MEC120301
Increases seat to floor height.

9. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

9.1. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired.

1. Q-Logic 3 Joystick

REQUIRED.

Q-Logic 3 EX Joystick Standard
Part: XRF000241

No Joystick, Setup for Q-Logic 3
Part: XRF000010
Requires selection of Enhanced Display.

Q-Logic 3e Joystick No Charge
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.



Q-Logic 3 EX Joystick

2. Q-Logic 3 Required Components

REQUIRED with 3SP-SS.

Q-Logic 3 EX Controller Module \$730
Part: XRF000002. HCPCS: E2377

Q-Logic 3 Expandable Controller Harness..... \$447
Part: XRF000264. HCPCS: E2313

9.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft Touch Standard
Part: DGN124307
Also matches Night Train shrouds.

Red Corvette..... \$63
Part: DGN124336-PRD1

Back in Black..... \$63
Part: DGN124336-PBK1

Grey Street..... \$63
Part: DGN124336-PSL3

Deep Purple..... \$63
Part: DGN124336-PPR5

Pretty in Pink..... \$63
Part: DGN124336-PPK3

Evolution Orange..... \$63
Part: DGN124336-POR7

Midnight Blue..... \$63
Part: DGN124336-PBU10

Electric Blue..... \$63
Part: DGN124336-PBU9

Lemon Crush..... \$63
Part: DGN124336-PYL4

Green Machine..... \$63
Part: DGN124336-PGN11

Military Green..... \$63
Part: DGN124336-PGN8

10. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323

Chin Cup (C)..... \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323

Medium Ball Handle (E)..... \$69
Part: ACCASMB1859. HCPCS: E2323

Soft Ball Handle (F)..... \$80
Part: ACCASMB1856. HCPCS: E2323

Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323

Bodypoint 3" Goal Post Handle \$104
Part: PLSKNCP1057. HCPCS: E2323

Bodypoint 4" Goal Post Handle \$104
Part: PLSKNCP1058. HCPCS: E2323

Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308



Joystick Handles

11. JOYSTICK MOUNTING BRACKETS


Optional. Select a mount to be used with a joystick selected in previous section.

11.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617

Left Side Mount
Part: PTOINDV1225

Right Side Mount
Part: PTOINDV1226



Swing Away Joystick

11.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$294
Part: ACC1908434. HCPCS: E1028

11.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1028

Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1028

Joystick mounting brackets section continued on following page.

11. JOYSTICK MOUNTING BRACKETS CONTINUED

11.1.d. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount \$368

Part: XRF000003. HCPCS: E1028

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount \$126


Part: XRF000004

Allows height adjustment of joystick.

Stealth Swing-Away/Flip Down Bracket \$250

Part: XRF000005. HCPCS: E1028

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.



Retract4



Multi-Axis Joystick Mount

11.2. Joystick Mount Location

Left Side Mount

Part: PTOINDV1225

Right Side Mount

Part: PTOINDV1226

12. SYNERGY SEATING SYSTEM

Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 22"W x 22"D.

12.1. Synergy Frame Selection

REQUIRED.

Synergy Seating System (For K0848 Only)..... Standard

Part: PTOINDV3746

Synergy Seating System (For K0856 Only)..... \$456

Part: XRF000244

12.2. Synergy Transit Kit

Optional. The Q6 Edge 2.0 is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

WC-19 Occupied Transit Option..... \$289

Part: KITASMB1830

12.3. Width

REQUIRED. Maximum seat width is 22". Part number is SETWDTH plus number following # below.

10" W (#1000)	11" W (#1001)	12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.4. Width Package

REQUIRED. Maximum seat width is 22".

10"-14" Width Package

Part: FRMASMB18011

14"-18" Width Package

Part: FRMASMB18012

16"-20" Width Package

Part: FRMASMB18013

20"-24" Width Package

Part: FRMASMB18014

12.5. Depth

REQUIRED. Maximum seat depth is 22". Part number is SETDPTH plus number following # below.

10" D (#1000)	11" D (#1001)	12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.6. Depth Package

REQUIRED. Maximum seat depth is 22".

10"-14" Depth Package

Part: FRMASMB6981

14"-18" Depth Package

Part: FRMASMB17823

15"-20" Depth Package

Part: FRMASMB7732

20"-24" Depth Package

Part: FRMASMB7733

12.7. Synergy Seat to Floor

REQUIRED.

Lowest Seat to Floor

Part: PTOINDV3558

Middle Seat to Floor

Part: PTOINDV3559

Highest Seat to Floor

Part: PTOINDV3560

Low	Med	High
Synergy Seat		
15.875"	16.875"	17.875"
Extended Towers		
Add 3 inches to above		
Seat to Floor Height		

Synergy seating section continued on following pages.

12. SYNERGY SEATING SYSTEM CONTINUED


12.8. Matching Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

Stealth TRU-Comfort 2 SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-TC2SPPXXYY. HCPCS: E2607/E2608


TRU-Comfort 2 SPP, Matching Seat Size
Reference part: XRF000116. MSRP: \$420/\$525



Stealth Solution SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-SOLSPPXXYY. HCPCS: E2607/E2608

Solution SPP, Matching Seat Size
Reference part: XRF000113. MSRP: \$394/\$452



Stealth Spectrum Foam SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-SPCFSPXXYY. HCPCS: E2607/E2608

Spectrum Foam SPP, Matching Seat Size
Reference part: XRF000112. MSRP: \$483/\$557



Stealth Premiere P Cushion

Positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-PREMPXXYY. HCPCS: E2605/E2606

Premiere P, Matching Seat Size
Reference part: XRF000122. MSRP: \$341/\$499



Stealth Zen SP Cushion

Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ZENSPXXYY. HCPCS: E2603/E2604

Zen SP, Matching Seat Size
Reference part: XRF000115. MSRP: \$331/\$446



Stealth Spectrum Gel SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-SPCGSPPXXYY. HCPCS: E2607/E2608

Spectrum Gel SPP, Matching Seat Size
Reference part: XRF000117. MSRP: \$483/\$557



Stealth Essence SPP Cushion

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ESSSPPXXYY. HCPCS: E2607/E2608

Essence SPP, No Wedge, Matching Seat Size
Reference part: XRF000118. MSRP: \$420/\$525

Essence SPP, Small Wedge, Matching Seat Size
Reference part: XRF000119. MSRP: \$457/\$562

Essence SPP, Medium Wedge, Matching Seat Size
Reference part: XRF000120. MSRP: \$488/\$593


Essence SPP, Large Wedge, Matching Seat Size
Reference part: XRF000121. MSRP: \$520/\$625



Stealth Glacial SP Cushion

Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-GLSPXXYY. HCPCS: E2603/E2604


Glacial SP, Matching Seat Size
Reference part: XRF000114. MSRP: \$347/\$462



Stealth Simplicity G Cushion

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.
Part series: ST-SIMGXXYY. HCPCS: E2601/E2602

Simplicity G, Matching Seat Size
Reference part: XRF000111. MSRP: \$116/\$173



12.9. Synergy Seating Back Cane Selection

REQUIRED.

Synergy Back Cane No Charge
Part: FRMASMB15844
Select a back cane height.

TRU-Comfort Back Canes No Charge
Part: XRF000133

Flip Forward Back Cane \$184
Part: FRMASMB15846
Select a back cane height. Not available with a vent tray, oxygen holder, or cane mounted and cantilever armrests.

2" Mesomorph (1" per side) Back Cane \$205
Part: FRMASMB12067
Select a cane height.

2" Endomorph (1" per side) Back Cane \$205
Part: FRMASMB12115
Select a cane height.

4" Endomorph (2" per side) Back Cane \$205
Part: FRMASMB12139
Select a cane height.

Manual Recline Back \$788
Part: FRMASMB12910. HCPCS: E1226
Select an overall back height. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, Shape Back, TRU-Comfort Back, rear accessories, power elevating seat, or articulating foot platform. Pivot point is 2.5" above the seat pan.

12.10. Synergy Seating Back Cane Height

REQUIRED EXCEPT if a TRU-Comfort back is desired. 16", 17", and 18" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles.

14" High Canes
Part: FRMASMB15832

15" High Canes
Part: FRMASMB15833

16" High Canes
Part: FRMASMB15834

17" High Canes
Part: FRMASMB15835

18" High Canes
Part: FRMASMB15836

19" High Canes
Part: FRMASMB15837

20" High Canes
Part: FRMASMB15838

21" High Canes
Part: FRMASMB15839

22" High Canes
Part: FRMASMB15840

23" High Canes
Part: FRMASMB15841

24" High Canes
Part: FRMASMB15842

25" High Canes
Part: FRMASMB15843

Synergy seating section continued on following pages.

12. SYNERGY SEATING SYSTEM CONTINUED

12.11. Synergy Seating Back Cushion

REQUIRED. Select one type of back.

1a. Synergy Back Selection

Synergy Back Cushion \$289

Select size below. Height is overall. Part number is SETASMB plus number following #.

12"Wx15"H Part: #6448	14"Wx21"H Part: #6491	15"Wx22"H Part: #6668	16"Wx23"H Part: #6512	18"Wx17"H Part: #6465	19"Wx18"H Part: #6642	20"Wx19"H Part: #6478	21"Wx22"H Part: #6674	22"Wx23"H Part: #6518
12"Wx16"H Part: #6624	14"Wx22"H Part: #6667	15"Wx23"H Part: #6511	16"Wx24"H Part: #6688	18"Wx18"H Part: #6641	19"Wx19"H Part: #6477	20"Wx20"H Part: #6654	21"Wx23"H Part: #6517	22"Wx24"H Part: #6694
13"Wx17"H Part: #6460	14"Wx23"H Part: #6510	15"Wx24"H Part: #6687	16"Wx25"H Part: #6531	18"Wx19"H Part: #6476	19"Wx20"H Part: #6653	20"Wx21"H Part: #6497	21"Wx24"H Part: #6693	22"Wx25"H Part: #6537
13"Wx18"H Part: #6636	14"Wx24"H Part: #6686	15"Wx25"H Part: #6530	17"Wx17"H Part: #6464	18"Wx20"H Part: #6652	19"Wx21"H Part: #6496	20"Wx22"H Part: #6673	21"Wx25"H Part: #6536	
14"Wx13"H Part: #6447	14"Wx25"H Part: #6529	16"Wx17"H Part: #6463	17"Wx18"H Part: #6640	18"Wx21"H Part: #6495	19"Wx22"H Part: #6672	20"Wx23"H Part: #6516	22"Wx17"H Part: #6469	
14"Wx14"H Part: #6623	15"Wx17"H Part: #6462	16"Wx18"H Part: #6639	17"Wx21"H Part: #6494	18"Wx22"H Part: #6671	19"Wx23"H Part: #6515	20"Wx24"H Part: #6692	22"Wx18"H Part: #6645	
14"Wx15"H Part: #6449	15"Wx18"H Part: #6638	16"Wx19"H Part: #6474	17"Wx22"H Part: #6670	18"Wx23"H Part: #6514	19"Wx24"H Part: #6691	20"Wx25"H Part: #6535	22"Wx19"H Part: #6480	
14"Wx16"H Part: #6625	15"Wx19"H Part: #6473	16"Wx20"H Part: #6650	17"Wx23"H Part: #6513	18"Wx24"H Part: #6690	19"Wx25"H Part: #6534	21"Wx19"H Part: #6479	22"Wx20"H Part: #6656	
14"Wx17"H Part: #6461	15"Wx20"H Part: #6649	16"Wx21"H Part: #6493	17"Wx24"H Part: #6689	18"Wx25"H Part: #6533	20"Wx17"H Part: #6467	21"Wx20"H Part: #6655	22"Wx21"H Part: #6499	
14"Wx18"H Part: #6637	15"Wx21"H Part: #6492	16"Wx22"H Part: #6669	17"Wx25"H Part: #6532	19"Wx17"H Part: #6466	20"Wx18"H Part: #6643	21"Wx21"H Part: #6498	22"Wx22"H Part: #6675	

1b. TRU-Comfort 4 Back Selection

See back availability chart for sizes. Back size is the product of seat width and overall back height selections.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width..... \$709 Part: TC4V10W21WBC. HCPCS: E2620	TRU-Comfort Back, Back Cane Option, 22" Width \$814 Part: TC4V22W32WBC. HCPCS: E2621
---	---

1c. ADI Back, AL Series, 4-Point, Active Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 14" Wide x 13" or 14" Overall Height \$592.25 Part: ALBR14LA4. HCPCS: E2613	Active Contour, 18" Wide x 21" or 22" Overall Height \$592.25 Part: ALBR18TA4. HCPCS: E2613
Active Contour, 14" Wide x 16" or 17" Overall Height \$592.25 Part: ALBR14MA4. HCPCS: E2613	Active Contour, 18" Wide x 23" or 24" Overall Height \$592.25 Part: ALBR18XA4. HCPCS: E2613
Active Contour, 16" Wide x 19" or 20" Overall Height \$592.25 Part: ALBR16SA4. HCPCS: E2613	Active Contour, 20" Wide x 19" or 20" Overall Height \$592.25 Part: ALBR20SA4. HCPCS: E2613
Active Contour, 16" Wide x 21" or 22" Overall Height \$592.25 Part: ALBR16TA4. HCPCS: E2613	Active Contour, 20" Wide x 21" or 22" Overall Height \$592.25 Part: ALBR20TA4. HCPCS: E2613
Active Contour, 16" Wide x 23" or 24" Overall Height \$592.25 Part: ALBR16XA4. HCPCS: E2613	Active Contour, 20" Wide x 23" or 24" Overall Height \$592.25 Part: ALBR20XA4. HCPCS: E2613
Active Contour, 18" Wide x 19" or 20" Overall Height \$592.25 Part: ALBR18SA4. HCPCS: E2613	

Synergy seating section continued on following pages.

12. SYNERGY SEATING SYSTEM CONTINUED

1d. ADI Back, AL Series, 4-Point, Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 14" Wide x 13" or 14" Overall Height..... \$700.40 Part: ALBR14LD4. HCPCS: E2615	Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$700.40 Part: ALBR18SD4. HCPCS: E2615
Deep Contour, 14" Wide x 16" or 17" Overall Height..... \$700.40 Part: ALBR14MD4. HCPCS: E2615	Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$700.40 Part: ALBR18TD4. HCPCS: E2615
Deep Contour, 16" Wide x 19" or 20" Overall Height..... \$700.40 Part: ALBR16SD4. HCPCS: E2615	Deep Contour, 20" Wide x 19" or 20" Overall Height \$700.40 Part: ALBR20SD4. HCPCS: E2615
Deep Contour, 16" Wide x 21" or 22" Overall Height..... \$700.40 Part: ALBR16TD4. HCPCS: E2615	Deep Contour, 20" Wide x 21" or 22" Overall Height..... \$700.40 Part: ALBR20TD4. HCPCS: E2615
Deep Contour, 16" Wide x 23" or 24" Overall Height..... \$700.40 Part: ALBR16XD4. HCPCS: E2615	Deep Contour, 20" Wide x 23" or 24" Overall Height \$700.40 Part: ALBR20XD4. HCPCS: E2615

1e. ADI Back, AL Series, 4-Point, Extra Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 16" Wide x 21" or 22" Overall Height.. \$782.80 Part: ALBR16T64. HCPCS: E2615	Extra Deep Contour, 18" Wide x 23" or 24" Overall Height.. \$782.80 Part: ALBR18X64. HCPCS: E2615
Extra Deep Contour, 16" Wide x 23" or 24" Overall Height.. \$782.80 Part: ALBR16X64. HCPCS: E2615	Extra Deep Contour, 20" Wide x 19" or 20" Overall Height.. \$782.80 Part: ALBR20S64. HCPCS: E2615
Extra Deep Contour, 18" Wide x 19" or 20" Overall Height.. \$782.80 Part: ALBR18S64. HCPCS: E2615	Extra Deep Contour, 20" Wide x 21" or 22" Overall Height.. \$782.80 Part: ALBR20T64. HCPCS: E2615
Extra Deep Contour, 18" Wide x 21" or 22" Overall Height.. \$782.80 Part: ALBR18T64. HCPCS: E2615	Extra Deep Contour, 20" Wide x 23" or 24" Overall Height.. \$782.80 Part: ALBR20X64. HCPCS: E2615

1f. Sling Back Selection

Available 18"-25" overall heights. Back cane height must be 3" taller than sling back overall back height. Not available with cantilever arms.

Sling Back..... \$103

Select size below. Height is overall and must match back cane height. Part number is SETASMB plus number following #.

13"Wx18"H Part: #3740	14"Wx21"H Part: #3751	15"Wx22"H Part: #3761	16"Wx23"H Part: #3770	17"Wx24"H Part: #3779	18"Wx25"H Part: #3788	20"Wx18"H Part: #3797	21"Wx19"H Part: #3806	22"Wx20"H Part: #3815
13"Wx19"H Part: #3741	14"Wx22"H Part: #3752	15"Wx23"H Part: #3762	16"Wx24"H Part: #3771	17"Wx25"H Part: #3780	19"Wx18"H Part: #3789	20"Wx19"H Part: #3798	21"Wx20"H Part: #3807	22"Wx21"H Part: #3816
13"Wx20"H Part: #3742	14"Wx23"H Part: #3753	15"Wx24"H Part: #3763	16"Wx25"H Part: #3772	18"Wx18"H Part: #3781	19"Wx19"H Part: #3790	20"Wx20"H Part: #3799	21"Wx21"H Part: #3808	22"Wx22"H Part: #3817
13"Wx21"H Part: #3743	14"Wx24"H Part: #3754	15"Wx25"H Part: #3764	17"Wx18"H Part: #3773	18"Wx19"H Part: #3782	19"Wx20"H Part: #3791	20"Wx21"H Part: #3800	21"Wx22"H Part: #3809	22"Wx23"H Part: #3818
13"Wx22"H Part: #3744	14"Wx25"H Part: #3764	16"Wx18"H Part: #3765	17"Wx19"H Part: #3774	18"Wx20"H Part: #3783	19"Wx21"H Part: #3792	20"Wx22"H Part: #3801	21"Wx23"H Part: #3810	22"Wx24"H Part: #3819
13"Wx23"H Part: #3745	15"Wx18"H Part: #3756	16"Wx19"H Part: #3766	17"Wx20"H Part: #3775	18"Wx21"H Part: #3784	19"Wx22"H Part: #3793	20"Wx23"H Part: #3802	21"Wx24"H Part: #3811	22"Wx25"H Part: #3820
13"Wx24"H Part: #3746	15"Wx19"H Part: #3757	16"Wx20"H Part: #3767	17"Wx21"H Part: #3776	18"Wx22"H Part: #3785	19"Wx23"H Part: #3794	20"Wx24"H Part: #3803	21"Wx25"H Part: #3812	
13"Wx25"H Part: #3747	15"Wx20"H Part: #3758	16"Wx21"H Part: #3768	17"Wx22"H Part: #3777	18"Wx23"H Part: #3786	19"Wx24"H Part: #3795	20"Wx25"H Part: #3804	22"Wx18"H Part: #3813	
14"Wx18"H Part: #3748	15"Wx21"H Part: #3760	16"Wx22"H Part: #3769	17"Wx23"H Part: #3778	18"Wx24"H Part: #3787	19"Wx25"H Part: #3796	21"Wx18"H Part: #3805	22"Wx19"H Part: #3814	

Synergy seating section continued on following pages.

12. SYNERGY SEATING SYSTEM CONTINUED

1g. Omit Back Selection

Omit Back..... No Charge
Part: PTOINDV1287

12.12. Synergy Seating Overall Back Height

REQUIRED.

- | | |
|--|--|
| 12" Overall Back Height
Part: XRF000134 | 19" Overall Back Height
Part: XRF000141 |
| 13" Overall Back Height
Part: XRF000135 | 20" Overall Back Height
Part: XRF000142 |
| 14" Overall Back Height
Part: XRF000136 | 21" Overall Back Height
Part: XRF000143 |
| 15" Overall Back Height
Part: XRF000137 | 22" Overall Back Height
Part: XRF000144 |
| 16" Overall Back Height
Part: XRF000138 | 23" Overall Back Height
Part: XRF000145 |
| 17" Overall Back Height
Part: XRF000139 | 24" Overall Back Height
Part: XRF000146 |
| 18" Overall Back Height
Part: XRF000140 | 25" Overall Back Height
Part: XRF000147 |

Overall Height	Width									
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"
15"	S	S	-	-	-	-	-	-	-	-
16"	S	S	-	-	-	-	-	-	-	TC
17"	-	S	S	S	S	S	S	S	-	S TC
18"	-	S	S	S	S	S	S	S	-	S TC
19"	-	S	S	S	-	S	S	S	S	S TC
20"	-	S	S	S	-	S TC	S	S	S	S TC
21"	-	S	S	S	S	S TC	S	S	S	S TC
22"	-	S	S	S TC	S TC	S TC	S	S TC	S	S TC
23"	-	S	S	S TC	S TC	S	S	S TC	S	S TC
24"	-	S	S	S	S	S TC	S	S TC	S	S
25"	-	S	S	S	S	S TC	S	S TC	S	S

Synergy TRU-Comfort 4

Synergy and TRU-Comfort Back Availability Chart

12.13. Synergy Armrests

REQUIRED with Synergy seating. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

- | | |
|--|---|
| <input type="checkbox"/> L <input type="checkbox"/> R Full Length, Flip Back Armrests..... Standard
Left part: FRMASMB7038. Right part: FRMASMB7039 | <input type="checkbox"/> L <input type="checkbox"/> R Desk Length, Flip Back Armrests Standard
Left part: FRMASMB7040. Right part: FRMASMB7041 |
|--|---|

1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

- | | |
|---|--|
| <input type="checkbox"/> L <input type="checkbox"/> R Full Length, Two-Post Armrests.....\$158 ea.
Left part: XRF000051. Right part: XRF000052. HCPCS: E0973 | <input type="checkbox"/> L <input type="checkbox"/> R Desk Length, Two-Post Armrests\$158 ea.
Left part: XRF000053. Right part: XRF000054. HCPCS: E0973 |
|---|--|

1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

- | | |
|--|--|
| <input type="checkbox"/> L <input type="checkbox"/> R Full Length, Heavy Duty Armrests, 8"-12" H \$90 ea.
Left Part: XRF000245. Right part: XRF000246. HCPCS: E0973 | <input type="checkbox"/> L <input type="checkbox"/> R Desk Length, Heavy Duty Armrests, 8"-12"H..... \$90 ea.
Left part: XRF000247. Right part: XRF000248. HCPCS: E0973 |
|--|--|

1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

- | | |
|---|--|
| <input type="checkbox"/> L <input type="checkbox"/> R Full Length, Heavy Duty Armrests, 10"-14"H..... \$90 ea.
Left Part: XRF000249. Right part: XRF000250. HCPCS: E0973 | <input type="checkbox"/> L <input type="checkbox"/> R Desk Length, Heavy Duty Armrests, 10"-14"H \$90 ea.
Left part: XRF000251. Right part: XRF000252. HCPCS: E0973 |
|---|--|

1e. Pediatric, Removable, Height Adjustable Armrests

Only available with desk length armpads.

- | |
|---|
| <input type="checkbox"/> L <input type="checkbox"/> R Pediatric, Removable Armrests \$79 ea.
Left Part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973 |
|---|

1f. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

- | | |
|--|---|
| <input type="checkbox"/> L <input type="checkbox"/> R Full Length, Cantilever Armrests\$103 ea.
Left part: FRMASMB16188. Right part: FRMASMB16189 | <input type="checkbox"/> L <input type="checkbox"/> R Desk Length, Cantilever Armrests.....\$103 ea.
Left part: FRMASMB16190. Right part: FRMASMB16191 |
|--|---|

1g. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

- | | |
|--|--|
| <input type="checkbox"/> L <input type="checkbox"/> R Full Length, Single Post Armrests\$79 ea.
Left Part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973 | <input type="checkbox"/> L <input type="checkbox"/> R Desk Length, Single Post Armrests\$79 ea.
Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973 |
|--|--|

12.14. Synergy Armpads

REQUIRED with armrest selection above. Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

- | | |
|--|---|
| Straight Armpads..... Standard
Part: XRF000148 | Waterfall Armpads..... No Charge
Part: XRF000149 |
| No Armpad Selection No Charge
Part: PTOINDV1373 | |

12.15. Joystick Bracket Receiver

REQUIRED.

- | | |
|---|--|
| <input type="checkbox"/> Standard Fixed Receiver Standard
Part: FRMASMB16062 | <input type="checkbox"/> Stealth Height Adjustable Receiver \$225
Part: XRF000064. HCPCS: K0108 |
|---|--|

13. HEADRESTS

Optional. Select a headrest pad and bracket.

13.1a. Headrests, Comfort Plus with Cool Core

1. Headrest Pads

- | | |
|--|--|
| Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$170.25
Part: ST-CP08CCFPM-Q, HCPCS: E0955 | Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) \$170.25
Part: ST-CP14CCFPM-Q, HCPCS: E0955 |
| Stealth 10" Comfort Plus, Cool Core (CP10CCFPM) \$170.25
Part: ST-CP10CCFPM-Q, HCPCS: E0955 | Stealth 6", Comfort Plus, Cool Core (CP06CCFPM) \$170.25
Part: ST-CP06CCFPM-Q, HCPCS: E0955 |

2. Headrest Brackets

- Comfort Plus Headrest Mounting Bracket, UniLink..... \$296.50
Part: ST-UHB2D33200-Q, HCPCS: E1028

13.1b. Headrests, Comfort Plus and Stealth Contoured

1. Headrest Pads

- | | |
|---|--|
| Stealth 8" Comfort Plus (CP450)..... \$170.25
Part: POS148484, HCPCS: E0955 | Stealth 14" Comfort Plus (CP150) \$170.25
Part: INDPART2327, HCPCS: E0955 |
| Stealth 10" Comfort Plus (CP250)..... \$170.25
Part: SETHEAD1057, HCPCS: E0955 | Pediatric Comfort Plus (6"Wx4.5"H) (CPN350)..... \$170.25
Part: SETHEAD1069, HCPCS: E0955 |

2. Headrest Brackets

- Multi-Axis Removable Headrest Bracket (HMO475P) \$165.32
Part: HDWASMB1084, HCPCS: E1028

13.2. Flip-down Headrest Bracket

- Stealth Flip-down Headrest Bracket (FDM350)..... \$162.23
Part: INDPART2219
The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

14. LAP BELTS

REQUIRED. Select one lap belt below.

14.1a. Standard Lap Belts

- 38" Lap Belt Part: ACCBELT1006 50" Lap Belt Part: ACCBELT1000 60" Lap Belt Part: ACCBELT1002 70" Lap Belt Part: ACCBELT1003 80" Lap Belt Part: ACCBELT1004 90" Lap Belt Part: ACCBELT1001

14.1b. Optional Lap Belts

- | | |
|--|---|
| <input type="checkbox"/> Stealth Positioning 60" Single Pull, Push Button, Medium..... \$47.38
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q | <input type="checkbox"/> Stealth Positioning 60" Single Pull, Push Button, Extra Large \$47.38
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q |
| <input type="checkbox"/> Stealth Positioning 60" Single Pull, Push Button, Large \$47.38
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q | <input type="checkbox"/> 60" Sure-Lok Lap Belt (Non-retractable)..... \$120
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016 |

15. POWER LEGREST CONTROLS

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

- | | |
|---|---|
| Independent Legrests (adds 2 actuator control requirement)
Part: XRF000108 | Articulating Foot Platform (adds 1 actuator control requirement)
Part: XRF000110 |
| Combined Legrests (adds 2 actuator control requirement)
Part: XRF000109 | |

16. BASE MOUNTED FOOT PLATFORM

REQUIRED. May interfere with seat mounted legrests. Not available with center mount foot platforms.

- | | |
|---|--|
| Base Mounted Foot Platform Standard
Part: FRMASMB16077 | Omit base mounted foot platform No Charge
Part: PTOINDV3674 |
| Bell Shaped, Base Mounted Foot Platform \$247
Part: RIGASMB7120041 | |

17. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

17.1a. Basic Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform Assembly

Basic Center Mount..... No Charge
 Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates \$180 Long Extension, 5"x9" Footplates \$180
 Part: XRF000312. HCPCS: K0040 Part: XRF000314. HCPCS: K0040

Medium Extension, 5"x9" Footplates \$180
 Part: XRF000313. HCPCS: K0040

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount..... \$158 Large Calf Pads, Basic Center Mount \$158
 Part: RIG2104369 Part: RIG2202997
 Only available with Long Length Extension.

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16
 Part: RIG2104240

17.1b. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform Assembly

Power Articulating Foot Platform \$3,145
 Part: PTOINDV3754. HCPCS: E1012+K0040

2. Lower Extension and Angle Adjustable Footplates

Medium Extension, 5"x6" Footplates Long Extension, 5"x6" Footplates
 Part: RIG150206 Part: RIG151505

Medium Extension, 5"x8" Footplates Long Extension, 5"x8" Footplates
 Part: RIG150203 Part: RIG151506

Medium Extension, 5"x11" Footplates Long Extension, 5"x11" Footplates
 Part: RIG150212 Part: RIG151508

Medium	Long
11.5"-18"	14"-20"
Power AFP Extension Range	

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
 Part: RIGASMB7120018 Part: RIGASMB7120019

3b. Center Mount Calf Panel

Small Calf Panel for AFP, 18"Wx8"H..... \$184 Large Calf Panel for AFP, 24"Wx8"H \$237
 Part: FRMASMB7110007 Part: RIG123824
 Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$69 4" Outside Edge (pair) \$69
 Part: RIG152517 Part: RIG152510
 Only available with 6"D footplates and 8"D footplate extension.

Width Extension for 11"D Footplates \$69 6" Outside Edge (pair) \$69
 Part: RIG152518 Part: RIG152516
 Available with 8"D and 11"D footplates and 11"D footplate extension.

2" Bell-Shaped Length Extension \$69
 Part: RIG152519

5. Additional Accessories

Small Reversed Dartex Calf Pad Cover (pair) \$158 Small Foam Calf Pad Cover (pair)..... \$63.74
 Part: RIG133717 Part: ST-CSCS0775-Q

Large Reversed Dartex Calf Pad Cover (pair) \$210 Large Foam Calf Pad Cover (pair) \$78.72
 Part: RIG133719 Part: ST-CSCS0985-Q

Small Reversed Dartex Panel Cover (pair)..... \$210 Small Foam Panel Cover (pair) \$167.57
 Part: RIG134673 Part: ST-CSCS0920-Q

Large Reversed Dartex Panel Cover (pair) \$210 Large Foam Panel Cover (pair)..... \$268.77
 Part: RIG134744 Part: ST-CSCS0927-Q

Footplate Connector \$16
 Part: RIGASMB7120023

Legrests section continued on following pages.

17. LEGRESTS CONTINUED

6. Stealth Gel Footplate Covers

- 5"x8" Tapered Gel Footplate Covers \$211.15 5"x8" Tapered Gel Footplate Covers for Width Extension ... \$231.75
Part: ST-FPG5008T-Q. HCPCS: K0108 *Part: ST-FPG6508T-Q. HCPCS: K0108*
- 5"x11" Tapered Gel Footplate Covers \$231.75 5"x11" Tapered Gel Footplate Cover for Width Extension ... \$257.50
Part: ST-FPG5011T-Q. HCPCS: K0108 *Part: ST-FPG7011T-Q. HCPCS: K0108*

17.1c. Pediatric, Swing-Away Legrests

Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Pediatric Legrests Upper

- P L R 90°, Pediatric, Swing-Away Legrest (7"-10" Length Extension) \$172 ea.
Pair part: FRMASMB8808. Left: FRMASMB8809
Right: FRMASMB8810

2. Angle Adjustable Footplate(s)

- P L R 4"Wx6"D Angle Adjustable \$90 ea.
Pair part: FRMASMB8812. Left: RIG133537. Right: RIG133560
HCPCS: K0040
- P L R 5"Wx6"D Angle Adjustable \$90 ea.
Pair part: FRMASMB8811. Left: RIG133538. Right: RIG133561
HCPCS: K0040

17.1d. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1a. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

- P L R 60°-90° Swing-Away Legrest, Static Seat \$86 ea.
Pair part: XRF000082. Left: XRF000083. Right: XRF000084
Not available with 11"-13" extensions.
- P L R 60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea.
Pair part: XRF000085. Left: XRF000086. Right: XRF000087
For use only with 11"-13" lower extensions.

1b. 70° Swing-Away Legrest Uppers

Only available with occupied transit WC-19 selections. Sold only in pairs. Select lower extension and footplates as pairs.

- 70° Swing-Away Legrest (pair) \$172
Pair part: FRMASMB7841
Not available with 11"-13" extensions.

Swing-Away Legrest Compatibility Matrix	Contoured/Synergy Seat			
	Low	Medium	High	
Extruded, 13" - 16.5"	!	✓	✓	<input checked="" type="checkbox"/> This rigging configuration causes no interference issues at this seat-to-floor height.
Extruded, 15.5" - 19.5"	!	!	!	<input type="checkbox"/> This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
Angle Adjustable, 13" - 16.5"	!	✓	✓	
Angle Adjustable, 15.5" - 19.5"	!	!	!	
Multi-Axis, 15.5" - 19.5"	!	!	!	<input checked="" type="checkbox"/> This rigging configuration does not work with this seat-to-floor height.

2. Lower Extension

- P L R 11"-13" Length
Pair part: XRF000088.
Left: XRF000089.
Right: XRF000090.
N/A with Multi-axis.
- P L R 11"-16.5" Length
Pair part: XRF000091.
Left: XRF000092.
Right: XRF000093.
Multi-Axis only.
- P L R 13"-16.5" Length
Pair part: XRF000094.
Left: XRF000095.
Right: XRF000096.
N/A with Multi-Axis.
- P L R 15.5"-19.5" Length
Pair part: XRF000097.
Left: XRF000098.
Right: XRF000099.

3a. Extruded Footplate(s)

- P L R 7.25"Wx6"D Extruded Footplate Standard
Part pair: FRMASMB7847. Left: FRMASMB7861
Right: FRMASMB7862
- P L R 8.5"Wx6"D Extruded Footplate Standard
Part pair: FRMASMB7848. Left: FRMASMB7863
Right: FRMASMB7864
May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- P L R 4"Wx6"D \$90 ea.
Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560
HCPCS: K0040+E0951
- P L R 6"Wx8"D¹ \$90 ea.
Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563
HCPCS: K0040+E0951
- P L R 5"Wx6"D \$90 ea.
Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561
HCPCS: K0040+E0951
- P L R 6"Wx11"D¹ \$90 ea.
Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564
HCPCS: K0040+E0951
- P L R 5"Wx8"D \$90 ea.
Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562
HCPCS: K0040+E0951
- P L R 8"Wx11"D² \$90 ea.
Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565
HCPCS: K0040+E0951

Legrests section continued on following pages.

17. LEGRESTS CONTINUED

3c. Multi-Axis Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- | | |
|---|---|
| <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 4"Wx6"D \$260 ea.
 Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552
 HCPCS: K0040+E0951</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx6"D \$260 ea.
 Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553
 HCPCS: K0040+E0951</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx8"D \$260 ea.
 Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554
 HCPCS: K0040+E0951</p> | <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx8"D¹ \$260 ea.
 Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555
 HCPCS: K0040+E0951</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx11"D¹ \$260 ea.
 Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558
 HCPCS: K0040+E0951</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8"Wx11"D² \$260 ea.
 Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559
 HCPCS: K0040+E0951</p> |
|---|---|

4. Swing-Away Legrest Accessories

- Gel Neoprene Footrest Sleeve (pair) \$80
 Part: INDPART2243. HCPCS: K0108

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

- | | |
|---|--|
| <p>Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 Part: INDPART2822. HCPCS: E1020</p> | <p>Therafin 32387 Amputee Support Gel Cover \$155 ea.
 Part: INDPART2823. HCPCS: K0108</p> |
|---|--|

17.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

- Elevating Legrest Assembly** \$171 ea.
 Pair part: FRMASMB6483. Left: RIG138012
 Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2. Lower Extension

- 11"-14" Length**
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only
- 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length**
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only
- 14"-17" Length**
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only
- 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length**
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

2b. Residual Limb Support

Do not select a footplate or accessory.

- | | |
|--|---|
| <p>Left Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5301. HCPCS: E1020</p> | <p>Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020</p> |
|--|---|

3a. Reinforced Cast Footplate(s)

- | | |
|---|---|
| <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 7.25"Wx6"D Reinforced Footplate..... Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056</p> | <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8.5"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."</p> |
|---|---|

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

- | | |
|---|---|
| <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 4"Wx6"D \$90 ea.
 Pair parts: FRMASMB8812 x2.
 HCPCS: K0040</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx6"D \$90 ea.
 Pair parts: FRMASMB8811 x2.
 HCPCS: K0040</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx8"D Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951</p> | <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951</p> <p><input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8"Wx11"D² Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951</p> |
|---|---|

4. Elevating Legrest Accessories

- Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

Legrests section continued on following pages.

Elevating Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Standard, 11"-14"	✓	✓	✓
Standard, 14"-17"	!	!	!
Angle Adjustable, 12"-16.5"	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!

This rigging configuration causes no interference issues at this seat-to-floor height.
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

17. LEGRESTS CONTINUED

17.2. Legrest Accessories

Heel Loops (pair)..... \$53 Wheel Bumper..... \$19 ea.
 Part: FRMASMB7873. HCPCS: E0951 Part: XRF000100
 Standard on angle adjustable and multi axis footplates when noted.

17.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

18. SYNERGY POWER OPTION ELECTRONICS

REQUIRED with Power Legrests or Elevating Seat. Select the desired control type for the seating system.

Single Actuator Through Controller Interface \$2,163 Multiple Actuator Through Controller and Interface..... \$3,040
 Part: XRF000335. HCPCS: E2310 Part: XRF000336. HCPCS: E2311

19. SEAT MOUNTED LED LIGHTS

Optional. Requires Q-Logic 3 Joystick for Full Light Package, or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045
 Part: XRF000130
 Not available with Base Mounted, LED, Full Light Package.

20. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Synergy Seat Cup Holder No Charge Part: ACCASMB2634	Quantum Backpack \$168 Part: ACC1905853
Oxygen Holder for Synergy Seat..... \$158 Not available with Synergy Manual Recline. Part: FRMASMB7030. HCPCS: E2208	Grip Lap Board (20005)..... \$38.95 Part: ACC152992
Walker Holder \$109 Synergy part: FRMASMB7031	Grip Drink Holder (20085) \$15.50 Part: ACC152993
Cane/Crutch Holder \$116 Synergy part: FRMASMB4822. HCPCS: E2207	Accessory Bar \$132 Part: XRF000128
Rear Basket \$124 Synergy part: FRMASMB4821	Large, Black, Weather Cover \$202 Part: ACCCOVR1021
Flag Holder \$55 Synergy/Flex part: ACCASMB2415.	Q-Logic EX/NE Programmer \$573 Part: ELEASMB6976
Saddle Bag..... \$40 Part: DWR1020D010	Q-Logic PC Programming Disk with Cable \$445 Part: ELEASMB5216
Stealth Cell Phone/Media Device Holder \$205 Part: ACC7120006	Personal Item Pouch \$84 Part: ACC2101017
Stealth Large Tablet Holder \$400 Part: ACC7120005	Transfer Handle \$105 ea. Part: MEC2205680 x2. HCPCS: E0705

21. BRANCH PATCHES

Optional.

Air Force Patch \$21 Part: ACC140240	Coast Guard Patch \$21 Part: ACC140239	Navy Patch \$21 Part: ACC140244
Army Patch \$21 Part: ACC140241	Marines Patch \$21 Part: ACC140243	Purple Heart Patch \$21 Part: ACC140238



22. POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

22.1. Stealth Thoracic Laterals

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, Rail Mounts, pair (TWBSTL)\$535
 Part: INDPART2896. HCPCS: E1028

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$238 Part: INDPART2663. HCPCS: E0956	Flat Pads, 3"x5", pair (TWBLP3X5) \$238 Part: INDPART2676. HCPCS: E0956
Curved Pads, 4"x6", pair (TWBL-4X6)..... \$238 Part: INDPART2664. HCPCS: E0956	Flat Pads, 4"x6", pair (TWBLP4X6) \$238 Part: INDPART2675. HCPCS: E0956

3. Installation

Install Laterals Part: PTOINDV3499	Do Not Install Part: PTOINDV3500
---------------------------------------	-------------------------------------

4. Lateral Extensions

REQUIRED with Synergy Back and thoracic lateral selections.

Lateral Extension for Synergy Back (pair)..... \$28 ea.
 Part: XRF000150

22.2. Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW)\$267.50 ea. Part: INDPART2804. HCPCS: E1028 QTY _____	Flip Down with LEFT Full Surface Contact\$267.75 ea. Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028
Removable Bracket (LP430)..... \$270.38 ea. Part: INDPART2805. HCPCS: E1028 QTY _____	Flip Down with RIGHT Full Surface Contact.....\$267.75 ea. Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028
Removable Bracket with 1" Offset \$363.61ea. (LP430 & LHW-121) Part: INDPART2854. HCPCS: E1028 QTY _____	

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122)..... \$27.04 Part: INDPART2821. HCPCS: K0108 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.	2 1/4" Adapter Plate (LHW-129)\$54.08 Part: POS140430. HCPCS: K0108 Only available with the flip down full surface bracket.
---	---

3. Pelvic/Thigh Guide Pad

Short Pad, 4"x4" (TWBADD-4X4)..... \$77.25 ea. Part: INDPART2801. HCPCS: E0953/E0956 QTY _____	Long Pad, 4"x12" (TWBADD-4X12) \$77.25 ea. Part: INDPART2803. HCPCS: E0953/E0956 QTY _____
Medium Pad, 4"x8" (TWBADD-4X8)..... \$77.25 ea. Part: INDPART2802. HCPCS: E0953/E0956 QTY _____	

4. Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

22.3. Stealth Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$146.26 ea. Part: ST-IFPS1103-Q. HCPCS: E0960	Upper Contour, Medium (IFPS2103)..... \$146.26 ea. Part: ST-IFPS2103-Q. HCPCS: E0960
Upper Classic, Large (IFPS1104) \$146.26 ea. Part: ST-IFPS1104-Q. HCPCS: E0960	Upper Contour, Large (IFPS2104) \$146.26 ea. Part: ST-IFPS2104-Q. HCPCS: E0960



22.4. Stealth Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor (IFCS4113)\$92.70 ea. Part: ST-IFCS4113-Q. HCPCS: E0960	Chest Strap, Large w/ Belt Anchor (IFCS4114).\$92.70 ea. Part: ST-IFCS4114-Q. HCPCS: E0960
---	---



Positioning components section continued on following page.

22. POSITIONING COMPONENTS CONTINUED

22.5. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

Facial Lateral for Comfort Series.....\$303.08 ea.
 Right Side (CPS9-R-HW)
 Part: INDPART2647. HCPCS: E1028

Facial Lateral for Comfort Series.....\$303.08 ea.
 Left Side (CPS9-L-HW)
 Part: INDPART2648. HCPCS: E1028

Facial Laterals for Comfort Series \$358.18 ea.
 Left and Right Side (CPS9-B-HW)
 Part: INDPART2649. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Right..\$303.08 ea.
 Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left ... \$303.08 ea.
 Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away ...\$468.39 ea.
 2 Upper & 1 Lower Right
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away ...\$468.39 ea.
 2 Upper & 1 Lower Left
 Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away ... \$413.29 ea.
 2 Upper
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028



Left and Right Facial Lateral with Pads

2. Facial Lateral Pads and Accessories


Egg Switch Mount (ESM) (A) \$54.08 ea.
 Part: INDPART2815

Figure 8 Bracket with 13" Rod.....\$275.53 ea.
 and Egg Switch Mount (B)
 Part: ST-F8675-13-Q

Facial Pad, 4".5x2", Cool Core (900M-Q) . \$57.68 ea.
 Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2", Cool Core (920M-Q)....\$57.68 ea.
 Part: ST-920M-Q. HCPCS: K0108

Swing-Away Joystick Mount,.....\$262.50 ea.
 18" Rod (C)
 Part: ST-SUS18-Q
 Allows mounting of facial pads and attaches to the headrest



Headrest Accessories

22.6. Stealth Elbow Support

1. Elbow Support Bracket

Left, Swing-Away, Elbow Support Bracket\$189.26 ea.
 Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket.....\$189.26 ea.
 Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$54.08 ea.
 Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4).....\$54.08 ea.
 Part: POS124067. HCPCS: K0108

22.7. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$513.71
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$513.71
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$513.71
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954



Stealth Foot Boxes

22.8. Stealth Positioning Shoe Holder

Small Shoe Holder,\$82.40 ea.
 3/4" Webbing (IFSH4211)
 Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Medium Shoe Holder, \$82.40 ea.
 3/4" Webbing (IFSH4212)
 Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Large Shoe Holder, \$82.40 ea.
 3/4" Webbing (IFSH4213)
 Part: ST-IFSH4213-Q. HCPCS: E0951+E0952



Stealth Shoe Holders

23. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Select a type of specialty controls below or to refer to the specialty controls order form.

23.1. Enhanced Display

Q-Logic 3 Enhanced Display Kit..... \$1,980

Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft TouchStandard Part: DGN137728	Pretty in Pink \$63 Part: DGN137729-PPK3	Lemon Crush \$63 Part: DGN137729-PYL4
Red Corvette \$63 Part: DGN137729-PRD1	Evolution Orange \$63 Part: DGN137729-POR7	Green Machine \$63 Part: DGN137729-PGN11
Electric Blue \$63 Part: DGN137729-PBU9	Deep Purple \$63 Part: DGN137729-PPR5	Military Green \$63 Part: DGN137729-PGN8
Back in Black \$63 Part: DGN137729-PBK1	Grey Street \$63 Part: DGN137729-PSL3	
Night Train \$63 Part: DGN137729-PBK3	Midnight Blue \$63 Part: DGN137729-PBU10	

23.2. Q-Logic 3 SAJ

Q-Logic 3 EX SAJ \$840

Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

23.3. Head Array

i-Drive 4.0 Bluetooth Tri-array Series, 3 Switch with Egg Switch \$4,546.42

Part: ST-IDHBT200-U-Q. HCPCS: E2330

Junior, Curved, i-Drive Tri-Array with Mounting Hardware \$4,546.42

Part: ST-IDHBT450-U-Q. HCPCS: E2330

i-Drive 4.0 Bluetooth Pro Series, 3 Switch with Egg Switch \$4,594.83

Part: ST-IDHBT300-U-Q. HCPCS: E2330

Adult, i-Drive Pro Series Head Array (IDHBT300-1) \$4,594.83

Part: ST-IDHBT300-1-Q. HCPCS: E2330

Adult, Curved, i-Drive Tri-Array with Mounting Hardware \$4,546.42

Part: ST-IDHBT400-U-Q. HCPCS: E2330

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

23.4. Egg Switches

Red Egg Switch (ESRED)..... \$83.28 ea.

Part: SWTMCRO1077

Green Egg Switch (ESGRN) \$83.28 ea.

Part: SWTMCRO1078

Blue Egg Switch (ESBLU) \$83.28 ea.

Part: SWTMCRO1079

Yellow Egg Switch (ESYEL) \$83.28 ea.

Part: SWTMCRO1080

Black Egg Switch (ESBLK) \$83.28 ea.

Part: SWTMCRO1081



Egg Switch

23.5. Switch Adapter

2 to 1 Switch Adapter \$15

Part: DWR1361H037

Connects 2 single switches to a stereo jack to act as a dual switch.

23.6. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520

Part: CTL168057